



CATALOGUE



Circuit Components		Manual Resuscitators Accessories	
Reusable Polysulfone Connectors	1, 2	Reusable PEEP Valves	60
Reusable Polycarbonate Connectors	3, 4, 5	Disposable PEEP Valves	60
Disposable Single Patient Use Connectors	6, 7, 8, 9	PEEP Valve Diverters	60
Reusable Anaesthetic Components	10, 11		
Catheter Mounts		Oxygen Therapy Devices	
Reusable Catheter Mounts	12	Oxygen Therapy Masks	61
Disposable Catheter Mounts	13	Aerosol Masks	62
		Face Tent Masks	62
		Nebulizer Bottles	63
		Tracheostomy Masks	64
		Tracheostomy Mask Sets	64
		Oxygen Masks with Concentration Jets	65
		Venturi Diluter	65
		Child Venturi Masks	65
		Adult Venturi Masks	65
		Adjustable Diluter	65
		Nebulizer Kit	65
Nasal Cannulas		Laryngoscopes	
Oxygen Cannulas	14, 15	Duoscope - World's First Dual-blade Disposable Emergency Laryngoscope	66
Oxygen Cannulas with CO ₂ Sampling Line	14	Onescope - New Low-cost Disposable Emergency Laryngoscope	67
		Disposable Fibreoptic / LED Macintosh Blades	68
		Disposable Fibreoptic / LED Miller Blades	68
		Disposable Fibreoptic / LED Robertshaw Blades	68
		Disposable Fibreoptic / LED Kessell Blades	68
		Disposable Fibreoptic / LED Seward Blades	68
		Disposable Standard Macintosh Blades	69
		Disposable Standard Miller Blades	69
		Reusable Macintosh Blades	70
		Reusable Miller Blades	70
		Reusable Wisconsin (Foregger) Blades	70
		Reusable Seward Blades	71
		Reusable Soper Blades	71
		Reusable Robertshaw Blades	71
		Reusable Miller/Port Blades	71
		Reusable Flagg Blades	71
		Reusable American Profile Macintosh Blades	71
		Reusable American Profile Miller Blades	71
		Reusable Oxford Blade	71
		Reusable Sunflex Macintosh Blade (with Flexible Tip)	72
		Reusable Kessell Blades	72
Oxygen Tubing		Laryngoscope Accessories	
Oxygen Tubing	16	Laryngoscope Handles - Standard, Fibreoptic and LED Styles	72
		Laryngoscope Bulbs - Standard, Fibreoptic and LED	73
		Laryngoscope Spare Parts	73
CO₂ Sample Lines		Otoscope	
CO ₂ Sample Lines	16	Otoscope	73
Suction Accessories		Magill Forceps	
Suction Tubing	17	Magill Forceps	73
Yankaur Suction Handpieces	17		
		Secur-A-Tag Security Ties	
		Secur-A-Tag Security Ties	74
Respiratory Humidification		Blood Warmers	
Humidifiers	18	Biegler BW685S	75
Disposable Humidifier Chambers	18	Thermal Grade Extension Sets	76
Disposable Heater Wire Breathing Circuits	19, 20	Biegler Protherm II	77
Disposable Non-Heated Breathing Circuits	21	Biegler Autopress	78
Heated Wire Adaptors	22	Biegler Emergency Transport Warmer	78
Temperature Probes	22		
Disposable Water Traps	22	Oxygen Therapy Hardware	
Respiratory Filters (FHME)	23	High/Low Flowmeters	79
High Performance Bacterial Filter	24	Low Flowmeters	79
Respiratory Filters (HME)	24	Perflow High Precision Flowmeter	79
		Gauge Type Flowmeter	79
		O-Series Regulators	80
		Piston Regulators	80
Anaesthetic Breathing Circuits			
Disposable 3 Limb Adult	25-32		
Disposable 3 Limb Paediatric	33		
Disposable Adult and Paediatric Circuit (with ADU Connector)	34, 35		
Disposable 2 Limb Adult Circuits	36-41		
Disposable 2 Limb Paediatric Circuits	42		
Coaxial Circuits	43		
Disposable Resuscitation Circuits	44, 45		
Reusable Anaesthetic / Resuscitation Circuits	46		
Disposable Mid-O-Gas/Analgesia Breathing Circuits	47, 48		
Rebreathing Bags			
Rebreathing Bags	49		
Plastic Tubing			
Disposable Corrugated Tubing	50		
Disposable Expandable Tubing	51		
Reusable Smooth-Bor Tubing	52		
Masks			
Reusable Silicone Masks	53		
Disposable PVC Masks	53		
Reusable Mask Hooks	53		
Disposable Air Cushion Masks	54, 55		
Disposable Air Cushion Masks (MRI Compatible)	56		
Disposable Twin-port Air Cushion Masks	57		
Headgear	57, 59		
Teeth Opener	57		
CPAP Masks			
Non-Vented Full-face Masks	58, 59		
Non-Vented Nasal Face Mask	58, 59		



RESPIRATORY CARE PRODUCTS

Contents

Suction Attachments and Regulators

Suction Attachments - with tap	81
Suction Attachments - without tap	81
Suction Regulators	81
Suction Filters	82
Suction Filter Housing	82
Suction Filter Kits	82
Twin-O-Vac Suction/Oxygen Device	83
Central Vacuum Systems	83
Disposable Twin-O-Vac Suction Jar	83
3-Way Adaptors	84
Entonox Analgesic System	84
Humidifier	84

Oxygen Therapy Hardware Accessories

Handwheel Kits	85
Tubing	85
Hose Assembly	85
Gauges	85
Cylinder Key, Wingscrew and Spanner	85
Flowmeter Nipple and Wingnut	85
Bodok Seals	85
Twin-O-Vac Filter Spare	85

Non-contact Thermometers

	86
--	----

Infrared Ear Thermometer

	86
--	----

Blood Pressure and Glucose Monitors

Compact Blood Pressure Monitor	87
Desktop Blood Pressure Monitor	87
Desktop Blood Pressure Monitor with Cuff Storage	87
Blood Pressure and Glucose Monitors	87
Glucose Monitor	87

SpO₂ Monitors

Handheld and Finger SpO ₂ Monitors	88
Disposable SpO ₂ Sensors	89
Reusable SpO ₂ Sensors	89
Wrist Style SpO ₂ /PR Oximeters	90

Infant Care Products

Neoshades	91
Neosmile	91
Hold-A-Line	91
NeoBar	92
NeoBridge	92
Little Suckers	92
Meconium Aspirator	93
Suction Caddy	93
NeoBond	93
EZ-Hold	93
NeoHold	93
Duotube	94
Lead-Saver	94
NeoPulse	94
NeoLead	94
Micro NeoLead	94
NeoGrip	95
EZCare	95
NeoFlex	95
Nipple Shields	95
Comepa Neonatal Electrodes	95
BabyIVBoards	96
Neonatal SpO ₂ Probe Cover / Securing Strap	96
Infant Resuscitation Circuit	96
CPAP Kits	97

Electric Scooters and Mobility Equipment

4 Wheel Extra Large Scooter	98
4 Wheel Large Scooter	98
3 Wheel Large Scooter	98
4 Wheel Large Scooter with Closed Canopy	99
4 Wheel Large Electric Chair	99
Mini-size 4-Wheel Foldable Scooter	99
Electric Bicycles	99

Circuit Components - Polysulfone (PS)

(Autoclavable up to 134°C)



A2000Y(PS)
Autoclavable 'Y' Piece without Ports
22mm M/15mm F to 22mm M/22mm M
Reusable Long Life Polysulfone (PS)

GTIN: 9330783008780



A2000YT(PS) 'Y' with Temperature Port
Autoclavable 'Y' Piece with 8mm F Temperature Port
22mm M/15mm F to 22mm M/22mm M
Reusable Long Life Polysulfone (PS)
(Port Plug Available)

GTIN: 9330783008810



A2000YPT(PS) 'Y' with Temperature & Pressure Port GTIN: 9330783008803
Autoclavable 'Y' Piece with 8mm F Temperature & 6mm F Pressure Port
22mm M/15mm F to 22mm M/22mm M
Reusable Long Life Polysulfone (PS)
(Port Plug Available)



A2000YP(PS) 'Y' with Pressure Port GTIN: 9330783008797
Autoclavable 'Y' Piece with 6mm F Pressure Port
22mm M/15mm F to 22mm M/22mm M
Reusable Long Life Polysulfone (PS)
(Port Plug Available)



A2100(PS) GTIN: 9330783008858
Autoclavable Elbow Connector
22mm M/15mm F to 22mm F
Reusable Long Life Polysulfone (PS)



APM240150005 GTIN: 9330783009244
Autoclavable Straight Adaptor
Tube Connector 15mm M x 15mm ID Tube
Reusable Long Life Polysulfone (PS)



APM240640005 GTIN: 9330783009282
Autoclavable Hose Connector Tube Connector 12.5mm M x 15mm M
Reusable Long Life Polysulfone (PS)



A90001R(PS) GTIN: 9330783009091
Autoclavable Straight Adaptor 22mm M/15mm F to 22mm M/15mm F
Reusable Long Life Polysulfone (PS)



A2200(PS) GTIN: 9330783008896
Autoclavable Hose Adaptor 22mm M/15mm F to 15mm M/10mm F
Reusable Long Life Polysulfone (PS)




A2200T(PS) Adaptor with Temperature Port GTIN: 9330783008926
Autoclavable Adaptor with 8mm F Temperature Port
22mm M/15mm F to 15mm M/10mm F
Reusable Long Life Polysulfone (PS)
(Port Plug Available)

Legend: M = Male (External) F = Female (Internal) T = Temperature Port P = Pressure Port (ID) = Internal Diameter (PS) = Polysulfone

Circuit Components - Polysulfone (PS)


(Autoclavable up to 134°C)

22mm F  22mm M /15mm F


A2300(PS) GTIN: 9330783008971
Autoclavable Hose Adaptor 22mm F to 22mm M/15mm F
Reusable Long Life Polysulfone (PS)

8mm F Port
22mm F  22mm M /15mm F


A2300T(PS) Adaptor with Temperature Port GTIN: 9330783009008
Autoclavable Adaptor with 8mm F Temperature Port
22mm F to 22mm M/15mm F
Reusable Long Life Polysulfone (PS)
(Port Plug Available)

6mm ID Tube
22mm F  22mm M /15mm F


A2300TOF(PS) GTIN: 9330783009039
Autoclavable Hose Adaptor with O₂ Fitting
22mm F to 22mm M/15mm F
Reusable Long Life Polysulfone (PS)

Tube Size: 4mm
8mm
8.5mm  15mm M


A53504 Straight Adaptor 15mm M to 4.0mm ID Tube (PS) GTIN: 9330783009046
A53510 Straight Adaptor 15mm M to 8.0mm ID Tube (PS) GTIN: 9330783033447
A53511 Straight Adaptor 15mm M to 8.5mm ID Tube (PS) GTIN: 9330783033454
Reusable Long Life Polysulfone
NOTE: Other sizes available upon request

Suit 15mm ID Tube
with ridge  22mm F


APM240050002 GTIN: 9330783009213
Autoclavable Tube Connector 15mm ID Tube to 22mm F
Reusable Long Life Polysulfone (PS)

12.5mm M  15mm F


APM240740006(PS) GTIN: 9330783009343
Autoclavable Tube Connector 12.5mm M to 15mm F
Reusable Long Life Polysulfone (PS)

12.5mm M  22mm M /15mm F


APM240740003(PS) GTIN: 9330783009329
Autoclavable Tube Connector 12.5mm M to 22mm M/15mm F
Reusable Long Life Polysulfone (PS)

12.5mm M  22mm F


APM240840002 GTIN: 9330783009350
Autoclavable Tube Connector 12.5mm M to 22mm F
Reusable Long Life Polysulfone (PS)

Suit 15mm ID Tube
with ridge  15mm F Swivel


APM6060003(PS) GTIN: 9330783009534
Autoclavable Straight Swivel Connector 15mm ID Tube to 15mm F
Reusable Long Life Polysulfone (PS)

15mm M  15mm F Swivel

APM6120 GTIN: 9330783009565
Autoclavable 90° Swivel Elbow Connector 15mm M to 15mm F
Reusable Long Life Polysulfone (PS)

5mm ID Tube  6mm M

YPP500(PS) GTIN: 9330783024339
90° Pressure Port Adaptor 5mm tube to 6mm M
Suitable for use with A2000P/A2000TP adaptors
Reusable Long Life Polysulfone (PS)

6mm ID Tube  8mm M

NIPA2300T GTIN: 9330783018949
Oxygen Nipple 6mm tube to 8mm M, Reusable Long Life Polysulfone (PS)
Suitable for use with A2000T/A2000TP/A2200T/A2300T adaptors.

Legend: M = Male (External) F = Female (Internal) T = Temperature Port P = Pressure Port OF = Oxygen Fitting ID = Inside Diameter (PS) = Polysulfone

Circuit Components - Polycarbonate (PC)

(Autoclavable up to 110°C)



A2000(PC)
Autoclavable 'Y' Piece without Ports
22mm M/15mm F to 22mm M/22mm M
Reusable Clear Polycarbonate (PC)

GTIN: 9330783008681



A2000T(PC) 'Y' with Temperature Port
Autoclavable 'Y' Piece with 8mm F Temperature Port
22mm M/15mm F to 22mm M/22mm M
Reusable Clear Polycarbonate (PC)
(Port Plug Available)

GTIN: 9330783008735



A2000TP(PC) 'Y' with Temperature & Pressure Port GTIN: 9330783008759
Autoclavable 'Y' Piece with 8mm F Temperature & 6mm F Pressure Port
22mm M/15mm F to 22mm M/22mm M
Reusable Clear Polycarbonate (PC)
(Port Plug Available)



A2000P(PC) 'Y' with Pressure Port GTIN: 9330783008704
Autoclavable 'Y' Piece with 6mm F Pressure Port
22mm M/15mm F to 22mm M/22mm M
Reusable Clear Polycarbonate (PC)
(Port Plug Available)



A2100(PC) GTIN: 9330783008834
Autoclavable Elbow Connector
22mm M/15mm F to 22mm F
Reusable Clear Polycarbonate (PC)



APM240150005(PC) GTIN: 9330783009268
Autoclavable Straight Adaptor
Tube Connector 15mm M x 15mm ID Tube
Reusable Clear Polycarbonate (PC)



APM240640005(PC) GTIN: 9330783009305
Autoclavable Hose Connector Tube Connector 12.5mm M x 15mm M
Reusable Clear Polycarbonate (PC)



A90001(PC) GTIN: 9330783009077
Autoclavable Straight Adaptor 22mm M/15mm F to 22mm M/15mm F
Reusable Clear Polycarbonate (PC)



A2200(PC) GTIN: 9330783008889
Autoclavable Hose Adaptor 22mm M/15mm F to 15mm M/10mm F
Reusable Clear Polycarbonate (PC)




A2200T(PC) Adaptor with Temperature Port GTIN: 9330783008919
Autoclavable Adaptor with 8mm F Temperature Port
22mm M/15mm F to 15mm M/10mm F
Reusable Clear Polycarbonate (PC)
(Port Plug Available)


Legend: M = Male (External) F = Female (Internal) T = Temperature Port P = Pressure Port (PC) = Polycarbonate


Circuit Components - Polycarbonate (PC)


(Autoclavable up to 110°C)


22mm F  22mm M /15mm F **A2300(PC)** Autoclavable Hose Adaptor 22mm F to 22mm M/15mm F Reusable Clear Polycarbonate (PC) GTIN: 9330783008964


8mm F Port
22mm F  22mm M /15mm F **A2300T(PC) Adaptor with Temperature Port** Autoclavable Adaptor with 8mm F Temperature Port 22mm F to 22mm M/15mm F Reusable Clear Polycarbonate (PC) (Port Plug Available) GTIN: 9330783008995


6mm ID Tube
22mm F  22mm M /15mm F **A2300TOF(PC)** Autoclavable Hose Adaptor with O₂ Fitting 22mm F to 22mm M/15mm F Reusable Clear Polycarbonate (PC) GTIN: 9330783009022


Suit 15mm ID Tube with ridge
 22mm F **APM240050002(PC)** Autoclavable Tube Connector 15mm ID Tube to 22mm F Reusable Clear Polycarbonate (PC) GTIN: 9330783009237

12.5mm M  15mm F **APM240740006** Autoclavable Tube Connector 12.5mm M to 15mm F Reusable Clear Polycarbonate (PC) GTIN: 9330783009336

12.5mm M  22mm M /15mm F **APM240740003** Autoclavable Tube Connector 12.5mm M to 22mm M/15mm F Reusable Clear Polycarbonate (PC) GTIN: 9330783009312

12.5mm M  22mm F **APM240840002(PC)** Autoclavable Tube Connector 12.5mm M to 22mm F Reusable Clear Polycarbonate (PC) GTIN: 9330783009367

Suit 15mm ID Tube with ridge
 15mm F Swivel **APM6060003(PC)** Straight Swivel Connector 15mm ID Tube to 15mm F (PC) Reusable Clear Polycarbonate (PC) GTIN: 9330783009527

15mm M  15mm F Swivel **APM6120(PC)** 90° Swivel Elbow Connector 15mm M to 15mm F Reusable Clear Polycarbonate (PC) GTIN: 9330783009572

Legend: M = Male (External) F = Female (Internal) T = Temperature Port P = Pressure Port OF = Oxygen Fitting ID = Inside Diameter (PC) = Polycarbonate

RESPIRATORY CARE PRODUCTS



Circuit Components - Polycarbonate (PC)


(Autoclavable up to 110°C)



APM6125/M

Autoclavable 90° Elbow Connector 15mm M to 22mm M/15mm F
Reusable Clear Polycarbonate (PC)

GTIN: 9330783009589



YPP500(PC) Pressure Port Adaptor

90° Pressure Port Adaptor 5mm tube to 6mm M
Reusable Clear Polycarbonate (PC)
Suitable for use with A2000P/A2000TP adaptors

GTIN: 9330783024322

NIPA2300T(PC) Oxygen Fitting

Oxygen Nipple 6mm tube to 8mm M
Reusable Clear Polycarbonate (PC)
Suitable for use with A2000T/A2000TP/A2200T/A2300T adaptors.


GTIN: 9330783018963



APM6125/M/1/CO₂

Autoclavable Paediatric 90° Elbow Connector
with Port & Silicone Cap
15mm F/15mm M
Reusable Clear Polycarbonate (PC)

GTIN: 9330783009596



APM6103/22/PG

Autoclavable Paediatric 'Y' Connector with Port & Silicone Cap
22mm M/15mm F to 15mm M/15mm M (suit 15mm Smoothbore Tubing)
Reusable Clear Polycarbonate (PC)

GTIN: 9330783009541



APM6103/MF/22

Autoclavable Paediatric 'Y' Connector
22mm M/15mm F to 15mm M/15mm M (suit 15mm Smoothbore Tubing)
Reusable Clear Polycarbonate (PC)


GTIN: 9330783009558



DHD236

Autoclavable 'T' Piece Connector 22mm M/M/M
Reusable Clear Polycarbonate (PC)

GTIN: 9330783012619



DHD235/A

Autoclavable 'T' Piece Connector with One Way Valves 22mm M/M/M
Reusable Clear Polycarbonate (PC) (S)

GTIN: 9330783012602



7116 Catheter Mount Tubing (S)

Autoclavable Corrugated Silicone Smoothbore
Tubing 10mm ID - 15mm ID
Reusable and Latex Free

GTIN: 9330783007417

Legend: M = Male (External) F = Female (Internal) T = Temperature Port P = Pressure Port ID = Internal Diameter (PC) = Polycarbonate (S) = Silicone

Circuit Components - Single Patient Use (D)

22mm M
/15mm F



APM51100
Fixed 'Y' Piece
22mm M/22mm M/22mm M (PP)
Disposable

GTIN: 9330783036578

22mm M
/15mm F



APM51100L
Fixed 'Y' Piece with Luer Lock Port
22mm M/22mm M/22mm M (PP)
Disposable

GTIN: 9330783036578

22mm M
/15mm F



APM51100TP
Fixed 'Y' Piece with Temperature Port and Pressure Port (Including Port Plugs)
22mm M/22mm M/22mm M (PP)
Disposable

GTIN: 9330783036578

22mm M
/15mm F



MP1000
Split 'Y' Piece
22mm M/15mm F to 22mm ID Tube/22mm ID Tube (PP)
Disposable

GTIN: 9330783036943

22mm M
/15mm F



MP1000L
Split 'Y' Piece with Luer Lock Port and Plug
22mm M/15mm F to 22mm ID Tube/22mm ID Tube (PP)
Disposable

GTIN: 9330783036950

22mm M
/15mm F



MP1000TP
Split 'Y' Piece with Temperature Port and Pressure Port (Including Port Plugs)
22mm M/15mm F to 22mm ID Tube/22mm ID Tube (PP)
Disposable

GTIN: 9330783036967

22mm M

22mm M



APM51093
Swivel 'Y' Connector - No Port
22mm M X 22mm M to 22mm M / 15mm F (HDPE)
Disposable

GTIN: 9330783009442

22mm M

22mm M



APM51098
Swivel 'Y' Connector with CO₂ Sampling Port
22mm M X 22mm M to 22mm M / 15mm F (HDPE)
with Sampling Port
Disposable

GTIN: 9330783009459

Legend: M = Male (External) F = Female (Internal) T = Temperature Port (S) = Silicone (HDPE) = High Density Polyethylene (PP) = Polypropylene (D) = Disposable (PE) = Polyethylene (SY) = Styrene

Circuit Components - Single Patient Use (D)



5112
Paediatric 'Y' Connector (PP)
22mm M/15mm F to 15mm M/15mm M (PP)
Disposable

GTIN: 9330783006168



PH30015
Oxygen 'T' Piece
22mm M/22mm M/15mm F (SY)
Disposable

GTIN: 9330783028474



APM51464
One Way Valve 22mm F to 22mm M (PP)(S)
Disposable

GTIN: 9330783009497



APM51465
One Way Valve 22mm F to 22mm M (PP)(S)
Disposable

GTIN: 9330783009503



TB2000-1
Adaptor 10mm Tube to 15mm M/10mm F (SY)
Clear Styrene

GTIN: 9330783036929



IS146
Bleed Valve Assembly 15mm M Adjustable Bleed (PE)
Disposable

GTIN: 9330783017102

6mm ID Tube

1. 22mm M / 15mm F
2. 22mm F



1. 15mm M / 10mm F
2. 22mm M / 15mm F

Hose Adaptor with O₂ Fitting

1. A2200TOF(D) 22mm M/15mm F to 15mm M/10mm F (PP) Disposable GTIN: 9330783008933

2. A2300TOF(D) 22mm F to 22mm M / 15mm F (PP) Disposable GTIN: 9330783009015



56000
Plastic Mouth Piece 15mm F / 18mm M (HDPE)
Disposable

GTIN: 9330783006885



56001
Plastic Mouth Piece 22mm F (HDPE)
Disposable

GTIN: 9330783006892

Legend: M = Male (External) F = Female (Internal) T = Temperature Port (S) = Silicone (HDPE) = High Density Polyethylene (PP) = Polypropylene (D) = Disposable (PE) = Polyethylene (SY) = Styrene

Circuit Components - Single Patient Use (D)

22mm Tube		22mm F	PH3000(D) Hose Adaptor 22mm ID Tube to 22mm F (PP) Disposable	GTIN: 9330783019762
22mm M /15mm F		22mm M /15mm F	PH3100(D) Straight Connector 22mm M/15mm F to 22mm M/15mm F (PP) Disposable	GTIN: 9330783019779
22mm M /15mm F		19mm M /15mm F	PH3200(D) Straight Scavenge Connector 22mm M/15mm F to 19mm M/15mm F (PP) Disposable	GTIN: 9330783019786
22mm F		6mm Tube	PH3300(D) Hose Adaptor 22mm F to 6mm Tube (PP) Disposable	GTIN: 9330783019793
22mm F		10mm Tube	APM51450 Hose Adaptor 22mm F to 10mm Tube (PP) Disposable	GTIN: 9330783009480
22mm F		22mm F	APM51487 Adaptor 22mm F to 22mm F (PP) Disposable	GTIN: 9330783009510
8mm Tube		8mm Tube 8mm Tube	APM88409 'Y' Hose Adaptor 5mm to 8mm ID Tube/Tube/Tube (PP) Disposable	GTIN: 9330783009626
22mm M		22mm F	PH70055 MDI Adaptor 22mm M/F with Port and Cap (PP) Disposable	GTIN: 9330783032686
22mm M		22mm F 22mm M	PH70068 MDI 'T' Adaptor 22mm M/F/F with Port and Cap (PP) Disposable	GTIN: 9330783033478

Legend: M = Male (External) F = Female (PP) = Polypropylene (S) = Silicone

Circuit Components - Single Patient Use (D)



APM51046
Mask Elbow 15mm M to 22mm M / 15mm F (PP)
Disposable

GTIN: 9330783009435



APM51033
Mask 90° Elbow with Luer Lock Port and Plug (PP)
15mm M to 22mm M / 15mm F
Disposable

GTIN: 9330783009411



5126
Norman Mask Fresh Gas 90° Elbow (PP)
15mm M to 22mm M / 15mm F with 6mm Outlet
Disposable

GTIN: 9330783006229



5114
Washington 'T' Piece (PP)
15mm M to 15mm F with 6mm Outlet
Disposable

GTIN: 9330783006175



PH221211
Multi-size Tubing Joiner
Accepts 8mm ID to 15mm ID Tubing (PP)
Disposable

GTIN: 9330783036585



60011
Tubing Joiner 6mm tube to 6mm tube (PP)
Disposable

GTIN: 9330783031795



APM51035
Straight Adaptor with Luer Lock & Cap (PP)
15mm M to 22mm M / 15mm F (PP)
Disposable

GTIN: 9330783009428



ADU1016
ADU Connector (PP)

GTIN: 9330783036936



1968
Reducing Adaptor 6mm tube to 22mm M (PP)
Disposable

GTIN: 9330783004270



APM8821
Plastic Cap (Dust Cover)
Fits 22mm M Connections (PP)

GTIN: 9330783041985

Legend: M = Male (External) F = Female (Internal) (PP) = Polypropylene (D) = Disposable

Anaesthetic Components - Reusable



15mm M

22mm M
/15mm F

00-233 Reusable Mask Elbow
22mm M/15mm F X 15mm M

GTIN: 9330783000937

Material: Light Weight Metal

6mm ID Tube Connection



15mm M

15mm F

00-240 Reusable Washington T-Piece
15mm M X 15mm F with 6mm ID Tube Connection

GTIN: 9330783000968

Material: Light Weight Metal

6mm ID Tube Connection



15mm M

22mm M
/15mm F

00-234 Reusable Norman Mask Elbow
22mm M/15mm F X 15mm M with 6mm ID Tube Connection

GTIN: 9330783000944

Material: Light Weight Metal



22mm M

00-97 Reusable Bag Cage
22mm M

GTIN: 9330783001071

Material: Light Weight Metal



1/4" FPT

1/4" (6mm) ID Tube
Connection

00-489 Reusable Threaded Tubing Adaptor
1/4" FPT X 1/4" (6mm) ID Tubing Attachment
(For use with APM00199 Manometer Gauge)

GTIN: 9330783001064

Material: Metal



1/4" FPT

1/4" (6mm) ID Tube
Connection

1/2"-20 Thread

P-159 Reusable 'T' Adaptor
1/4" FPT to 1/2"-20 X 1/4" (6mm) ID Tubing
(For use with APM00199 Manometer Gauge)

GTIN: 9330783019533

Material: Metal



1/2"-20 Thread

00-273 Reusable Pressure Relief Valve
Adjustable 5-55 cmH2O. Mounting Thread - 1/2"-20
(To suit APMP-159 Adaptor)

GTIN: 9330783000982

Material: Metal

Legend: M = Male (External) F = Female (Internal) (FPT) = Female Pipe Thread (NPT)

Anaesthetic Components - Reusable



00-379 In-line Scavenge Valve
with 40cmH₂O Relief Valve

Material: Metal

GTIN: 9330783001033



00-389 Bag Tail Scavenge Valve
with 40cmH₂O Relief Valve

Material: Metal

GTIN: 9330783001057



00-199 Manometer Gauge

GTIN: 9330783000913

00-150 Gauge Guard (not shown)

GTIN: 9330783000906



00-274 Manometer Gauge Assembly

GTIN: 9330783000999

Consisting of:
00-199 Manometer Gauge
P-159 'T' Adaptor
00-273 Pressure Relief Valve



DF655 Reusable Expiratory Valve with Gas Outlet
22mm M to 22mm F SmoothBor Tubing Adaptor
with 6mm ID Tube Outlet

GTIN: 9310574006249

Material: Metal

Seal/Service Kit

DF655-SK01 Seal/Service Kit to suit DF655

GTIN: 9330783037506



4418
Reusable T-Piece for Oxygen Sensor Cell
22mm F x 22mm M x 22mm M/15mm F

GTIN: 9330783006069

Legend: M = Male (External) F = Female (Internal) (FPT) = Female Pipe Thread (NPT)


Reusable Catheter Mounts - LATEX FREE

(Autoclavable up to 134°C*)

	A. 10mm ID	B. 15mm ID	C. 15mm F Swivel	6000-1 GTIN: 9330783007042 Comprising of: A. 7116 Corrugated Silicone Tube 10-15mm (S) B. APM6120 Right Angle Swivel Elbow Connector 15mm M to 15mm F (PS)
15mm M	A. 10mm ID	B. 15mm ID	C. 15mm F Swivel	6000-3 GTIN: 9330783007226 Comprising of: A. APM240640005 Straight Connector 15mm M to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM6060002(PS) Straight Swivel Connector 15mm M to 15mm F (PS)
22mm F	A. 10mm ID	B. 15mm ID	C. 15mm F	6000-5 GTIN: 9330783007301 Comprising of: A. APM240840002 Straight Connector 22mm F to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM6125/M* Right Angle Elbow Connector 15mm M to 15mm F (PC) <small>* Polycarbonate only autoclavable to 110°C</small>
15mm F	A. 10mm ID	B. 15mm ID	C. 22mm F	6000-7 GTIN: 9330783007318 Comprising of: A. APM240740006(PS) Straight Connector 15mm F to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM240050002 Straight Connector 15mm M to 22mm F (PS)
22mm F	A. 10mm ID	B. 15mm ID	C. 15mm F Swivel	6000-8 GTIN: 9330783000319 Comprising of: A. APM240840002 Straight Connector 22mm F to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM6060002(PS) Straight Swivel Connector 15mm M to 15mm F (PS)
22mm M / 15mm F	A. 10mm ID	B. 15mm ID	C. 22mm F	6000-9 GTIN: 9330783007325 Comprising of: A. APM240740006(PS) Straight Connector 22mm M / 15mm F to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM240050002 Straight Connector 15mm M to 22mm F (PS)
15mm M	A. 10mm ID	B. 15mm ID	C. 22mm M / 15mm F	6000-12 GTIN: 9330783007073 Comprising of: A. APM240640005 Straight Connector 15mm M to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. A2200(PS) Straight Connector 15mm M to 22mm M / 15mm F (PS)
15mm M	A. 10mm ID	B. 15mm ID	C. 22mm F	6000-14 GTIN: 9330783007080 Comprising of: A. APM240640005 Straight Connector 15mm M to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM240050002 Straight Connector 15mm M to 22mm F (PS)
15mm M	A. 10mm ID	B. 15mm ID	C. 15mm F Swivel	6000-15 GTIN: 9330783007097 Comprising of: A. APM240640005 Straight Connector 15mm M to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM6120 Right Angle Swivel Elbow Connector 15mm M to 15mm F (PS)
22mm F	A. 10mm ID	B. 15mm ID	C. 15mm F Swivel	6000-16 GTIN: 9330783007110 Comprising of: A. APM240840002 Straight Connector 22mm F to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM6120 Right Angle Swivel Elbow Connector 15mm M to 15mm F (PS)
15mm F	A. 10mm ID	B. 15mm ID	C. 15mm F	6000-17 GTIN: 9330783007127 Comprising of: A. APM240740006(PS) Straight Connector 15mm F to 12.5mm M (PS) B. 7116 Corrugated Silicone Tube 10-15mm (S) C. APM6125/M* Right Angle Connector 15mm M to 15mm F (PC) <small>* Polycarbonate only autoclavable to 110°C</small>

Legend: M = Male (External) F = Female (Internal) ID = Inside Diameter (PC) = Polycarbonate (PS) = Polysulfone (S) = Silicone

Disposable Catheter Mounts - LATEX FREE




6000-25(D)

Comprising of:

- A. APM6060003(PC) Straight Swivel 15mm F (PC)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240050002(PC) Straight Connector 22mm F (PC)

GTIN: 9330783000272




6000-26(D)

Comprising of:

- A. APM6120(PC) Swivel 90° Elbow 15mm F (PC)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240050002(PC) Straight Connector 22mm F (PC)

GTIN: 9330783000289




6000-28(D)

Comprising of:

- A. APM6120(PC) Swivel 90° Elbow 15mm F (PC)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240150005(D) Straight Connector 15mm M (D)

GTIN: 9330783000296




6000-29(D)

Comprising of:

- A. APM6060003(PC) Straight Swivel 15mm F (PC)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240150005(D) Straight Connector 15mm M (D)

GTIN: 9330783000302




6000-30(D)

Comprising of:

- A. A2200(PC) Fixed Straight Adaptor 22mmM/15mm F (D)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240050002(PC) Straight Connector 22mm F (PC)

GTIN: 9330783007233




6000-33(D)

Comprising of:

- A. APM51046 90° Fixed Elbow 22mmM/15mm F (D)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240050002(PC) Straight Connector 22mm F (PC)

GTIN: 9330783007264




6000-34(D)

Comprising of:

- A. APM51046 90° Fixed Elbow 22mmM/15mm F (D)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240150005(D) Straight Connector 15mm M (D)

GTIN: 9330783007264




95003203 with Suction Port and Cap

Comprising of:

- A. I143 Double Swivel 22mm M/15mm F (PP)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240150005(D) Straight Connector 15mm M (D)

GTIN: 9330783008407



95003208 with Suction Port and Cap

Comprising of:

- A. I143 Double Swivel 22mm M/15mm F (PP)
- B. MP-EX15N6 Collapsible Tubing 15mm ID/15mm ID (PP)
- C. APM240050002(PC) Straight Connector 22mm F (PC)

GTIN: 9330783008421

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (PC) = Polycarbonate

"Softi" Smoothflow® Cannulas & Tubing



"Softi" Smoothflow® Oxygen Nasal Curved Tips Cannula (PVC)

- Contoured "Softi" Smoothflow® cannula provides total patient comfort
- Kink resistant tubing with inner star lumen design
- Over the ear lit secures position on the patient

Part No.	Cannula Size	Tube Length	
PH1161	Adult	2.1 metre	GTIN: 9330783019670
PH1163	Paediatric	2.1 metre	GTIN: 9330783019700
PH1164	Infant	2.1 metre	GTIN: 9330783019724
PH1165	Neonate	2.1 metre	GTIN: 9330783019748

"Softi" Smoothflow® Cannula



"Softi" Smoothflow® Oxygen Nasal CO₂ Sampling Cannula (PVC)(PC)(PP)

- Contoured "Softi" Smoothflow® cannula provides total patient comfort
- The design allows delivery of oxygen to one nasal prong and sampling of expired CO₂ gases from the other prong simultaneously
- Kink resistant tubing with inner star lumen design
- Over the ear lit secures position on the patient

Part No.	Cannula Size	Tube Length	Connection	
PHC02-261	Adult	2.1 metre	Male Luer	GTIN: 9330783020171
PHC02-261/F	Adult	2.1 metre	Fem. Luer	GTIN: 9330783020188
PHC02-261-4	Adult	4.0 metre	Male Luer	GTIN: 9330783020195
PHC02-261-4/F	Adult	4.0 metre	Fem. Luer	GTIN: 9330783020201
PHC02-263	Paediatric	2.1 metre	Male Luer	GTIN: 9330783033485
PHC02-263/F	Paediatric	2.1 metre	Fem. Luer	GTIN: 9330783033492
PHC02-263-4	Paediatric	4.0 metre	Male Luer	GTIN: 9330783033508
PHC02-263-4/F	Paediatric	4.0 metre	Fem. Luer	GTIN: 9330783033515

"Softi" Smoothflow® Cannula



Legend: M = Male (External) F = Female (Internal) ID = Internal Diameter (PVC) = Polyvinyl Chloride (PC) = Polycarbonate

“Softi” Smoothflow® Cannulas & Tubing

“Softi” Smoothflow® Oxygen Nasal Curved Tips Cannula (PVC)



DR-0165

- Contoured “Softi” Smoothflow® cannula provides total patient comfort
- Over the ear lit secures position on the patient

Part No.	Cannula Size	Connection	
DR-0161	Adult	6mm M Barbed Nipple	GTIN: 9330783034215
DR-0163	Paediatric	6mm M Barbed Nipple	GTIN: 9330783034222
DR-0164	Infant	6mm M Barbed Nipple	GTIN: 9330783034239
DR-0165	Neonate	6mm M Barbed Nipple	GTIN: 9330783034246

“Softi” Smoothflow® Oxygen Nasal Curved Tips Cannula (PVC)



NC-05550

- Contoured “Softi” Smoothflow® cannula provides total patient comfort
- Over the ear lit secures position on the patient

Part No.	Cannula Size	Connection	
NC-01110	Adult	Female	GTIN: 9330783040452
NC-03330	Paediatric	Female	GTIN: 9330783040469
NC-04440	Infant	Female	GTIN: 9330783040476
NC-05550	Neonate	Female	GTIN: 9330783040483

“Softi” Smoothflow® Diagnostic Cannula (PVC)



- Can be used in conjunction with pressure transducers or thermistors in Polysomnography (Sleep Studies)
- Optimizes pressure/snore waveform signals
- Inexpensive
- Contoured “Softi” Smoothflow® cannula provides total patient comfort
- Over the ear lit secures position on the patient
- Available in various lengths
- Male luer lock connection

Part No.	Size	Length	Connection	
LD50ML61	Adult	50cm	Male Luer	GTIN: 9330783040490
LD50ML63	Paediatric	50cm	Male Luer	GTIN: 9330783040506
LD100ML61	Adult	100cm	Male Luer	GTIN: 9330783040513
LD100ML63	Paediatric	100cm	Male Luer	GTIN: 9330783040520
LD130ML61	Adult	130cm	Male Luer	GTIN: 9330783040537
LD130ML63	Paediatric	130cm	Male Luer	GTIN: 9330783040544

Extension Line with Stopcock



Part No.	Length	I.D.	O.D.
DC-10HP	100mm	1.0mm	2.0mm
DC-15HP	150mm	1.0mm	2.0mm
DC-10-15HP	100mm	1.5mm	3.0mm
DC-15-15HP	150mm	1.5mm	3.0mm

Pressure Monitoring Line



Part No.	Length	I.D.	O.D.
DC-0150HP	150mm	1.0mm	2.0mm
DC-0200HP	200mm	1.0mm	2.0mm

Legend: M = Male (External) F = Female (Internal) T = Temperature Port P = Pressure Port ID = Internal Diameter (PC) = Polycarbonate (S) = Silicone

Oxygen Tubing and CO₂ Sample Lines



Non-Kink Star Lumen Disposable Oxygen Tubing (PVC)

- Kink resistant tubing with inner star lumen design
- Flexible PVC tubing with standard connectors
- Lengths available in 2.1, 3.0, 6.0 and 10.0 metres

34502.1/1	2.1m, 6mmM, 4mmF Oxygen Tube	GTIN: 9330783005444
34502/3	3.0m, 6mmM, 4mmF Oxygen Tube	GTIN: 9330783005451
34502/6	6.0m, 6mmM, 4mmF Oxygen Tube	GTIN: 9330783031771
34502/10	10.0m, 6mmM, 4mmF Oxygen Tube	GTIN: 9330783032440

Other Lengths Available By Request



CO₂ Sample Lines with Inner "Star Lumen" Design (PVC)(PC)

Male to Female Luer Connection

MP0300	3.0 metre (10ft.) CO ₂ Sample Line	GTIN: 9330783024841
	M/F Luer connection	
MP0315	4.5 metre (15ft.) CO ₂ Sample Line	GTIN: 9330783030910
	M/F Luer connection	
MP0330	9.0 metre (30ft.) CO ₂ Sample Line	GTIN: 9330783030903
	M/F Luer connection	



Male to Male Luer Connection

MP0500	3.0 metre (10ft.) CO ₂ Sample Line	GTIN: 9330783025060
	M/M Luer connection	
MP0515	4.5 metre (15ft.) CO ₂ Sample Line	GTIN: 9330783030927
	M/M Luer connection	
MP0530	9.0 metre (30ft.) CO ₂ Sample Line	GTIN: 9330783030897
	M/M Luer connection	

Legend: M = Male (External) F = Female (Internal) T = Temperature Port P = Pressure Port ID = Internal Diameter (PC) = Polycarbonate (S) = Silicone

Suction Tubing and Accessories

Suction Connection Tubing (PVC)

- Minimal coil memory
- Clear PVC formulation with striated walls
- Flexible tubing with Yellow cuffed ends
- 6mm lumen



Part No. Length

Single Packed - Non-Sterilised

DC6015	1.5m	GTIN: 9330783026333
DC6020	2.0m	GTIN: 9330783026357
DC6030	3.0m	GTIN: 9330783026371
DC6045	4.5m	GTIN: 9330783026395
DC6060	6.0m	GTIN: 9330783026418

Double Packed - Sterilised

DC6015(DPS)	1.5m	GTIN: 9330783026630
DC6020(DPS)	2.0m	GTIN: 9330783031757
DC6030(DPS)	3.0m	GTIN: 9330783031733
DC6045(DPS)	4.5m	GTIN: 9330783026661
DC6060(DPS)	6.0m	GTIN: 9330783026678



Yankaur Suction Handpiece (PVC) - Sterile Packed

PH/YSH/001 Suction Handpiece Crown tip without Vent GTIN: 9330783032471

PH/YSH/002 Suction Handpiece Crown tip with Vent GTIN: 9330783032457



Yankaur Suction Handpiece (PVC) - Sterile Packed

PH/YSH/003 Suction Handpiece Standard Tip without Vent GTIN: 9330783032488

PH/YSH/004 Suction Handpiece Standard tip with Vent GTIN: 9330783032464



Suction Tubing with Yankaur Suction Handpiece (PVC) - Sterile Packed

DC6015/YS(DPS) 1.5m Suction Tubing and Crown Tip Yankaur Suction Handpiece without Vent - Double Packed Sterile GTIN: 9330783032495

DC6020/YS(DPS) 2.0m Suction Tubing and Crown Tip Yankaur Suction Handpiece without Vent - Double Packed Sterile GTIN: 9330783031764

DC6030/YS(DPS) 3.0m Suction Tubing and Crown Tip Yankaur Suction Handpiece without Vent - Double Packed Sterile GTIN: 9330783031740

DC6045/YS(DPS) 4.5m Suction Tubing and Crown Tip Yankaur Suction Handpiece without Vent - Double Packed Sterile GTIN: 9330783031726

DC6060/YS(DPS) 6.0m Suction Tubing and Crown Tip Yankaur Suction Handpiece without Vent - Double Packed Sterile GTIN: 9330783031702



Suction Tubing (Yellow) - 10m Roll (PVC)

Designed for use behind bed areas or where there is a risk of beds, furniture or foot traffic cutting off suction.

- For all areas with High Suction outlets
- Will not collapse under high suction

MP-YS0010 Yellow Suction Tubing (10m. Roll) 8mm ID/13.5mmOD GTIN: 9330783031573

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (PVC) = Polyvinyl Carbonate (DPS) = Double Packed Sterile

Respiratory Humidification



Respiratory Humidifiers

Respiratory humidification is a method of the artificial warming and humidification of respiratory gas for the patient during mechanical ventilation.

When a patient requires mechanical ventilation over a longer period of time, it is imperative to take measures to compensate for the losses of heat and moisture in order to avoid complications.

[V-2600 Assembled Unit Pictured]



VH2100



VH2600

Respiratory Humidifiers

VH2100
Respiratory Humidifier
[General Purpose]

GTIN: 9330783040209

VH2600
Respiratory Humidifier
Microprocessor Servo Control and Heater Wire Circuit

GTIN: 9330783034512



22mm M

22mm M

VH3140 Disposable Humidifier Chamber

GTIN: 9330783034826

Inlet / Outlet size : 22M / 22M
Compressible volume-empty : 325 ml
Compressible volume-full : 175 ml
Maximum water capacity : 150 ml
Heat conduction diameter : 12.2 cm
Maximum peak flow : 180 L/min (in 40sec)
Resistance to flow : 0.15 cmH₂O @30 LPM
Compliance : 0.38 ml / cmH₂O
Chassis diameter : 12.2 cm
Material : Single use material



22mm M

22mm M

VH3141 Disposable Autofeed Humidifier Chamber

GTIN: 9330783034833


Inlet / Outlet size : 22M / 22M
Compressible volume-empty : 325 ml
Compressible volume-full : 175 ml
Maximum water capacity : 150 ml
Heat conduction diameter : 12.2 cm
Maximum peak flow : 180 L/min (in 40sec)
Resistance to flow : 0.15 cmH₂O @30 LPM
Compliance : 0.38 ml / cmH₂O
Chassis diameter : 12.2 cm
IV line length : 130 cm
Material : Single use material

RESPIRATORY CARE PRODUCTS



Disposable Corrugated Heater Wire Breathing Circuit Adult

Heated



CIRCUIT		Type	A	B	C	D	E	F	G	H	Equivalent
VA-4213-21	Adult	F&P	150cm, 1 pc	75cm, 2pcs	60cm, 1pc	•	22mm M, 1pc	Single	•		
VA-4213-22	Adult	F&P	150cm, 1 pc	150cm, 1 pc	60cm, 1pc	•		Dual	•		
VA-4213-23	Adult	F&P	180cm, 1 pc	90cm, 2pcs	60cm, 1pc	•	22mm M, 1pc	Single	•		
VA-4213-24	Adult	F&P	150cm, 1 pc	150cm, 1 pc	60cm, 1pc	•		Dual	•	•	RT200
VA-4213-25	Adult	F&P	150cm, 1 pc	75cm, 2pcs	60cm, 1pc	•	22mm M, 1pc	Single	•	•	RT205
VA-4213-26	Adult	F&P	150cm, 1 pc	75cm, 2 pcs	60cm, 1pc	•	22mm M, 1pc	Single	•	•	RT206
VA-4213-11N			Optional Nebulizer Adaptor Kit								

GTIN: 9330783040292

GTIN: 9330783034734

GTIN: 9330783034727


GTIN: 9330783040308

GTIN: 9330783040315

GTIN: 9330783040322

GTIN: 9330783040377


Heated



CIRCUIT		Type	A	B	C	D	E	F
VA-4213-29	Adult	F&P	22mmOD/22mmID, 1pc with Temp Port	150cm, 1pc	1pc	1pc	1pc	Single

GTIN: 9330783042043


Heated



CIRCUIT		Type	A	B	C	D	E	F
VA-4213-27	Adult	F&P	22mmOD/22mmID, 1pc with Temp Port	150cm, 1pc	1pc	Single	60cm, 1pc	22mm M, 2pc

GTIN: 9330783040254


Disposable Corrugated Heater Wire Breathing Circuit Paediatric



CIRCUIT		Type	A	B	C	D	E	F	G	H	Equivalent
VA-4223-21	Paediatric	F&P	130cm, 1 pc	60cm, 2pcs	60cm, 1pc	•	22mm M, 1pc	Single	•		
VA-4223-22	Paediatric	F&P	130cm, 1 pc	130cm, 1 pc	60cm, 1pc	•		Dual	•		
VA-4223-23	Paediatric	F&P	130cm, 1 pc	60cm, 2pcs	60cm, 1pc	•	22mm M, 1pc	Single	•	•	
VA-4223-24	Paediatric	F&P	130cm, 1 pc	130cm, 1 pc	60cm, 1pc	•		Dual	•	•	
VA-4213-11N	Optional Nebulizer Adaptor Kit										

GTIN: 9330783034710
 GTIN: 9330783034703
 GTIN: 9330783040353
 GTIN: 9330783040360
 GTIN: 9330783040377

Neonatal



CIRCUIT		Type	A	B	B1	C	D	E	F	G	H	Equivalent
VA-4233-21	Neonate	F&P	130cm, 1 pc	60cm, 2pc		60cm, 1pc	•	22mm M, 1pc	Single	•		
VA-4233-22	Neonate	F&P	130cm, 1 pc	130cm, 1 pc		60cm, 1pc	•		Dual	•		
VA-4233-23	Neonate	F&P	110cm, 1 pc	60cm, 2pc	30cm, 1pc	60cm, 1pc	•	22mm M, 1pc	Single	•		
VA-4233-24	Neonate	F&P	130cm, 1 pc	110cm, 1 pc	30cm, 1pc	60cm, 1pc	•		Dual	•		
VA-4233-25	Neonate	F&P	130cm, 1 pc	60cm, 2pc		60cm, 1pc	•	22mm M, 1pc	Single	•	•	
VA-4233-26	Neonate	F&P	130cm, 1 pc	130cm, 1 pc		60cm, 1pc	•		Dual	•	•	
VA-4233-27	Neonate	F&P	110cm, 1 pc	60cm, 2pc	30cm, 1pc	60cm, 1pc	•	22mm M, 1pc	Single	•	•	RT225/RT226
VA-4233-28	Neonate	F&P	130cm, 1 pc	110cm, 1 pc	30cm, 1pc	60cm, 1pc	•		Dual	•	•	

GTIN: 9330783034697
 GTIN: 9330783034680
 GTIN: 9330783040216
 GTIN: 9330783040223
 GTIN: 9330783040230
 GTIN: 9330783040247
 GTIN: 9330783040254
 GTIN: 9330783040261


RESPIRATORY CARE PRODUCTS



Disposable Corrugated Non-heated Breathing Circuits

Adult


Non-Heated



CIRCUIT		A	B	C	Cuff	D	E	F	G	H
VA-4212-01	Adult	75cm, 2pcs	75cm, 2pcs	60cm, 1pc	22mm F	•	22mm M, 2pcs	•	•	•
VA-4212-02	Adult	75cm, 2pcs	75cm, 2pcs	60cm, 1 pc	22mm F	•	22mm M, 2pcs	•	•	•
VA-4222-02	Paediatric	75cm, 2pcs	75cm, 2pcs	60cm, 1pc	15mm F	•	22mm M, 2pcs	•	•	•
VA-4212-01N	Optional Nebulizer Adaptor Kit									

GTIN: 9330783040834
 GTIN: 9330783040841
 GTIN: 9330783040902
 GTIN: 9330783042487


Non-Heated



CIRCUIT		Type	A	B	C	D	E
VA-4432-04	Adult	F&P	22mmOD/22mmID, 1pc with Temp Port	150cm, 1pc	60cm, 1pc	22mmOD / 22mmID, 3pc	Single

GTIN: 9330783042050

Non-Heated



CIRCUIT		Type	A	B	C	D	E
VA-4432-06	Adult	F&P	22mmOD/22mmID, 2pc	150cm, 1pc	1pc	1pc	Single

GTIN: 9330783042067

Respiratory Humidification



HWA-730F01



HWA-730F02



HWA-850F02



HWA-850F01

Heated Wire Adaptors

HWA-730F01 Disposable Single Heated Wire F & P Connector (Compatible with MR901) GTIN: 9330783041602

HWA-730F02 Disposable Dual Heated Wire F & P Connector (Compatible with MR900) GTIN: 9330783041619

HWA-850F01 Disposable Single Heated Wire F & P Connector (Compatible with MR801) GTIN: 9330783041626

HWA-850F02 Disposable Dual Heated Wire F & P Connector (Compatible with MR806) GTIN: 9330783041633



Temperature Probes

HWP73011 Reusable Temperature Probe (110cm) (Compatible with MR568) GTIN: 9330783040179

HWP73013 Reusable Temperature Probe (130cm) (Compatible with MR560) GTIN: 9330783034673

HWP73015 Reusable Temperature Probe (150cm) (Compatible with MR569) GTIN: 9330783034666

HWP73018 Reusable Temperature Probe (180cm) (Compatible with MR561) GTIN: 9330783034659



PH70520



PH70521

Disposable Water Traps

PH70520 Adult Disposable Water Trap 22mm M (PP) GTIN: 9330783029020

PH70521 Child Disposable Water Trap 15mm M (PP) GTIN: 9330783029037



FHME Filters (Microwell)

VH3110 FHME Filter - Adult	GTIN: 9330783034642
VH3210 FHME Filter - Adult/Child	GTIN: 9330783034628
VH3310 FHME Filter - Paediatric	GTIN: 9330783034604



FHME Angled Filters (Microwell)

VH3115 FHME Angled Filter - Adult	GTIN: 9330783040148
VH3215 FHME Angled Filter - Adult/Child	GTIN: 9330783040155
VH3315 FHME Angled Filter - Paediatric	GTIN: 9330783040162



FHME Filters (Pleated)

VH3100 FHME Filter - Adult	GTIN: 9330783034635
VH3200 FHME Filter - Adult/Child	GTIN: 9330783034611
VH3300 FHME Filter - Paediatric	GTIN: 9330783034598



FHME Angled Filters (Pleated)

VH3105 FHME Angled Filter - Adult	GTIN: 9330783040117
VH3205 FHME Angled Filter - Adult/Child	GTIN: 9330783040124
VH3305 FHME Angled Filter - Paediatric	GTIN: 9330783040131

Respiratory Filtration



HME-P Filters

VH1130	HME-P Filter - Adult Filter (Pleated)	GTIN: 9330783034574
VH1131	HME-P Filter - Adult Filter (Microwell)	GTIN: 9330783034581
VH1132	HME-P Filter - Paediatric (Microwell)	GTIN: 9330783034567
VH1133	HME-P Filter - Neonate (Microwell)	GTIN: 9330783034550



HME-T Filters

VH1160	HME-T Filter with Oxygen Port (Microwell)	GTIN: 9330783040186
VH1161	HME-T Filter (Microwell)	GTIN: 9330783034543
VH1161-1	Oxygen Kit for VH1161	GTIN: 9330783040193



Disposable High Performance Bacterial Filters


VF2150	Adult Bacterial Filter with Port	GTIN: 9330783034529
VF2151	Paediatric Bacterial Filter with Port	GTIN: 9330783040988
VF2152	Neonate Bacterial Filter with Port	GTIN: 9330783040995



Disposable High Performance / Efficiency Bacterial Filters

VF2160	Adult Bacterial Filter without Port	GTIN: 9330783040957
VF2160-1	Paediatric Bacterial Filter with Port	GTIN: 9330783040964
VF2170	Adult / Paediatric High Efficiency Bacterial Filter without Port	GTIN: 9330783034536

Corrugated Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' 22mm M/M/M Adaptor (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-100 Circuit - Length = 1.5m (60")

GTIN: 9330783038039

A3C72-100 Circuit - Length = 1.8m (72")


GTIN: 9330783038046

A3C96-100 Circuit - Length = 2.4m (96")

GTIN: 9330783038053

A3C120-100 Circuit - Length = 3.0m (120")

GTIN: 9330783038060



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51100L Fixed 'Y' 22mm M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-200 Circuit - Length = 1.5m (60")

GTIN: 9330783038077

A3C72-200 Circuit - Length = 1.8m (72")


GTIN: 9330783038084

A3C96-200 Circuit - Length = 2.4m (96")

GTIN: 9330783038091

A3C120-200 Circuit - Length = 3.0m (120")

GTIN: 9330783038107



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-300 Circuit - Length = 1.5m (60")

GTIN: 9330783038114

A3C72-300 Circuit - Length = 1.8m (72")

GTIN: 9330783038121

A3C96-300 Circuit - Length = 2.4m (96")


GTIN: 9330783038138

A3C120-300 Circuit - Length = 3.0m (120")

GTIN: 9330783038145

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Corrugated Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-400 Circuit - Length = 1.5m (60")

GTIN: 9330783038152

A3C72-400 Circuit - Length = 1.8m (72")


GTIN: 9330783038169

A3C96-400 Circuit - Length = 2.4m (96")

GTIN: 9330783038176

A3C120-400 Circuit - Length = 3.0m (120")

GTIN: 9330783038183



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x MP1000L Split 'Y' M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-500 Circuit - Length = 1.5m (60")

GTIN: 9330783038657

A3C72-500 Circuit - Length = 1.8m (72")


GTIN: 9330783038664

A3C96-500 Circuit - Length = 2.4m (96")

GTIN: 9330783038671

A3C120-500 Circuit - Length = 3.0m (120")

GTIN: 9330783038688



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-600 Circuit - Length = 1.5m (60")

GTIN: 9330783038695

A3C72-600 Circuit - Length = 1.8m (72")

GTIN: 9330783038701

A3C96-600 Circuit - Length = 2.4m (96")


GTIN: 9330783038718

A3C120-600 Circuit - Length = 3.0m (120")

GTIN: 9330783038725

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Corrugated Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' 22mm M/M/M Adaptor (HDPE)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-700 Circuit - Length = 1.5m (60")

GTIN: 9330783038732

A3C72-700 Circuit - Length = 1.8m (72")


GTIN: 9330783038749

A3C96-700 Circuit - Length = 2.4m (96")

GTIN: 9330783038756

A3C120-700 Circuit - Length = 3.0m (120")

GTIN: 9330783038763



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x 2MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51098 Swivel 'Y' M/M/M Adaptor with Capped Luer Port (HDPE)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3C60-800 Circuit - Length = 1.5m (60")

GTIN: 9330783038770

A3C72-800 Circuit - Length = 1.8m (72")


GTIN: 9330783038787

A3C96-800 Circuit - Length = 2.4m (96")

GTIN: 9330783038794

A3C120-800 Circuit - Length = 3.0m (120")

GTIN: 9330783038800



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3C60-102 Circuit - Length = 1.5m (60")

GTIN: 9330783038190

A3C72-102 Circuit - Length = 1.8m (72")

GTIN: 9330783038206

A3C96-102 Circuit - Length = 2.4m (96")


GTIN: 9330783038213

A3C120-102 Circuit - Length = 3.0m (120")

GTIN: 9330783038220

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Corrugated Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3C60-302 Circuit - Length = 1.5m (60")

GTIN: 9330783038237

A3C72-302 Circuit - Length = 1.8m (72")


GTIN: 9330783038244

A3C96-302 Circuit - Length = 2.4m (96")

GTIN: 9330783038251

A3C120-302 Circuit - Length = 3.0m (120")

GTIN: 9330783038268



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3C60-402 Circuit - Length = 1.5m (60")

GTIN: 9330783038275

A3C72-402 Circuit - Length = 1.8m (72")


GTIN: 9330783038282

A3C96-402 Circuit - Length = 2.4m (96")

GTIN: 9330783038299

A3C120-402 Circuit - Length = 3.0m (120")

GTIN: 9330783038305



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3C60-602 Circuit - Length = 1.5m (60")

GTIN: 9330783038312

A3C72-602 Circuit - Length = 1.8m (72")


GTIN: 9330783038329

A3C96-602 Circuit - Length = 2.4m (96")

GTIN: 9330783038336

A3C120-602 Circuit - Length = 3.0m (120")

GTIN: 9330783038343



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' M/M/M Adaptor (HDPE)
- D. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3C60-702 Circuit - Length = 1.5m (60")

GTIN: 9330783038350

A3C72-702 Circuit - Length = 1.8m (72")

GTIN: 9330783038367

A3C96-702 Circuit - Length = 2.4m (96")


GTIN: 9330783038374

A3C120-702 Circuit - Length = 3.0m (120")

GTIN: 9330783038381

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Expandable Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' 22mm M/M/M Adaptor (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-100 Circuit - Length = 1.8m (72")

GTIN: 9330783038572

A3E120-100 Circuit - Length = 3.0m (120")

GTIN: 9330783038589



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51100L Fixed 'Y' 22mm M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-200 Circuit - Length = 1.8m (72")

GTIN: 9330783038596

A3E120-200 Circuit - Length = 3.0m (120")

GTIN: 9330783038602



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-300 Circuit - Length = 1.8m (72")


GTIN: 9330783038619

A3E120-300 Circuit - Length = 3.0m (120")

GTIN: 9330783038626

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Expandable Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-400 Circuit - Length = 1.8m (72")

GTIN: 9330783038633

A3E120-400 Circuit - Length = 3.0m (120")

GTIN: 9330783038640



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x MP1000L Split 'Y' M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-500 Circuit - Length = 1.8m (72")

GTIN: 9330783038497

A3E120-500 Circuit - Length = 3.0m (120")

GTIN: 9330783038503



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-600 Circuit - Length = 1.8m (72")


GTIN: 9330783038510

A3E120-600 Circuit - Length = 3.0m (120")

GTIN: 9330783038527

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Expandable Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' 22mm M/M/M Adaptor (HDPE)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-700 Circuit - Length = 1.8m (72")

GTIN: 9330783038534

A3E120-700 Circuit - Length = 3.0m (120")

GTIN: 9330783038541



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51098 Swivel 'Y' M/M/M Adaptor with Capped Luer Port (HDPE)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)

A3E72-800 Circuit - Length = 1.8m (72")

GTIN: 9330783038558

A3E120-800 Circuit - Length = 3.0m (120")

GTIN: 9330783038565



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3E72-102 Circuit - Length = 1.8m (72")


GTIN: 9330783038398

A3E120-102 Circuit - Length = 3.0m (120")

GTIN: 9330783038404

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Expandable Anaesthetic Breathing Circuits - Disposable



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3E72-302 Circuit - Length = 1.8m (72")

GTIN: 9330783038411

A3E120-302 Circuit - Length = 3.0m (120")

GTIN: 9330783038428



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3E72-402 Circuit - Length = 1.8m (72")

GTIN: 9330783038435

A3E120-402 Circuit - Length = 3.0m (120")

GTIN: 9330783038442



Adult Disposable Anaesthetic Circuit

Consisting of:


- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3E72-602 Circuit - Length = 1.8m (72")

GTIN: 9330783038459

A3E120-602 Circuit - Length = 3.0m (120")

GTIN: 9330783038466



Adult Disposable Anaesthetic Circuit

Consisting of:

- A. 3 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' M/M/M Adaptor (HDPE)
- D. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- E. 2 x PH3100(D) 22mm M/M Adaptor (PP)
- F. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- G. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A3E72-702 Circuit - Length = 1.8m (72")


GTIN: 9330783038473

A3E120-702 Circuit - Length = 3.0m (120")

GTIN: 9330783038480

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (HDPE) = High Density Polyethylene (***) = Tubing Length

Paediatric Anaesthetic Breathing Circuits - Disposable



Disposable Paediatric Anaesthetic Circuit - Corrugated

Consisting of:


- A. 3 x APM240050002(PC) 15mm ID tube to 22mm F (PC)
- B. 2 x MP-CO15N*** 15mm ID Natural Corrugated Disposable Tube (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F/22mm M (PP)
- D. 1 x A2200(PC) 15mm ID tube to 22mm M (PC)
- E. 1 x MP-CO15B*** 15mm ID Blue Corrugated Disposable Tube (PP)
- F. 1 x PH3100(D) 15mm F/22mm M to 15mm F/22mm M Adaptor
- G. 1 x 54936 1.0 Litre Latex Free Disposable Rebreathing Bag (CP)

P3C72-100 Circuit - Length = 1.8m (72")

GTIN: 9330783038817

P3C120-100 Circuit - Length = 3.0m (120")

GTIN: 9330783038824



Disposable Paediatric Anaesthetic Circuit - Corrugated

Consisting of:


- A. 3 x APM240050002(PC) 15mm ID tube to 22mm F (PC)
- B. 2 x MP-CO15N*** 15mm ID Natural Corrugated Disposable Tube (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F/22mm M (PP)
- D. 1 x A2200(PC) 15mm ID tube to 22mm M (PC)
- E. 1 x MP-CO15B*** 15mm ID Blue Corrugated Disposable Tube (PP)
- F. 1 x PH3100(D) 15mm F/22mm M to 15mm F/22mm M Adaptor
- G. 1 x 54936 1.0 Litre Latex Free Disposable Rebreathing Bag (CP)
- H. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

P3C72-102 Circuit - Length = 1.8m (72")

GTIN: 9330783038831

P3C120-102 Circuit - Length = 3.0m (120")

GTIN: 9330783038848



Disposable Paediatric Anaesthetic Circuit - Expandable

Consisting of:


- A. 3 x APM240050002(PC) 15mm ID tube to 22mm F (PC)
- B. 2 x MP-EX15N*** 15mm ID Natural Expandable Disposable Tube (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F/22mm M (PP)
- D. 1 x A2200(PC) 15mm ID tube to 22mm M (PC)
- E. 1 x MP-EX15B*** 15mm ID Blue Expandable Disposable Tube (PP)
- F. 1 x PH3100(D) 15mm F/22mm M to 15mm F/22mm M Adaptor
- G. 1 x 54936 1.0 Litre Latex Free Disposable Rebreathing Bag (CP)

P3E72-100 Circuit - Length = 1.8m (72")

GTIN: 9330783038855

P3E120-100 Circuit - Length = 3.0m (120")

GTIN: 9330783038862



Disposable Paediatric Anaesthetic Circuit - Expandable

Consisting of:

- A. 3 x APM240050002(PC) 15mm ID tube to 22mm F (PC)
- B. 2 x MP-EX15N*** 15mm ID Natural Expandable Disposable Tube (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F/22mm M (PP)
- D. 1 x A2200(PC) 15mm ID tube to 22mm M (PC)
- E. 1 x MP-EX15B*** 15mm ID Blue Expandable Disposable Tube (PP)
- F. 1 x PH3100(D) 15mm F/22mm M to 15mm F/22mm M Adaptor
- G. 1 x 54936 1.0 Litre Latex Free Disposable Rebreathing Bag (CP)
- H. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

P3E72-102 Circuit - Length = 1.8m (72")


GTIN: 9330783038879

P3E120-102 Circuit - Length = 3.0m (120")

GTIN: 9330783038886

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PC) = Polycarbonate (PP) = Polypropylene (CP) = Chloroprene (***) = Tubing Length

Adult Anaesthetic Corrugated Circuit with ADU Connector - Disposable



Disposable Adult Anaesthetic Circuit with ADU Connector

Consisting of:

- A. 4 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x APM51046 Fixed Elbow 15mm M to 22mm M/15mm F (PP)
- E. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- F. 1 x PH3100(D) 22mm M/M Adaptor (PP)
- G. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- H. 1 x ADU1016 ADU Connector (PP)
- I. 1 x MP-CO22N48 22mm ID Natural Disposable Corrugated Tubing x 1.2m (48") (PP)


A4C60-108-A Circuit - Length = 1.5m (60")
Ventilation Limb Length = 1.2m (48")

GTIN: 9330783039630

A4C120-108-A Circuit - Length = 3.0m (120")
Ventilation Limb Length = 1.2m (48")

GTIN: 9330783039647

Paediatric Anaesthetic Corrugated Circuit with ADU Connector - Disposable



Disposable Paediatric Anaesthetic Circuit with ADU Connector

Consisting of:

- A. 3 x APM240050002(PC) 22mm F to 15mm Tube Adaptors (PP)
- B. 2 x MP-CO15N*** 15mm ID Natural Disposable Corrugated Tubing (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F Adaptor (PP)
- D. 1 x APM51046 Fixed Elbow 15mm M to 22mm M/15mm F (PP)
- E. 1 x MP-CO22B*** 22mm ID Blue Disposable Corrugated Tubing (PP)
- F. 1 x PH3100(D) 22mm M/M Adaptor (PP)
- G. 1 x 54936 1 Litre Latex Free Disposable Rebreathing Bag (CP)
- H. 1 x ADU1016 ADU Connector (PP)
- I. 1 x MP-CO15N48 15mm ID Natural Disposable Corrugated Tubing x 1.2m (48") (PP)
- J. 1 x A2200(PC) 15mm ID tube to 22mm M (PC)

P4C60-108-A Circuit - Length = 1.5m (60")
Ventilation Limb Length = 1.2m (48")


GTIN: 9330783039654

P4C120-108-A Circuit - Length = 3.0m (120")
Ventilation Limb Length = 1.2m (48")

GTIN: 9330783040391

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (PC) = Polycarbonate (***) = Tubing Length

Adult Anaesthetic Expandable Circuit with ADU Connector Disposable



Disposable Adult Anaesthetic Circuit with ADU Connector

Consisting of:

- A. 4 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x APM51046 Fixed Elbow 15mm M to 22mm M/15mm F (PP)
- E. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- F. 1 x PH3100(D) 22mm M/M Adaptor (PP)
- G. 1 x 54937 2 Litre Latex Free Disposable Rebreathing Bag (CP)
- H. 1 x ADU1016 ADU Connector (PP)
- I. 1 x MP-EX22N48 22mm ID Natural Disposable Expandable Tubing x 1.2m (48") (PP)


A4E60-108-A Circuit - Length = 1.5m (60")
Ventilation Limb Length = 1.2m (48")

GTIN: 9330783040087

A4E120-108-A Circuit - Length = 3.0m (120")
Ventilation Limb Length = 1.2m (48")

GTIN: 9330783040094

Paediatric Anaesthetic Expandable Circuit with ADU Connector Disposable



Disposable Paediatric Anaesthetic Circuit with ADU Connector

Consisting of:

- A. 3 x APM240050002(PC) 22mm F to 15mm Tube Adaptors (PP)
- B. 2 x MP-EX15N*** 15mm ID Natural Disposable Expandable Tubing (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F Adaptor (PP)
- D. 1 x APM51046 Fixed Elbow 15mm M to 22mm M/15mm F (PP)
- E. 1 x MP-EX22B*** 22mm ID Blue Disposable Expandable Tubing (PP)
- F. 1 x PH3100(D) 22mm M/M Adaptor (PP)
- G. 1 x 54936 1 Litre Latex Free Disposable Rebreathing Bag (CP)
- H. 1 x ADU1016 ADU Connector (PP)
- I. 1 x MP-EX15N48 15mm ID Natural Disposable Expandable Tubing x 1.2m (48") (PP)
- J. 1 x A2200(PC) 15mm ID tube to 22mm M (PC)

P4E60-108-A Circuit - Length = 1.5m (60")
Ventilation Limb Length = 1.2m (48")


GTIN: 9330783040100

P4E120-108-A Circuit - Length = 3.0m (120")
Ventilation Limb Length = 1.2m (48")

GTIN: 9330783040384

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (CP) = Chloroprene (PC) = Polycarbonate (***) = Tubing Length

Corrugated Breathing Circuits - Disposable



Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2C60-101 Circuit - Length = 1.5m (60")

GTIN: 9330783037551

A2C72-101 Circuit - Length = 1.8m (72")


GTIN: 9330783037568

A2C96-101 Circuit - Length = 2.4m (96")

GTIN: 9330783037575

A2C120-101 Circuit - Length = 3.0m (120")

GTIN: 9330783037582



Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51100L Fixed 'Y' M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2C60-201 Circuit - Length = 1.5m (60")

GTIN: 9330783037599

A2C72-201 Circuit - Length = 1.8m (72")


GTIN: 9330783037605

A2C96-201 Circuit - Length = 2.4m (96")

GTIN: 9330783037612

A2C120-201 Circuit - Length = 3.0m (120")

GTIN: 9330783037629



Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2C60-301 Circuit - Length = 1.5m (60")

GTIN: 9330783037636

A2C72-301 Circuit - Length = 1.8m (72")


GTIN: 9330783037643

A2C96-301 Circuit - Length = 2.4m (96")

GTIN: 9330783037650

A2C120-301 Circuit - Length = 3.0m (120")

GTIN: 9330783037667



Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2C60-401 Circuit - Length = 1.5m (60")

GTIN: 9330783037674

A2C72-401 Circuit - Length = 1.8m (72")

GTIN: 9330783037681

A2C96-401 Circuit - Length = 2.4m (96")

GTIN: 9330783037698

A2C120-401 Circuit - Length = 3.0m (120")

GTIN: 9330783037704


Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (HDPE) = High Density Polyethylene (***) = Tubing Length

Corrugated Breathing Circuits - Disposable

Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x MP1000L Split 'Y' M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)



A2C60-501 Circuit - Length = 1.5m (60")

GTIN: 9330783037711

A2C72-501 Circuit - Length = 1.8m (72")

GTIN: 9330783037728

A2C96-501 Circuit - Length = 2.4m (96")

GTIN: 9330783037735


A2C120-501 Circuit - Length = 3.0m (120")

GTIN: 9330783037742

Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)



A2C60-601 Circuit - Length = 1.5m (60")

GTIN: 9330783037759

A2C72-601 Circuit - Length = 1.8m (72")

GTIN: 9330783037766

A2C96-601 Circuit - Length = 2.4m (96")

GTIN: 9330783037773


A2C120-601 Circuit - Length = 3.0m (120")

GTIN: 9330783037780

Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' M/M/M Adaptor (HDPE)
- D. 1 x APM8821 22mm F Plastic Cap (PP)



A2C60-701 Circuit - Length = 1.5m (60")

GTIN: 9330783037797

A2C72-701 Circuit - Length = 1.8m (72")

GTIN: 9330783037803

A2C96-701 Circuit - Length = 2.4m (96")

GTIN: 9330783037810


A2C120-701 Circuit - Length = 3.0m (120")

GTIN: 9330783037827

Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51098 Swivel 'Y' M/M/M Adaptor with Luer Port and Cap (HDPE)
- D. 1 x APM8821 22mm F Plastic Cap (PP)



A2C60-801 Circuit - Length = 1.5m (60")

GTIN: 9330783037834

A2C72-801 Circuit - Length = 1.8m (72")

GTIN: 9330783037841

A2C96-801 Circuit - Length = 2.4m (96")


GTIN: 9330783037858

A2C120-801 Circuit - Length = 3.0m (120")

GTIN: 9330783037865

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (HDPE) = High Density Polyethylene (***) = Tubing Length

Corrugated Breathing Circuits - Disposable




Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2C60-102 Circuit - Length = 1.5m (60")	GTIN: 9330783039715
A2C72-102 Circuit - Length = 1.8m (72")	GTIN: 9330783039722
A2C96-102 Circuit - Length = 2.4m (96")	GTIN: 9330783039739
A2C120-102 Circuit - Length = 3.0m (120")	GTIN: 9330783039746




Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2C60-302 Circuit - Length = 1.5m (60")	GTIN: 9330783039753
A2C72-302 Circuit - Length = 1.8m (72")	GTIN: 9330783039760
A2C96-302 Circuit - Length = 2.4m (96")	GTIN: 9330783039777
A2C120-302 Circuit - Length = 3.0m (120")	GTIN: 9330783039784




Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2C60-402 Circuit - Length = 1.5m (60")	GTIN: 9330783039791
A2C72-402 Circuit - Length = 1.8m (72")	GTIN: 9330783039807
A2C96-402 Circuit - Length = 2.4m (96")	GTIN: 9330783039814
A2C120-402 Circuit - Length = 3.0m (120")	GTIN: 9330783039821




Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2C60-602 Circuit - Length = 1.5m (60")	GTIN: 9330783039838
A2C72-602 Circuit - Length = 1.8m (72")	GTIN: 9330783039845
A2C96-602 Circuit - Length = 2.4m (96")	GTIN: 9330783039852
A2C120-602 Circuit - Length = 3.0m (120")	GTIN: 9330783039869



Adult Disposable Circuit


Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-CO22N*** 22mm ID Natural Corrugated Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' M/M/M Adaptor (HDPE)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2C60-702 Circuit - Length = 1.5m (60")	GTIN: 9330783039876
A2C72-702 Circuit - Length = 1.8m (72")	GTIN: 9330783039883
A2C96-702 Circuit - Length = 2.4m (96")	GTIN: 9330783039890
A2C120-702 Circuit - Length = 3.0m (120")	GTIN: 9330783039906

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (HDPE) = High Density Polyethylene (***) = Tubing Length

Expandable Breathing Circuits - Disposable



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-101 Circuit - Length = 1.8m (72")

GTIN: 9330783037872

A2E120-101 Circuit - Length = 3.0m (120")

GTIN: 9330783037889



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51100L Fixed 'Y' M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-201 Circuit - Length = 1.8m (72")

GTIN: 9330783037896

A2E120-201 Circuit - Length = 3.0m (120")

GTIN: 9330783037902



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-301 Circuit - Length = 1.8m (72")

GTIN: 9330783037919

A2E120-301 Circuit - Length = 3.0m (120")

GTIN: 9330783037926



Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-401 Circuit - Length = 1.8m (72")


GTIN: 9330783037933

A2E120-401 Circuit - Length = 3.0m (120")

GTIN: 9330783037940

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (HDPE) = High Density Polyethylene (***) = Tubing Length

Expandable Breathing Circuits - Disposable



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x MP1000L Split 'Y' M/M/M Adaptor with Luer Port and Plug (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-501 Circuit - Length = 1.8m (72")

GTIN: 9330783037957

A2E120-501 Circuit - Length = 3.0m (120")

GTIN: 9330783037964



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-601 Circuit - Length = 1.8m (72")

GTIN: 9330783037971

A2E120-601 Circuit - Length = 3.0m (120")

GTIN: 9330783037988



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' M/M/M Adaptor (HDPE)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-701 Circuit - Length = 1.8m (72")

GTIN: 9330783037995

A2E120-701 Circuit - Length = 3.0m (120")

GTIN: 9330783038008



Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51098 Swivel 'Y' M/M/M Adaptor with Luer Port and Cap (HDPE)
- D. 1 x APM8821 22mm F Plastic Cap (PP)

A2E72-801 Circuit - Length = 1.8m (72")


GTIN: 9330783038015

A2E120-801 Circuit - Length = 3.0m (120")

GTIN: 9330783038022

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (HDPE) = High Density Polyethylene (***) = Tubing Length

Expandable Breathing Circuits - Disposable



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51100 Fixed 'Y' M/M/M Adaptor (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2E72-102 Circuit - Length = 1.8m (72")

GTIN: 9330783039913

A2E120-102 Circuit - Length = 3.0m (120")

GTIN: 9330783039920



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51100TP Fixed 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2E72-302 Circuit - Length = 1.8m (72")

GTIN: 9330783039937

A2E120-302 Circuit - Length = 3.0m (120")

GTIN: 9330783039944



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x MP1000 Split 'Y' M/M/M Adaptor (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2E72-402 Circuit - Length = 1.8m (72")

GTIN: 9330783039951

A2E120-402 Circuit - Length = 3.0m (120")

GTIN: 9330783039968



Adult Disposable Circuit

Consisting of:


- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x MP1000TP Split 'Y' M/M/M Adaptor with Plugged Temp. and Pres. Ports (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2E72-602 Circuit - Length = 1.8m (72")

GTIN: 9330783039975

A2E120-602 Circuit - Length = 3.0m (120")

GTIN: 9330783039982



Adult Disposable Circuit

Consisting of:

- A. 2 x PH3000(D) 22mm M to 22mm F Adaptors (PP)
- B. 2 x MP-EX22N*** 22mm ID Natural Expandable Tubing (PP)
- C. 1 x APM51093 Swivel 'Y' M/M/M Adaptor (HDPE)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)

A2E72-702 Circuit - Length = 1.8m (72")

GTIN: 9330783039999

A2E120-702 Circuit - Length = 3.0m (120")

GTIN: 9330783040001


Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (SY) = Styrene (HDPE) = High Density Polyethylene (***) = Tubing Length

Paediatric Breathing Circuits - Disposable

Disposable Paediatric Anaesthetic Circuit - Corrugated

Consisting of:

- A. 2 x APM240050002(PC) 15mm ID tube to 22mm F (PC)
- B. 2 x MP-CO15N*** 15mm ID Natural Corrugated Disposable Tube (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F/22mm M (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)



P2C60-102 Circuit - Length = 1.5m (60")

GTIN: 9330783040018

P2C72-102 Circuit - Length = 1.8m (72")

GTIN: 9330783040025

P2C96-102 Circuit - Length = 2.4m (96")

GTIN: 9330783040032


P2C120-102 Circuit - Length = 3.0m (120")

GTIN: 9330783040049

Disposable Paediatric Anaesthetic Circuit - Expandable

Consisting of:

- A. 2 x APM240050002(PC) 15mm ID tube to 22mm F (PC)
- B. 2 x MP-EX15N*** 15mm ID Natural Expandable Disposable Tube (PP)
- C. 1 x 5112 Fixed 'Y' 15mm M/15mm M to 15mm F/22mm M (PP)
- D. 1 x APM51033 Disposable Elbow 15mm M to 22mm M/15mm F with Luer Port (PP)



P2E72-102 Circuit - Length = 1.8m (72")


GTIN: 9330783040056

P2E120-102 Circuit - Length = 3.0m (120")

GTIN: 9330783040063

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PC) = Polycarbonate (PP) = Polypropylene (***) = Tubing Length

Coaxial Circuits



Coaxial Circuits

The coaxial breathing circuit is supplied with a detachable proximal fitting for providing respiratory gases to and receiving expiratory gases from a patient. The circuit contains adaptor fittings for connection to patient devices or assisted ventilation systems.



Coaxial Circuits

PH70159 Disposable Paediatric Coaxial Circuit, Length = 1.5m with Mask Elbow and Luer Port (PE)(PP) GTIN: 9330783036943

PH70159-1.8 Disposable Paediatric Coaxial Circuit, Length = 1.8m with Mask Elbow and Luer Port (PE)(PP) GTIN: 9330783036950


PH70125 Disposable Adult Coaxial Circuit, Length = 1.5m with Mask Elbow and Luer Port (PE)(PP) GTIN: 9330783036592

PH70125-1.8 Disposable Adult Coaxial Circuit, Length = 1.8m with Mask Elbow and Luer Port (PE)(PP) GTIN: 9330783036608

PH70125-2.4 Disposable Adult Coaxial Circuit, Length = 2.4m with Mask Elbow and Luer Port (PE)(PP) GTIN: 9330783030934

PH70125-3.0 Disposable Adult Coaxial Circuit, Length = 3.0m with Mask Elbow and Luer Port (PE)(PP) GTIN: 9330783030101

PH70125-3.3 Disposable Adult Coaxial Circuit, Length = 3.3m with Mask Elbow and Luer Port (PE)(PP) GTIN: 9330783028979



Coaxial Circuits with Bagging Limb

PH70159-B1L Paediatric Coaxial Circuit, Length = 1.5m with 1.5m Blue Limb and 1Lit. Breathing Bag (PE)(PP)(CP) GTIN: 9330783036967

PH70159-1.8B1L Paediatric Coaxial Circuit, Length = 1.8m with 1.8m Blue Limb and 1Lit. Breathing Bag (PE)(PP)(CP) GTIN: 9330783036974


PH70125-B2L Adult Coaxial Circuit, Length = 1.5m with 1.5m Blue Limb and 2Lit. Breathing Bag (PE)(PP)(CP) GTIN: 9330783031023

PH70125-1.8B2L Adult Coaxial Circuit, Length = 1.8m with 1.8m Blue Limb and 2Lit. Breathing Bag (PE)(PP)(CP) GTIN: 9330783030095

PH70125-2.4B2L Adult Coaxial Circuit, Length = 2.4m with 2.4m Blue Limb and 2Lit. Breathing Bag (PE)(PP)(CP) GTIN: 9330783030941


PH70125-3B2L Adult Coaxial Circuit, Length = 3.0m with 3.0m Blue Limb and 2Lit. Breathing Bag (PE)(PP)(CP) GTIN: 9330783030118

PH70125-3.3B2L Adult Coaxial Circuit, Length = 3.3m with 3.3m Blue Limb and 2Lit. Breathing Bag (PE)(PP)(CP) GTIN: 9330783036981



Legend: (PVC) = Polyvinyl Chloride (PE) = Polyethylene (PP) = Polypropylene (CP) = Chlorprene

Resuscitation Breathing Circuits - Disposable




PHC/3500 Disposable Adult Resuscitation Circuit

GTIN: 9330783019908

Consisting of:

- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x 5126 Norman Mask Elbow 15mm M to 22mm M/15mm F (PP)
- C. 2 x PH3100(D) Adaptor 22mm M/15mm F to 22mm M/15mm F (PP)
- D. 1 x MP-CO22N12 22mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- E. 1 x 54937 2.0 lit. Disposable Rebreathing Bag (CP)



PHC/3520 Disposable Paediatric Resuscitation Circuit

GTIN: 9330783019939


Consisting of:

- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x 5126 Norman Mask Elbow 15mm M to 22mm M/15mm F (PP)
- C. 2 x PH3100(D) Adaptor 22mm M/15mm F to 22mm M/15mm F (PP)
- D. 1 x MP-CO22N12 22mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- E. 1 x 54935 0.5 lit. Disposable Rebreathing Bag (CP)

PHC/3520/1L Disposable Paediatric Resus. Circuit

GTIN: 9330783019953

(As above circuit with 1.0 lit. bag in lieu of 0.5 lit. bag) (CP)




PHC/3550/1L/SNIP Disposable Paediatric Resus. Circuit

GTIN: 9330783025145

Consisting of:

- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x 5126 Norman Mask Elbow 15mm M to 22mm M/15mm F (PP)
- C. 1 x 5114 Adaptor 15mm M/F with 6mm Outlet (PP)
- D. 1 x MP-CO15N12 15mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- E. 1 x A2200(PC) Adaptor 22mm M/15mm F to 15mm M/10mm F (PP)
- F. 1 x 54936 1.0 lit. Disposable Rebreathing Bag (Endtail Snipped) (CP)



PHC/3560 Disposable Paediatric Resuscitation Circuit


GTIN: 9330783020072

Consisting of:

- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x PH3300(D) Adaptor 22mm F to 6mm Tube (PP)
- C. 1 x APM51033 Mask Elbow with Luer Port and Plug 15mm M to 22mm M/15mm F (PP)
- D. 1 x 5114 Adaptor 15mm M/F with 6mm Outlet (PP)
- E. 1 x MP-CO15N12 15mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- F. 1 x A2200(PC) Adaptor 22mm M/15mm F to 15mm M/10mm F (PP)
- G. 1 x 54935 0.5 lit. Disposable Rebreathing Bag (CP)

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (CP) = Chloroprene (HDPE) = High Density Polyethylene (PC) = Polycarbonate


Resuscitation Breathing Circuits - Disposable



PHC/3540/A/SNIP Disposable Paediatric Resus. Circuit GTIN: 9330783020003

Consisting of:


- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x PH3300(D) Adaptor 22mm F to 6mm Tube (PP)
- C. 1 x 5126 Norman Mask Elbow 15mm M to 22mm M/15mm F (PP)
- D. 1 x APM51035 Adaptor 15mm M/F with Luer Port and Cap (PP)
- E. 1 x MP-CO15N24 15mm ID Natural Corrugated Tubing x 600mm (24") Length (PP)
- F. 1 x A2200(PC) Adaptor 22mm M/15mm F to 15mm M/10mm F (PP)
- G. 1 x 54935 0.5 lit. Disposable Rebreathing Bag (Endtail Snipped) (CP)



PHC/3540 Disposable Neonatal Resuscitation Circuit GTIN: 9330783019984

Consisting of:


- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x 5126 Norman Mask Elbow 15mm M to 22mm M/15mm F (PP)
- C. 1 x MP-CO15N12 15mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- D. 1 x A2200(PC) Adaptor 22mm M/15mm F to 15mm M/10mm F (PP)
- E. 1 x 54935 0.5 lit. Disposable Rebreathing Bag (CP)



PHC/3540/WT Disposable Neonatal Resus. Circuit GTIN: 9330783020034

Consisting of:


- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x 5114 Washington 'T' 15mm M to 15mm F (PP)
- C. 1 x MP-CO15N12 15mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- D. 1 x A2200(PC) Adaptor 22mm M/15mm F to 15mm M/10mm F (PP)
- E. 1 x 54935 0.5 lit. Disposable Rebreathing Bag (CP)
- F. 1 x PH3300(D) Adaptor 22mm F to 6mm Tube (PP)



PHC/3540/M-1 Disposable Neonatal Resus. Circuit GTIN: 9330783020027

Consisting of:

- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x 5126 Norman Mask Elbow 15mm M to 22mm M/15mm F (PP)
- C. 1 x MP-CO15N12 15mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- D. 1 x A2200(PC) Adaptor 22mm M/15mm F to 15mm M/10mm F (PP)
- E. 1 x PH20117 Inflatable Mask Size #1 (PVC) (HDPE)
- F. 1 x 54935 0.5 lit. Disposable Rebreathing Bag (CP)




PHC/3550 Disposable Neonatal Resuscitation Circuit GTIN: 9330783020058

Consisting of:

- A. 1 x 34502.1/1 Oxygen Tubing x 2.1m (PVC)
- B. 1 x 5126 Norman Mask Elbow 15mm M to 22mm M/15mm F (PP)
- C. 1 x 5114 Adaptor 15mm M/F with 6mm Outlet (PP)
- D. 1 x MP-CO15N12 15mm ID Natural Corrugated Tubing x 300mm (12") Length (PP)
- E. 1 x A2200(PC) Adaptor 22mm M/15mm F to 15mm M/10mm F (PP)
- F. 1 x 54935 0.5 lit. Disposable Rebreathing Bag (CP)

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (PP) = Polypropylene (CP) = Chloroprene (HDPE) = High Density Polyethylene (PC) = Polycarbonate

Anaesthetic Breathing Circuits - Reusable




5200P.R. Paediatric Anaesthetic Circuit

GTIN: 9330783006496

Consisting of:

- A. 3 x APM240050002(PS) Adaptor 22mm F to 15mm Tube (PP)
- B. 3 x SB405-15-72 15mm ID x 1800mm Smoothbor Tubing (S)
- C. 1 x APM6103/MF/22 Adaptor 22mm M/15mm F/M Tube/M Tube Adaptor (PC)
- D. 1 x APM50258 0.5 Litre Latex Free Rebreathing Bag (CP)
- E. 1 x A2200(PS) Adaptor 22mm M to 15mm M (PS)

Tubing Length = 1.8m (72")




5200R Adult Anaesthetic Circuit

GTIN: 9330783006519

Consisting of:

- A. 3 x A2300(PS) Adaptor 22mm F/15mm F X 22mm M (PS)
- B. 2 x SB403-60 Smoothbore Tubing 22mm ID X 1500mm (S)
- C. 1 x A2000Y(PS) 'Y' Piece 22mm M/15mm F X 22mm M X 22mm M (PS)
- D. 1 x SB403-60BC Blue Cuff Smoothbor Tubing 22mm ID X 1500mm (S)
- E. 1 x A90001R(PS) Adaptor 22mm M/15mm F X 22mm M/15mm F (PS)
- F. 1 x APM50260 2.0 lit. Rebreathing Bag (CP)

Tubing Length = 1.5m (60")




5205 PAED Paediatric Resuscitation Circuit

GTIN: 9330783006564

Consisting of:

- A. 1 x 00-234 Norman Mask Elbow 15mm F/22mm M X 15mm M (LWM)
- B. 1 x 00-240 Washington 'T' Piece 15mm M/F (LWM)
- C. 1 x SB405-15-12 Smoothbor Tubing 15mm ID X 300mm (S)
- D. 1 x A2200(PS) Adaptor 22mm M/15mm F X 15mm M/10mm F (PS)
- E. 1 x APM50258 0.5 lit. Rebreathing Bag (CP)
- F. 1 x 00-274 Gauge & Relief Assembly (M)




5205 P.RE Paediatric Resuscitation Circuit

GTIN: 9330783006571

Consisting of:

- A. 1 x 00-234 Norman Mask Elbow 15mm F/22mm M X 15mm M (LWM)
- B. 1 x 00-240 Washington 'T' Piece 15mm M/F (LWM)
- C. 1 x SB405-15-12 Smoothbor Tubing 15mm ID X 300mm (S)
- D. 1 x A2200(PS) Adaptor 22mm M/15mm F X 15mm M/10mm F (PS)
- E. 1 x APM50258 0.5 lit. Rebreathing Bag (CP)



5206 P.RE Paediatric Resuscitation Circuit


GTIN: 9330783006588

Consisting of:

- A. 1 x 00-234 Norman Mask Elbow 15mm F/22mm M X 15mm M (LWM)
- B. 1 x SB405-15-12 Smoothbor Tubing 15mm ID X 300mm (S)
- C. 1 x A2200(PS) Adaptor 22mm M/15mm F X 15mm M/10mm F (PS)
- D. 1 x APM50258 0.5 lit. Rebreathing Bag (CP)

Legend: M = Male (External) F = Female (Internal) ID = Inside Diameter (CP) = Chloroprene (PS) = Polysulfone (LWM) = Light Weight Metal (M) = Metal

Disposable Mid-O-Gas / Analgesia Breathing Circuits



MC/3200 Disposable Breathing Circuit

- A. MC/HP Handpiece 22mm M with one-way valve (SY)(S)
- B. 56001 Mouthpiece 22mm F (HDPE)
- C. MP-EX22B120 22mm ID Blue Expandable Tubing (PP)

GTIN: 9330783000210


Tube length = 3.0 metres

MC/3300 Disposable Breathing Circuit

- A. MC/HP Handpiece 22mm M with one-way valve (SY)(S)
- B. 56001 Mouthpiece 22mm F (HDPE)
- C. MP-EX22B60 22mm ID Blue Expandable Tubing (PP)

GTIN: 9330783000326

Tube length = 1.5 metres



MC/3200/F Disposable Breathing Circuit

- A. MC/HP Handpiece 22mm M with one-way valve (SY)(S)
- B. 55000104-500 In-line Filter 22mm M/F (PP)
- C. 56001 Mouthpiece 22mm F (HDPE)
- D. MP-EX22B120 22mm ID Blue Expandable Tubing (PP)

GTIN: 9330783000227


Tube length = 3.0 metres

MC/3300/F Disposable Breathing Circuit

- A. MC/HP Handpiece 22mm M with one-way valve (SY)(S)
- B. 55000104-500 In-line Filter 22mm M/F (PP)
- C. 56001 Mouthpiece 22mm F (HDPE)
- D. MP-EX22B60 22mm ID Blue Expandable Tubing (PP)

GTIN: 9330783000203

Tube length = 1.5 metres




MC/3200/F/A Disposable Breathing Circuit

- A. MC/HP Handpiece 22mm M with one-way valve (SY)(S)
- B. 55000104-500 In-line Filter 22mm M/F (PP)
- C. 56001 Mouthpiece 22mm F (HDPE)
- D. MP-EX22B120 22mm ID Blue Expandable Tubing (PP)
- E. PH3000(D) Adaptor 22mm M to 22mm F (PP)

GTIN: 9330783000234

Tube length = 3.0 metres



MC/3200/F/M* Disposable Breathing Circuits

- A. MC/HP Handpiece 22mm M with one-way valve (SY)(S)
- B. 55000104-500 In-line Filter 22mm M/F (PP)
- C. Inflatable Mask with Hook, Size: 2 or 4 (refer below) (PVC)(SS)(R)(PP)
- D. MP-EX22B120 22mm ID Blue Expandable Tubing (PP)

Tube length = 3.0 metres

MC/3200/F/M2 Circuit, Inflatable Mask with Hook, Size: 2


GTIN: 9330783017867

MC/3200/F/M4 Circuit, Inflatable Mask with Hook, Size: 4

GTIN: 9330783017874

Legend: ID = Internal Diameter (PP) = Polypropylene (SY) = Styrene (PVC) = Polyvinyl chloride (HDPE) = High Density Polyethylene (SS) = Stainless Steel (R) = Rubber (S) = Silicone


Disposable Mid-O-Gas / Analgesia Breathing Circuits



MC/3000 Disposable Breathing Circuit GTIN: 9330783000241
 A. MP-EX22P120 22mm ID Pink Expandable Tubing (PP)
 B. MC-HP Handpiece 22mm M/M/M with one-way valve (SY)(S)
 C. 56001 Mouthpiece 22mm F (HDPE)
 D. MP-EX22B120 22mm ID Blue Expandable Tubing (PP)

MC/3000/SA Disposable Breathing Circuit GTIN: 9330783017850
 (with Scavenge Adaptor)
 Above circuit including Scavenge Adaptor.
 E. PH3200(D) Adaptor 22mm M to 19mm M (PP)

Tube lengths = 3.0 metres



MC/3000/F Disposable Breathing Circuit GTIN: 9330783000258
 A. MP-EX22P120 22mm ID Pink Expandable Tubing (PP)
 B. MC-HP Handpiece 22mm M/M/M with one-way valve (SY)(S)
 C. 55000104-500 In-line Filter 22mm M/F (PP)
 D. 56001 Mouthpiece 22mm F (HDPE)
 E. MP-EX22B120 22mm ID Blue Expandable Tubing (PP)


MC/3000/F/SA Disposable Breathing Circuit GTIN: 9330783017843
 (with Scavenge Adaptor)
 Above circuit including Scavenge Adaptor.
 F. PH3200(D) Adaptor 22mm M to 19mm M (PP)

MC/3000/F/M* Disposable Breathing Circuit
 Above circuit including Mask in lieu of mouthpiece.
 * Mask size must be stated at ordering (refer below).
 G. Mask: (PVC)(SS)(R)(PP)


Mask Sizes: #2 - Child
 #3 - Large Child
 #4 - Adult
 #5 - Large Adult

MC/3000/F/M2 Circuit with Size #2 Mask GTIN: 9330783025954
MC/3000/F/M3 Circuit with Size #3 Mask GTIN: 9330783031375
MC/3000/F/M4 Circuit with Size #4 Mask GTIN: 9330783031382
MC/3000/F/M5 Circuit with Size #5 Mask GTIN: 9330783031399

Tube lengths = 3.0 metres



MC/3500 Disposable Circuit Components GTIN: 9330783017881
 A. MC-HP Handpiece 22mm M/M/M with one-way valve (SY)(S)
 B. 56001 Mouthpiece 22mm F (HDPE)



MC/3500/F Disposable Circuit Components GTIN: 9330783017898
 A. MC-HP Handpiece 22mm M/M/M with one-way valve (SY)(S)
 B. 55000104-500 In-line Filter 22mm M/F (PP)
 C. 56001 Mouthpiece 22mm F (HDPE)

Legend: (PP) = Polypropylene (SY) = Styrene (PVC) = Polyvinyl chloride (HDPE) = High Density Polyethylene (SS) = Stainless Steel (R) = Rubber (S) = Silicone

Rebreathing Bags



Disposable Rebreathing Bags

54935	0.5 Litre Capacity, Ferrule Neck (22mm F) (CP)	GTIN: 9330783006670
54936	1.0 Litre Capacity, Ferrule Neck (22mm F) (CP)	GTIN: 9330783006687
54937	2.0 Litre Capacity, Ferrule Neck (22mm F) (CP)	GTIN: 9330783006694
54938	3.0 Litre Capacity, Ferrule Neck (22mm F) (CP)	GTIN: 9330783029136



Reusable Rebreathing Bags

APM50258	0.5 Litre Capacity, Ferrule Neck (22mm F) with tail and loop (CP)	GTIN: 9330783009374
APM50259	1.0 Litre Capacity, Ferrule Neck (22mm F) with tail and loop (CP)	GTIN: 9330783009381
APM50260	2.0 Litre Capacity, Ferrule Neck (22mm F) with tail and loop (CP)	GTIN: 9330783009398
APM50261	3.0 Litre Capacity, Ferrule Neck (22mm F) with tail and loop (CP)	GTIN: 9330783009404




Reusable Silicone Rebreathing Bags

PH70501	0.5 Litre Capacity, Ferrule Neck (22mm F) (S)	GTIN: 9330783028986
PH70502	1.0 Litre Capacity, Ferrule Neck (22mm F) (S)	GTIN: 9330783028993
PH70503	2.0 Litre Capacity, Ferrule Neck (22mm F) (S)	GTIN: 9330783029006
PH70504	3.0 Litre Capacity, Ferrule Neck (22mm F) (S)	GTIN: 9330783029013

Legend: M = Male (Outside) F = Female (Inside) (D) = Disposable ID = Inside Diameter (CP) = Chloroprene (S) = Silicone


Corrugated Disposable Plastic Tubing



10mm ID (PP)


Part No.	Colour	Length	GTIN
MP-C010B30M-12	Blue	30m Roll	9330783038893
MP-C010N30M-12	Natural	30m Roll	9330783038909
MP-C010P30M-12	Pink	30m Roll	9330783038916

NOTE: Cuff Every 300mm (12")



15mm ID (PP)

Part No.	Colour	Length	GTIN
MP-C015B12	Blue	300mm (12")	9330783038923
MP-C015B60	Blue	1.5m (60")	9330783038930
MP-C015B72	Blue	1.8m (72")	9330783038947
MP-C015B96	Blue	2.4m (96")	9330783038954
MP-C015B120	Blue	3.0m (120")	9330783038961
MP-C015B30M-12	Blue	30m Roll (12" cuff)	9330783038978
MP-C015N12	Natural	300mm (12")	9330783038985
MP-C015N60	Natural	1.5m (60")	9330783038992
MP-C015N72	Natural	1.8m (72")	9330783039005
MP-C015N96	Natural	2.4m (96")	9330783039012
MP-C015N120	Natural	3.0m (120")	9330783039029
MP-C015N30M-12	Natural	30m Roll (12" cuff)	9330783039036
MP-C015P12	Pink	300mm (12")	9330783039043
MP-C015P60	Pink	1.5m (60")	9330783039050
MP-C015P72	Pink	1.8m (72")	9330783039067
MP-C015P96	Pink	2.4m (96")	9330783039074
MP-C015P120	Pink	3.0m (120")	9330783039081
MP-C015P30M-12	Pink	30m Roll (12" cuff)	9330783039098



22mm ID (PP)

Part No.	Colour	Length	GTIN
MP-C022B12	Blue	300mm (12")	9330783039104
MP-C022B60	Blue	1.5m (60")	9330783039111
MP-C022B72	Blue	1.8m (72")	9330783039128
MP-C022B96	Blue	2.4m (96")	9330783039135
MP-C022B120	Blue	3.0m (120")	9330783039142
MP-C022B30M-12	Blue	30m Roll (12" cuff)	9330783039159
MP-C022N12	Natural	300mm (12")	9330783039166
MP-C022N60	Natural	1.5m (60")	9330783039173
MP-C022N72	Natural	1.8m (72")	9330783039180
MP-C022N96	Natural	2.4m (96")	9330783039197
MP-C022N120	Natural	3.0m (120")	9330783039203
MP-C022N30M-12	Natural	30m Roll (12" cuff)	9330783039210
MP-C022P12	Pink	300mm (12")	9330783039227
MP-C022P60	Pink	1.5m (60")	9330783039234
MP-C022P72	Pink	1.8m (72")	9330783039241
MP-C022P96	Pink	2.4m (96")	9330783039258
MP-C022P120	Pink	3.0m (120")	9330783039265
MP-C022P30M-12	Pink	30m Roll (12" cuff)	9330783039272

Mediplas	Corrugated	Tubing ID (mm)	Tubing Colour	Tubing Length (inches)
MP	CO	22	B	72
B=Blue, N=Natural, P=Pink				
12", 60", 72", 96", 120", 30m Roll				

Other lengths available upon request

Legend: (D) = Disposable ID = Inside Diameter (PP) = Polypropylene

Expandable Disposable Plastic Tubing


15mm ID Expandable Tubing (PP)



MP-EX15N6

Part No.	Colour	Length	GTIN
MP-EX15B60	Blue	1.5m (60")	9330783037179
MP-EX15B72	Blue	1.8m (72")	9330783037186
MP-EX15B96	Blue	2.4m (96")	9330783037193
MP-EX15B120	Blue	3.0m (120")	9330783037209
MP-EX15N6	Natural	150mm (6")	9330783037254
MP-EX15N60	Natural	1.5m (60")	9330783037216
MP-EX15N72	Natural	1.8m (72")	9330783037223
MP-EX15N96	Natural	2.4m (96")	9330783037230
MP-EX15N120	Natural	3.0m (120")	9330783037247
MP-EX15P60	Pink	1.5m (60")	9330783037261
MP-EX15P72	Pink	1.8m (72")	9330783037278
MP-EX15P96	Pink	2.4m (96")	9330783037285
MP-EX15P120	Pink	3.0m (120")	9330783037292
MP-EX15G60	Green	1.5m (60")	9330783037308
MP-EX15G72	Green	1.8m (72")	9330783037315
MP-EX15G96	Green	2.4m (96")	9330783037322
MP-EX15G120	Green	3.0m (120")	9330783037339

22mm ID Expandable Tubing (PP)




Part No.	Colour	Length	GTIN
MP-EX22B60	Blue	1.5m (60")	9330783037346
MP-EX22B72	Blue	1.8m (72")	9330783037353
MP-EX22B96	Blue	2.4m (96")	9330783037360
MP-EX22B120	Blue	3.0m (120")	9330783037377
MP-EX22N60	Natural	1.5m (60")	9330783037384
MP-EX22N72	Natural	1.8m (72")	9330783037391
MP-EX22N96	Natural	2.4m (96")	9330783037407
MP-EX22N120	Natural	3.0m (120")	9330783037414
MP-EX22P60	Pink	1.5m (60")	9330783037421
MP-EX22P72	Pink	1.8m (72")	9330783037438
MP-EX22P96	Pink	2.4m (96")	9330783037445
MP-EX22P120	Pink	3.0m (120")	9330783037452
MP-EX22G60	Green	1.5m (60")	9330783037469
MP-EX22G72	Green	1.8m (72")	9330783037476
MP-EX22G96	Green	2.4m (96")	9330783037483
MP-EX22G120	Green	3.0m (120")	9330783037490

	Mediplas	Expandable	Tubing ID (mm)	Tubing Colour - G=Green, B=Blue, N=Natural, P=Pink	Tubing Length (inches)
Legend:	MP	EX	15	N	120

Legend: (D) = Disposable ID = Inside Diameter (PP) = Polypropylene


Smooth-Bor Tubing

22mm Adult Cuff - Reusable



Part Number	Description	GTIN
SB403-6	150mm (6") Length	9330783021994
SB403-12	300mm (12") Length	9330783021901
SB403-16	400mm (16") Length	9330783021932
SB403-24	600mm (24") Length	9330783021949
SB403-36	900mm (36") Length	9330783021956
SB403-48	1.2m (48") Length	9330783021963
SB403-48BC	1.2m (48") Length	9330783021970
SB403-52	1.3m (52") Length	9330783021987
SB403-60	1.5m (60") Length	9330783022007
SB403-60BC	1.5m (60") Length	9330783022014
SB403-72	1.8m (72") Length	9330783022021
SB403-72-BC	1.8m (72") Length	9330783022038
SB403-84	2.1m (84") Length	9330783022045
SB403-84BC	2.1m (84") Length	9330783022052
SB403-96	2.4m (96") Length	9330783022069
SB403-96-BC	2.4m (96") Length	9330783022076
SB403-108	2.7m (108") Length	9330783021895
SB403-120	3.0m (120") Length	9330783021918
SB103-120-BC	3.0m (120") Length	9330783021925

15mm Child Cuff - Reusable



Part number	Description	GTIN
SB405-15-12	300mm (12") Length	9330783022113
SB405-15-12HE	300mm (12") Length	9330783022120
SB405-15-18	400mm (18") Length	9330783022137
SB405-15-24	600mm (24") Length	9330783022144
SB405-15-24HE	600mm (24") Length	9330783022151
SB405-15-28	700mm (28") Length	9330783022168
SB405-15-42HE	1.05m (42") Length	9330783034208
SB405-15-48	1.2m (48") Length	9330783022175
SB405-15-60	1.5m (60") Length	9330783022199
SB405-15-72	1.8m (72") Length	9330783022205

10mm Paediatric Cuff - Reusable



Part number	Description	GTIN
SB405-10-48	1.2m (48") Length	9330783022083
SB405-10-60	1.5m (60") Length	9330783022090
SB405-10-72	1.8m (72") Length	9330783022106

Legend: M = Male (External) F = Female (Internal) (D) = Disposable ID = Inside Diameter (S) = Silicone

Masks



Reusable Silicone Donut Masks (S)

PH10110 Silicone Donut Mask, Size #0, Small Infant, 15mm OD GTIN: 9330783029334

PH10111 Silicone Donut Mask, Size #1, Large Infant, 15mm OD GTIN: 9330783029341

PH10112 Silicone Donut Mask, Size #2, Small Child, 22mm ID GTIN: 9330783029358

Reusable and Autoclavable up to 134°C



Reusable Rendell Baker Transparent Silicone Masks (S)

1550 Silicone Mask, Size #0, Small Infant, 22mm ID GTIN: 9330783003907

1551 Silicone Mask, Size #1, Small Child, 22mm ID GTIN: 9330783003914

1552 Silicone Mask, Size #2, Medium Child, 22mm ID GTIN: 9330783003921

1553 Silicone Mask, Size #3, Large Child, 22mm ID GTIN: 9330783003938

Reusable and Autoclavable up to 134°C



Reusable Silicone Masks (S)

PH10113 Silicone Mask, Size #3, Child, 22mm ID GTIN: 9330783029365

PH10114 Silicone Mask, Size #4, Small Adult, 22mm ID GTIN: 9330783029372

PH10115 Silicone Mask, Size #5, Large Adult, 22mm ID GTIN: 9330783029389

Reusable and Autoclavable up to 134°C



Silicone Cushion Masks

PH10133 Silicone Mask, Size #3, Child, 22mm ID (S) (PC) GTIN: 9330783029402

PH10134 Silicone Mask, Size #4, Small Adult, 22mm ID (S) (PC) GTIN: 9330783029419

PH10135 Silicone Mask, Size #5, Large Adult, 22mm ID (S) (PC) GTIN: 9330783029426

Autoclavable up to 110°C



PVC Masks (PVC)

PH10143 PVC Mask, Size #3, Child, 22mm ID GTIN: 9330783029303

PH10144 PVC Mask, Size #4, Small Adult, 22mm ID GTIN: 9330783029310

PH10145 PVC Mask, Size #5, Large Adult, 22mm ID GTIN: 9330783029327

Autoclavable up to 80°C



Reusable Mask Hooks

PH10001 Mask Hook, Polysulfone (PS)

GTIN: 9330783034000

Reusable and Autoclavable up to 134°C

Legend: (PC) = Polycarbonate (PVC) = Polyvinyl Chloride (PS) = Polysulfone (S) = Silicone

Air Cushion Masks (with Inflate and Deflate Port)



Disposable Air Cushion Mask (PVC)(HDPE)

PH20116 Air Cushion Mask with Valve, Size #0, Neonate,
15mm OD, White Hook

GTIN: 9330783028412



Disposable Air Cushion Mask (PVC)(HDPE)

PH20117 Air Cushion Mask with Valve, Size #1, Infant,
15mm OD, Pink Hook

GTIN: 9330783028429



Disposable Air Cushion Mask (PVC)(HDPE)

PH20118 Air Cushion Mask with Valve, Size #2, Paediatric,
22mm ID, Yellow Hook

GTIN: 9330783028436



Disposable Air Cushion Mask (PVC)(HDPE)

PH20119 Air Cushion Mask with Valve, Size #3, Child,
22mm ID, Green Hook

GTIN: 9330783028443



Disposable Air Cushion Mask (PVC)(HDPE)

PH20120 Air Cushion Mask with Valve, Size #4, Small Adult,
22mm ID, Red Hook

GTIN: 9330783028450



Disposable Air Cushion Mask (PVC)(HDPE)

PH20121 Air Cushion Mask with Valve, Size #5, Large Adult,
22mm ID, Blue Hook

GTIN: 9330783028467

Legend: (PS) = Polysulfone (S) = Silicone (PVC) = Polyvinyl Chloride

RESPIRATORY CARE PRODUCTS

Scented Air Cushion Masks (with Inflate and Deflate Port)



Disposable Air Cushion Mask (PVC)(HDPE)

PH20140-B Blueberry Scented Air Cushion Mask with Valve, Size #0, Neonate, 15mm OD, White Hook GTIN: 9330783029150

PH20140-S Strawberry Scented Air Cushion Mask with Valve, Size #0, Neonate, 15mm OD, White Hook GTIN: 9330783029259

PH20140-V Vanilla Scented Air Cushion Mask with Valve, Size #0, Neonate, 15mm OD, White Hook GTIN: 9330783029204



Disposable Air Cushion Mask (PVC)(HDPE)

PH20141-B Blueberry Scented Air Cushion Mask with Valve, Size #1, Infant, 15mm OD, Pink Hook GTIN: 9330783029167

PH20141-S Strawberry Scented Air Cushion Mask with Valve, Size #1, Infant, 15mm OD, Pink Hook GTIN: 9330783029143

PH20141-V Vanilla Scented Air Cushion Mask with Valve, Size #1, Infant, 15mm OD, Pink Hook GTIN: 9330783029211



Disposable Air Cushion Mask (PVC)(HDPE)

PH20142-B Blueberry Scented Air Cushion Mask with Valve, Size #2, Paediatric, 22mm ID, Yellow Hook GTIN: 9330783029174

PH20142-S Strawberry Scented Air Cushion Mask with Valve, Size #2, Paediatric, 22mm ID, Yellow Hook GTIN: 9330783029266

PH20142-V Vanilla Scented Air Cushion Mask with Valve, Size #2, Paediatric, 22mm ID, Yellow Hook GTIN: 9330783029228



Disposable Air Cushion Mask (PVC)(HDPE)

PH20143-B Blueberry Scented Air Cushion Mask with Valve, Size #3, Child, 22mm ID, Green Hook GTIN: 9330783029181

PH20143-S Strawberry Scented Air Cushion Mask with Valve, Size #3, Child, 22mm ID, Green Hook GTIN: 9330783029273

PH20143-V Vanilla Scented Air Cushion Mask with Valve, Size #3, Child, 22mm ID, Green Hook GTIN: 9330783029235



Disposable Air Cushion Mask (PVC)(HDPE)

PH20144-B Blueberry Scented Air Cushion Mask with Valve, Size #4, Small Adult, 22mm ID, Red Hook GTIN: 9330783029280

PH20144-S Strawberry Scented Air Cushion Mask with Valve, Size #4, Small Adult, 22mm ID, Red Hook GTIN: 9330783029280

PH20144-V Vanilla Scented Air Cushion Mask with Valve, Size #4, Small Adult, 22mm ID, Red Hook GTIN: 9330783029242

Legend: (PVC) = Polyvinyl Chloride (HDPE) = High Density Polyethylene

Air Cushion Masks (MRI compatible with Inflate and Deflate Port)



Disposable Air Cushion Mask (PVC)(HDPE)(S)

PH20160 Air Cushion Mask with Valve, Size #0, Neonate, 15mm OD, White Hook and Non-metal valve - can be used during MRI process GTIN: 9330783036738



Disposable Air Cushion Mask (PVC)(HDPE)(S)

PH20161 Air Cushion Mask with Valve, Size #1, Infant, 15mm OD, Pink Hook and Non-metal valve - can be used during MRI process GTIN: 9330783036745



Disposable Air Cushion Mask (PVC)(HDPE)(S)

PH20162 Air Cushion Mask with Valve, Size #2, Paediatric, 22mm ID, Yellow Hook and Non-metal valve - can be used during MRI process GTIN: 9330783036752



Disposable Air Cushion Mask (PVC)(HDPE)(S)

PH20163 Air Cushion Mask with Valve, Size #3, Child, 22mm ID, Green Hook and Non-metal valve - can be used during MRI process GTIN: 9330783036769



Disposable Air Cushion Mask (PVC)(HDPE)(S)

PH20164 Air Cushion Mask with Valve, Size #4, Small Adult, 22mm ID, Red Hook and Non-metal valve - can be used during MRI process GTIN: 9330783036776



Disposable Air Cushion Mask (PVC)(HDPE)(S)

PH20165 Air Cushion Mask, Size #5, Large Adult, 22mm ID, Blue Hook and Non-metal valve - can be used during MRI process GTIN: 9330783036783



Disposable Air Cushion Mask (PVC)(HDPE)(S)

PH20166 Air Cushion Mask with Valve, Size #6, Extra Large Adult, 22mm ID, Orange Hook and Non-metal valve - can be used during MRI process GTIN: 9330783036790

Legend: (PVC) = Polyvinyl Chloride (HDPE) = High Density Polyethylene (S) = Silicone



Disposable Twin Port Air Cushion Mask - Child (PVC)(PC)

PH10283 Air Cushion Mask with Valve, Size: Child,
15mm OD Inspiratory / 22mm OD Expiratory

GTIN: 9330783036615



Disposable Twin Port Air Cushion Mask - Adult (PVC)(PC)

PH10285 Air Cushion Mask with Valve, Size: Adult,
15mm OD Inspiratory / 22mm OD Expiratory

GTIN: 9330783036622

Headgear (S)



PH20001 Silicone Head Strap, Length: 700mm
Autoclavable to 134°C

GTIN: 9330783028399



PH20003 Silicone Head Strap, Length: 440mm
Autoclavable to 134°C

GTIN: 9330783028405



PH60025
Teeth Opener
Colour: Blue (PC)

GTIN: 9330783036554

Legend: (PVC) = Polyvinyl Chloride (PC) = Polycarbonate (HDPE) = High Density Polyethylene

"Non-Vented" CPAP Full Face and Nasal Masks

'Comfort Seal, Soft & Flexible'



Reusable 'COMFO' Full Face CPAP Mask (PC) (S)

- Standard Head Gear
- Transparent Blue Colour Elbow - Non Vented

PH10217 'COMFO' CPAP Full Face Mask - Small GTIN: 9330783042296

PH10214 'COMFO' CPAP Full Face Mask - Medium GTIN: 9330783042302

PH10211 'COMFO' CPAP Full Face Mask - Large GTIN: 9330783042319



Reusable 'COMFO' Nasal CPAP Mask (PC) (S)

- Standard Head Gear
- Transparent Blue Colour Elbow - Non Vented

PH10208 'COMFO' CPAP Nasal Mask - Small GTIN: 9330783042326

PH10205 'COMFO' CPAP Nasal Mask - Medium GTIN: 9330783042333

PH10202 'COMFO' CPAP Nasal Mask - Large GTIN: 9330783042340



Reusable 'EASE' Full Face CPAP Mask (PC) (S)

- Standard Head Gear
- Transparent Blue Colour Elbow - Non Vented

PH10311 'EASE' CPAP Full Face Mask - Small GTIN: 9330783042357

PH10312 'EASE' CPAP Full Face Mask - Medium GTIN: 9330783042364

PH10313 'EASE' CPAP Full Face Mask - Large GTIN: 9330783042371



Reusable 'EASE' Nasal CPAP Mask (PC) (S)

- Standard Head Gear
- Transparent Blue Colour Elbow - Non Vented

PH10301 'EASE' CPAP Nasal Mask - Small GTIN: 9330783042388

PH10302 'EASE' CPAP Nasal Mask - Medium GTIN: 9330783042395

PH10303 'EASE' CPAP Nasal Mask - Large GTIN: 9330783042401

Legend: (PC) = Polycarbonate (S) = Silicone

RESPIRATORY CARE PRODUCTS

"Non-Vented" CPAP Full Face and Nasal Masks and Headgear

'Comfort Seal, Soft & Flexible'



Disposable '810 Type' Full Face CPAP Mask (PVC) (PC)

- Polypropylene Frame
- White Colour Head Gear
- Transparent Blue Colour Elbow - Non Vented
- Economic

PH11111 '810 Type' CPAP Full Face Mask - Small GTIN: 9330783042234

PH11112 '810 Type' CPAP Full Face Mask - Medium GTIN: 9330783042241

PH11113 '810 Type' CPAP Full Face Mask - Large GTIN: 9330783042258



Disposable '810 Type' Nasal CPAP Mask (PVC) (PC)

- Polypropylene Frame
- White Colour Head Gear
- Transparent Blue Colour Elbow - Non Vented
- Economic

PH11101 '810 Type' CPAP Nasal Mask - Small GTIN: 9330783042265

PH11102 '810 Type' CPAP Nasal Mask - Medium GTIN: 9330783042272

PH11103 '810 Type' CPAP Nasal Mask - Large GTIN: 9330783042289



PH10228



PH10229



PH10230

'EASE' & 'COMFO' Headgear

Suit all 'EASE' and 'COMFO' masks and available in 3 styles to suit your budget.

PH10228 Premium Quality Headgear (Blue) GTIN: 9330783042494

PH10229 Standard Quality Headgear (Black) GTIN: 9330783042500

PH10230 Silicone Headgear (S) GTIN: 9330783042517

Legend: (PVC) = Polyvinyl Chloride (PC) = Polycarbonate (S) = Silicone



Mediplas Respiratory Products Sdn Bhd
No 7 Jalan KIP 3, Taman Perindustrian KIP, 52200 Kuala Lumpur, Malaysia
Ph: +60 3 6272 7625 Fax: +60 3 6273 6452
E-mail: sales@mediplas.com.my Web: www.mediplas.com.my

Manual Resuscitator Accessories



Reusable PEEP Valves - Polysulfone

PH60014 PEEP Valve, 2-10cm H₂O, 22mm OD (PS)

GTIN: 9330783028726

PH60018 PEEP Valve, 5-20cm H₂O, 22mm OD (PS)

GTIN: 9330783028740



Reusable PEEP Valves - Polysulfone

PH60016 PEEP Valve, 2-10cm H₂O, 30mm ID (PS)

GTIN: 9330783028733

PH60020 PEEP Valve, 5-20cm H₂O, 30mm ID (PS)

GTIN: 9330783028757



Disposable PEEP Valves

PH60013 PEEP Valve, 2-10cm H₂O, 22mm OD (PC)

GTIN: 9330783037513

PH60017 PEEP Valve, 5-20cm H₂O, 22mm OD (PC)

GTIN: 9330783037537



Disposable PEEP Valves

PH60015 PEEP Valve, 2-10cm H₂O, 30mm ID (PC)

GTIN: 9330783037520

PH60019 PEEP Valve, 5-20cm H₂O, 30mm ID (PC)

GTIN: 9330783037544



PH60032

GTIN: 9330783028801

Reusable PEEP Valve Diverter (PS) (S)



PHSA60100046

GTIN: 9330783032549

Disposable PEEP Valve Diverter (PC)

Legend: OD = Outside Diameter ID = Inside Diameter (S) = Silicone (PC) = Polycarbonate (PS) = Polysulfone (PP) = Polypropylene

Oxygen Therapy



Oxygen Therapy Adult Mask

PH30130

Disposable Adult O₂ Therapy Mask with Elastic Strap (PE)(PVC)
Colour: Clear/Transparent

GTIN: 9330783028481



Oxygen Therapy Child Mask

PH30131

Disposable Child O₂ Therapy Mask with Elastic Strap (PE)(PVC)
Colour: Clear/Transparent

GTIN: 9330783028498



Oxygen Therapy Adult Mask

PH30132-2.1

Disposable Adult O₂ Therapy Mask with Elastic Strap and
2.1m Oxygen Tubing (PE)(PVC)
Mask Colour: Clear/Transparent

GTIN: 9330783028504



Oxygen Therapy Child Mask

PH30133-2.1

Disposable Child O₂ Therapy Mask with Elastic Strap and
2.1m Oxygen Tubing (PE)(PVC)
Mask Colour: Clear/Transparent

GTIN: 9330783028511



High Concentration Oxygen Therapy Adult Mask

PH30134-2.1

Disposable Adult O₂ Therapy Mask with Elastic Strap,
2.1m Oxygen Tubing and Reservoir Bag (PE)(PVC)
Mask Colour: Clear/Transparent

GTIN: 9330783028528



High Concentration Oxygen Therapy Child Mask

PH30135-2.1

Disposable Child O₂ Therapy Mask with Elastic Strap,
2.1m Oxygen Tubing and Reservoir Bag (PE)(PVC)
Mask Colour: Clear/Transparent

GTIN: 9330783028535

Legend: (PVC) = Polyvinyl Chloride (PE) = Polyethylene

Oxygen Therapy



Aerosol Therapy Adult Mask

PH30150

Disposable Adult Aerosol Mask with Elastic Strap (PVC)
Colour: Clear/Transparent

GTIN: 9330783028542



Aerosol Therapy Child Mask

PH30151

Disposable Child Aerosol Mask with Elastic Strap (PVC)
Colour: Clear/Transparent

GTIN: 9330783028559



Aerosol Therapy Adult Mask Set

PH30152-2.1

Disposable Adult Aerosol Mask with Elastic Strap, 2.1m Oxygen Tubing
and Jet Nebulizer (PVC)
Mask Colour: Clear/Transparent

GTIN: 99330783032297



Aerosol Therapy Child Mask Set

PH30153-2.1

Disposable Child Aerosol Mask with Elastic Strap, 2.1m Oxygen Tubing
and Jet Nebulizer (PVC)
Mask Colour: Clear/Transparent

GTIN: 9330783028566



PH30180



PH30181-2.1



PH30182



PH30183-2.1

Face Tent Mask

PH30180 Disposable Face Tent with Oxygen Adaptor
(PE)(PVC)

GTIN: 9330783028610

PH30181-2.1 Disposable Face Tent with Oxygen Adaptor and
2.1m Oxygen Tubing (PE)(PVC)

GTIN: 9330783028627

PH30182 Disposable Face Tent with Aerosol Adaptor
(PP)(PVC)

GTIN: 9330783028634

PH30183-2.1 Disposable Face Tent with Aerosol Adaptor,
Nebulizer Bottle and 2.1m Oxygen Tubing (PP)(PVC)

GTIN: 9330783033522

Legend: (PVC) = Polyvinyl Chloride (PE) = Polyethylene (PP) = Polypropylene



Nebulizer Bottle

PH30507 Reusable Nebulizer Bottle (PC)(PP)
18mm OD, 6cc Volume

GTIN: 9330783032358



Nebulizer Bottle

PH30508 Reusable Nebulizer Bottle (PC)(PP)
22mm ID / 18mm OD, 6cc Volume

GTIN: 9330783032372



Nebulizer Bottle

PH30505 Nebulizer Bottle (PP)
18mm OD, 20cc Volume

GTIN: 9330783033676

Legend: (PVC) = Polyvinyl Chloride (PE) = Polyethylene (PP) = Polypropylene

Oxygen Therapy



Tracheostomy Mask

PH30160 GTIN: 9330783028573
 Disposable Adult O₂ Tracheostomy Mask with Elastic Strap & 22mm Connector
 (PP)(PE)(PVC)

PH30161 GTIN: 9330783028580
 Disposable Child O₂ Tracheostomy Mask with Elastic Strap & 22mm Connector
 (PP)(PE)(PVC)



Tracheostomy Mask

PH30162 GTIN: 9330783028597
 Disposable Adult O₂ Tracheostomy Mask with Elastic Strap & 6mm Connector
 (PP)(PE)(PVC)

PH30163 GTIN: 9330783028603
 Disposable Child O₂ Tracheostomy Mask with Elastic Strap & 6mm Connector
 (PP)(PE)(PVC)



Tracheostomy Adult Mask Set

PH30164-2.1 GTIN: 9330783033577
 Disposable Adult O₂ Tracheostomy Mask with Elastic Strap and 2.1m
 Oxygen Tubing (PP)(PE)(PVC)



Tracheostomy Child Mask Set

PH30165-2.1 GTIN: 9330783033584
 Disposable Child O₂ Tracheostomy Mask with Elastic Strap and 2.1m
 Oxygen Tubing (PP)(PE)(PVC)



Tracheostomy Adult Mask Set with Concentration Jets

PH30166-2.1 GTIN: 9330783033591
 Disposable Adult O₂ Tracheostomy Mask with Elastic Strap, 6 Colour Coded
 Concentration Jets and 2.1m Oxygen Tubing (PP)(PE)(PVC)
 Diluters: Blue = 24%, Yellow = 28%, White = 31%, Green = 35%
 Red = 40%, Orange = 50%



Tracheostomy Child Mask Set with Concentration Jets

PH30167-2.1 GTIN: 9330783033607
 Disposable Child O₂ Tracheostomy Mask with Elastic Strap, 6 Colour Coded
 Concentration Jets and 2.1m Oxygen Tubing (PP)(PE)(PVC)
 Diluters: Blue = 24%, Yellow = 28%, White = 31%, Green = 35%
 Red = 40%, Orange = 50%

Legend: (PVC) = Polyvinyl Chloride (PE) = Polyethylene (PP) = Polypropylene



Adult Oxygen Mask with Concentration Jets

PH30154-2.1

GTIN: 9330783033539

Disposable Adult Oxygen Mask with Elastic Strap, 150mm Corugated Tubing, 2.1m Oxygen Tubing and 6 Diluters (PP)(PE)(PVC)

Diluters: Blue = 24%, Yellow = 28%, White = 31%, Green = 35%
Red = 40%, Orange = 50%

Mask Colour: Clear/Transparent



Child Oxygen Mask with Concentration Jets

PH30155-2.1

GTIN: 9330783033546

Disposable Child Oxygen Mask with Elastic Strap, 150mm Corugated Tubing, 2.1m Oxygen Tubing and 6 Diluters (PP)(PE)(PVC)

Diluters: Blue = 24%, Yellow = 28%, White = 31%, Green = 35%
Red = 40%, Orange = 50%

Mask Colour: Clear/Transparent



Venturi Diluter

PH30001 Diluter Blue - 24% / 4 LPM (PP)

GTIN: 9330783032044

PH30002 Diluter Yellow - 28% / 4 LPM (PP)

GTIN: 9330783032037

PH30003 Diluter White - 31% / 6 LPM (PP)

GTIN: 9330783032020

PH30004 Diluter Green - 35% / 8 LPM (PP)

GTIN: 9330783032013

PH30005 Diluter Red - 40% / 8 LPM (PP)

GTIN: 9330783032006

PH30006 Diluter Orange - 50% / 12 LPM (PP)

GTIN: 9330783031993

PH30007 Diluter Light Green - 60% / 12 LPM (PP)

GTIN: 9330783031986

PH30008 Venturi Hood 22mm ID (PP)

GTIN: 9330783031979



Child Venturi Mask

PH30156-2.1

GTIN: 9330783031962

Disposable Child Oxygen Mask with Elastic Strap, 2.1m Oxygen Tubing and Adjustable Diluter (PP)(PVC)(N)

Diluter: 28% - 60%

Mask Colour: Clear/Transparent



Adult Venturi Mask

PH30157-2.1

GTIN: 9330783031955

Disposable Adult Oxygen Mask with Elastic Strap, 2.1m Oxygen Tubing and Adjustable Diluter (PP)(PVC)(N)

Diluter: 28% - 60%

Mask Colour: Clear/Transparent



Adjustable Diluter

PH30020 Adjustable Diluter (PP)(N)

GTIN: 9330783033553

Diluter: 28% - 60%



Nebulizer Kit

PH30148-2.1

GTIN: 9330783033560

Disposable Nebulizer Kit including 150mm Corugated Tubing, Jet Nebulizer, 'T' Piece, Mouth piece and 2.1m Oxygen Tubing (PP)(PE)(PVC)(PS)

Legend: (PVC) = Polyvinyl Chloride (PE) = Polyethylene (PP) = Polypropylene (N) = Nylon (PS) = Polysulfone



RESPIRATORY CARE PRODUCTS

Laryngoscopes - Disposable

DUOSCOPE™

THE NEW DUAL BLADE DISPOSABLE EMERGENCY LARYNGOSCOPE



Available in both Adult/Child or Adult/Adult dual blade configurations, the DUOSCOPE™'s simplistic functional design facilitates immediate accessibility and instantaneous response in emergency situations. Produced in high quality plastic and equipped with dual illumination lamps (allowing over 24 hours continuous illumination), the disposable DUOSCOPE™ is a cost effective premium intubation for the emergency cart or kit.

- ✓ **Compact**
- ✓ **Easy to Use**
- ✓ **2 Blades in 1**
- ✓ **Light Weight**
- ✓ **Illumination Tester**
- ✓ **Batteries Easily Removed**
- ✓ **Long-Life Illumination**
(Over 24 hours continuous Illumination)
- ✓ **No Assembly - Quick Action**

RB2000	Miller Size #1 + MacIntosh Size #3 (Child)	GTIN: 9330783025664
RB3000	Miller Size #1 + MacIntosh Size #4 (Adult)	GTIN: 9330783025671
RB4000	MacIntosh Size #3 + MacIntosh Size #4 (Large Adult)	GTIN: 9330783025688

MULTIPLE APPLICATIONS:



Trauma Rescue



Ambulance



Medical Air Service



Emergency Carts



Military



RB2000

RB3000

RB4000



NOTE: Military Model available in "Desert Storm" Colour



Mediplas Respiratory Products Sdn Bhd
No 7 Jalan KIP 3, Taman Perindustrian KIP, 52200 Kuala Lumpur, Malaysia
Ph: +60 3 6272 7625 Fax: +60 3 6273 6452
E-mail: sales@mediplas.com.my Web: www.mediplas.com.my

Laryngoscopes - Disposable

ONESCOPE™

THE NEW LOW-COST DISPOSABLE EMERGENCY LARYNGOSCOPE



Handle can be used as a storage area for gloves etc...

- ✓ **Compact**
- ✓ **Easy to Use**
- ✓ **Light Tester**
- ✓ **Light Weight**
- ✓ **Glove Storage Area**
- ✓ **Battery is Easy to Remove**
- ✓ **Long-Life Illumination**
(Over 24 hours continuous illumination)
- ✓ **No Assembly Required**

Available in three popular blade configurations, the ONESCOPE™'s simplistic functional design facilitates immediate accessibility and instantaneous response in emergency situations. Produced in high quality plastic and equipped with a long-life illumination lamp (allowing over 24 hours continuous illumination), the disposable ONESCOPE™ is a cost effective premium intubation device for the emergency cart or kit.

GB1000	Miller Size #1	(Child)	GTIN: 9330783029433
GB3000	MacIntosh Size #3	(Adult)	GTIN: 9330783029440
GB4000	MacIntosh Size #4	(Large Adult)	GTIN: 9330783029457

MULTIPLE APPLICATIONS:



Trauma Rescue



Ambulance



Medical Air Service



Emergency Carts



Military



GB1000

GB4000

GB3000

Laryngoscopes - Disposable Fibreoptic / LED Lightweight Stainless Steel



MacIntosh - Fibreoptic / LED Blades

DLB/MAC/4	Large Adult Blade	Size: 4	GTIN: 9330783000357
DLB/MAC/3	Adult Blade	Size: 3	GTIN: 9330783000340
DLB/MAC/2	Paediatric Blade	Size: 2	GTIN: 9330783000333
DLB/MAC/1	Infant Blade	Size: 1	GTIN: 9330783000364
DLB/MAC/0	Neonate Blade	Size: 0	GTIN: 9330783031832



Miller - Fibreoptic / LED Blades

DLB/MIL/3	Adult Blade	Size: 3	GTIN: 9330783000401
DLB/MIL/2	Paediatric Blade	Size: 2	GTIN: 9330783000395
DLB/MIL/1	Infant Blade	Size: 1	GTIN: 9330783000388
DLB/MIL/0	Neonate Blade	Size: 0	GTIN: 9330783000371
DLB/MIL/00	Small Neonate Blade	Size: 00	GTIN: 9330783028153



Robertshaw - Fibreoptic / LED Blades

DLB/ROB/1	Infant Blade	Size: 1	GTIN: 9330783031856
DLB/ROB/0	Neonate Blade	Size: 0	GTIN: 9330783031849



Kessell - Fibreoptic / LED Blades

DLB/KES/4	Large Adult Blade	Size: 4	GTIN: 9330783030217
DLB/KES/3	Adult Blade	Size: 3	GTIN: 9330783030200



Seward - Fibreoptic / LED Blades

DLB/SEW/1	Infant Blade	Size: 1	GTIN: 9330783031863
-----------	--------------	---------	---------------------

RESPIRATORY CARE PRODUCTS

Laryngoscopes - Disposable LED Standard / Conventional

Lightweight Stainless Steel



Macintosh - LED Standard / Conventional Blades

MAC4/DSLB	Large Adult Blade	Size: 4	GTIN: 9330783033713
MAC3/DSLB	Adult Blade	Size: 3	GTIN: 9330783033706
MAC2/DSLB	Paediatric Blade	Size: 2	GTIN: 9330783033690
MAC1/DSLB	Infant Blade	Size: 1	GTIN: 9330783033683



Miller - LED Standard / Conventional Blades

MIL4/DSLB	Large Adult Blade	Size: 4	GTIN: 9330783033768
MIL3/DSLB	Adult Blade	Size: 3	GTIN: 9330783033751
MIL2/DSLB	Paediatric Blade	Size: 2	GTIN: 9330783033744
MIL1/DSLB	Infant Blade	Size: 1	GTIN: 9330783033737
MIL0/DSLB	Neonate Blade	Size: 0	GTIN: 9330783033720



Laryngoscope blades are available in a large variety of styles and sizes to suit every requirement.

Laryngoscopes - Reusable

MacIntosh Blades



STANDARD	WB02	Neonatal Blade	Size: 0	GTIN: 9330783023974
STANDARD	WB04	Infant Blade	Size: 1	GTIN: 9330783023981
STANDARD	WB06	Paediatric Blade	Size: 2	GTIN: 9330783024001
STANDARD	WB08	Adult Blade	Size: 3	GTIN: 9330783024018
STANDARD	WB10	Large Blade	Size: 4	GTIN: 9330783024025
STANDARD	WB05	Extra Large Blade	Size: 5	GTIN: 9330783023998
STANDARD	LHWB08	Left-hand Blade	Size: 3	GTIN: 9330783030002
STANDARD	HPWB08	Adult Blade (Hinge Pin Fitting)	Size: 3	GTIN: 9330783031047
STANDARD	WB39	Adult Blade (Polio)	Size: 3	GTIN: 9330783031030



FIBROPTIC	M702	Neonatal Blade	Size: 0	GTIN: 9330783017492
FIBROPTIC	M704	Infant Blade	Size: 1	GTIN: 9330783017508
FIBROPTIC	M706	Paediatric Blade	Size: 2	GTIN: 9330783017522
FIBROPTIC	M708	Adult Blade	Size: 3	GTIN: 9330783017539
FIBROPTIC	M710	Large Blade	Size: 4	GTIN: 9330783017553
FIBROPTIC	M705	Extra Large Blade	Size: 5	GTIN: 9330783017515
FIBROPTIC	LHM708	Left-hand Blade	Size: 3	GTIN: 9330783030224
FIBROPTIC	SM708	"Shiney" Blade	Size: 3	GTIN: 9330783030231



Fibre-Seal (German Style) - Macintosh Blades (enclosed fibreoptic light guide)				
FIBROPTIC	M756	Fibre-Seal - Infant Blade	Size: 1	GTIN: 9330783017621
FIBROPTIC	M757	Fibre-Seal - Paediatric Blade	Size: 2	GTIN: 9330783017638
FIBROPTIC	M758	Fibre-Seal - Adult Blade	Size: 3	GTIN: 9330783017645
FIBROPTIC	M759	Fibre-Seal - Large Adult Blade	Size: 4	GTIN: 9330783017652

Miller Blades



STANDARD	WB11	Neonatal Blade	Size: 00	GTIN: 9330783024032
STANDARD	WB12	Neonatal Blade	Size: 0	GTIN: 9330783024049
STANDARD	WB14	Infant Blade	Size: 1	GTIN: 9330783024056
STANDARD	WB16	Paediatric Blade	Size: 2	GTIN: 9330783024063
STANDARD	WB18	Adult Blade	Size: 3	GTIN: 9330783024070
STANDARD	WB20	Large Blade	Size: 4	GTIN: 9330783024087



FIBROPTIC	M724	Neonatal Blade	Size: 00	GTIN: 9330783017614
FIBROPTIC	M712	Neonatal Blade	Size: 0	GTIN: 9330783017560
FIBROPTIC	M714	Infant Blade	Size: 1	GTIN: 9330783017577
FIBROPTIC	M716	Paediatric Blade	Size: 2	GTIN: 9330783017584
FIBROPTIC	M718	Adult Blade	Size: 3	GTIN: 9330783017591
FIBROPTIC	M720	Large Blade	Size: 4	GTIN: 9330783017607



Fibre-Seal (German Style) - Miller Blades (enclosed fibreoptic light guide)				
FIBROPTIC	M767	Fibre-Seal - Infant Blade	Size: 00	GTIN: 9330783017690
FIBROPTIC	M769	Fibre-Seal - Paediatric Blade	Size: 0	GTIN: 9330783017706
FIBROPTIC	M760	Fibre-Seal (German Style)	Size: 1	GTIN: 9330783017669
FIBROPTIC	M762	Fibre-Seal (German Style)	Size: 2	GTIN: 9330783017676
FIBROPTIC	M763	Fibre-Seal (German Style)	Size: 3	GTIN: 9330783017683

Wisconsin Blades (Foregger)



STANDARD	WB22	Neonatal Blade	Size: 0	GTIN: 9330783024094
STANDARD	WB24	Neonatal Blade	Size: 1	GTIN: 9330783024100
STANDARD	WB26	Infant Blade	Size: 2	GTIN: 9330783024117
STANDARD	WB28	Paediatric Blade	Size: 3	GTIN: 9330783024131
STANDARD	WB30	Adult Blade	Size: 4	GTIN: 9330783024148
FIBROPTIC	W700	Neonatal Blade	Size: 0.5	GTIN: 9330783023882
FIBROPTIC	W702	Neonatal Blade	Size: 0	GTIN: 9330783023899
FIBROPTIC	W704	Infant Blade	Size: 1	GTIN: 9330783023905
FIBROPTIC	W706	Paediatric Blade	Size: 2	GTIN: 9330783023912
FIBROPTIC	W708	Adult Blade	Size: 3	GTIN: 9330783023936
FIBROPTIC	W710	Large Blade	Size: 4	GTIN: 9330783023950

Laryngoscopes - Reusable

		Seward Blades				
		STANDARD	WB47	Infant Blade	Size: 1	GTIN: 9330783024223
		STANDARD	WB49	Paediatric Blade	Size: 2	GTIN: 9330783024230
		FIBREOPTIC	S747	Infant Blade	Size: 1	GTIN: 9330783021826
		FIBREOPTIC	S749	Paediatric Blade	Size: 2	GTIN: 9330783021840
		<hr/>				
		Soper Blades				
		STANDARD	WB41	Infant Blade	Size: 1	GTIN: 9330783024162
		STANDARD	WB43	Paediatric Blade	Size: 2	GTIN: 9330783024179
		STANDARD	WB45	Adult Blade	Size: 3	GTIN: 9330783024186
 <p>Standard Blade Shown</p>		FIBREOPTIC	PS705	Infant Blade	Size: 1	GTIN: 9330783021048
		FIBREOPTIC	PS706	Paediatric Blade	Size: 2	GTIN: 9330783021055
		FIBREOPTIC	PS707	Adult Blade	Size: 3	GTIN: 9330783029853
<hr/>						
		Robertshaw Blades				
		STANDARD	WB53	Neonate Blade	Size: 0	GTIN: 9330783024247
		STANDARD	WB54	Infant Blade	Size: 1	GTIN: 9330783024254
		FIBREOPTIC	R703	Neonate Blade	Size: 0	GTIN: 9330783021345
		FIBREOPTIC	R704	Infant Blade	Size: 1	GTIN: 9330783021352
		<hr/>				
		Miller/Port Blades				
		STANDARD	WB65	Neonate Blade	Size: 1	GTIN: 9330783031818
		STANDARD	WB67	Infant Blade	Size: 2	GTIN: 9330783031825
		<hr/>				
		Flagg Blades				
		STANDARD	WB75	Infant Blade	Size: 1	GTIN: 9330783024261
		STANDARD	WB76	Paediatric Blade	Size: 2	GTIN: 9330783029907
		STANDARD	WB77	Adult Blade	Size: 3	GTIN: 9330783029914
 <p>Standard Blade Shown</p>		FIBREOPTIC	M775	Infant Blade	Size: 1	GTIN: 9330783029877
		FIBREOPTIC	M776	Paediatric Blade	Size: 2	GTIN: 9330783029860
		FIBREOPTIC	M777	Adult Blade	Size: 3	GTIN: 9330783029884
<hr/>						
		American Profile - Macintosh Blades				
		STANDARD	APWB02	Neonate Blade	Size: 1	GTIN: 9330783030873
		STANDARD	APWB04	Infant Blade	Size: 1	GTIN: 9330783030880
		STANDARD	APWB06	Paediatric Blade	Size: 2	GTIN: 9330783029921
		STANDARD	APWB08	Adult Blade	Size: 3	GTIN: 9330783029938
		STANDARD	APWB10	Large Adult Blade	Size: 4	GTIN: 9330783029945
 <p>Standard Blade Shown</p>		FIBREOPTIC	APM708	Adult Blade	Size: 3	GTIN: 9330783029891
		<hr/>				
		American Profile - Miller Blades				
		STANDARD	APWB12	Neonate Blade	Size: 1	GTIN: 9330783029952
		STANDARD	APWB14	Infant Blade	Size: 1	GTIN: 9330783029969
		STANDARD	APWB16	Paediatric Blade	Size: 2	GTIN: 9330783029976
		STANDARD	APWB18	Adult Blade	Size: 3	GTIN: 9330783029983
		STANDARD	APWB20	Large Adult Blade	Size: 4	GTIN: 9330783029990
<hr/>						
		Oxford Blade				
		STANDARD	WB27			GTIN: 9330783024124



RESPIRATORY CARE PRODUCTS

Laryngoscopes - Reusable and Laryngoscope Handles



Fibreoptic Blade Shown
(Removable Light Guide)

Sunflex Macintosh Blades with Flexible Tip

STANDARD	S0020	Adult Blade	Size: 3	GTIN: 9330783021581
STANDARD	S0021	Large Adult Blade	Size: 4	GTIN: 9330783021598
FIBROPTIC	G0072	Adult Blade	Size: 3	GTIN: 9330783013777
FIBROPTIC	G0073	Large Adult Blade	Size: 4	GTIN: 9330783013791



Fibreoptic Blade Shown
(Removable Light Guide)

Kessell Blades

STANDARD	K101	Adult Blade	Size: 3	GTIN: 9330783017256
STANDARD	K104	Large Adult Blade	Size: 4	GTIN: 9330783017263
FIBROPTIC	K202	Adult Blade	Size: 3	GTIN: 9330783017270
FIBROPTIC	K204	Large Adult Blade	Size: 4	GTIN: 9330783017294

Laryngoscope Handles



Standard / Conventional Laryngoscope Handles - Reusable

B906	Small	L:150mm	D: 20mm	GTIN: 9330783009749
B904	Medium	L:150mm	D: 30mm	GTIN: 9330783009725
B808	Stubby	L:125mm	D: 32mm	GTIN: 9330783009718



Fibreoptic Laryngoscope Handles - Reusable

B706	Small	L:150mm	D: 20mm	GTIN: 9330783009701
B704	Medium	L:150mm	D: 30mm	GTIN: 9330783009688
B608	Stubby	L:125mm	D: 32mm	GTIN: 9330783009671
GB704	Medium (for Fibre-Seal Style Blades (German Style))			GTIN: 9330783013951
GB706	Small (for Fibre-Seal Style Blades (German Style))			GTIN: 9330783013968



Assembled
Laryngoscope

LED Laryngoscope Handles - Reusable

B704/LED	Medium	L:150mm	D: 30mm	GTIN: 9330783031276
B706/LED	Small	L:150mm	D: 20mm	GTIN: 9330783031283
GB704/LED	Medium (for Fibre-Seal Style Blades (German Style))			GTIN: 9330783033775
GB706/LED	Small (for Fibre-Seal Style Blades (German Style))			GTIN: 9330783033782

Laryngoscopes Accessories



Laryngoscope Bulbs

STANDARD	B951	Large - Frosted	GTIN: 9330783009756
(Suit: All Parker Reusable Standard Blades Not Listed Below)			
STANDARD	B952	Small - Frosted	GTIN: 9330783009763
(Suit: WB02, WB11, WB12, WB14 Parker Reusable Blades)			
FIBREOPTIC	B968	Halogen, Fine Thread (Silver)	GTIN: 9330783009794
(Suit: B704, B706, B608, GB704, GB706 Parker Reusable Fiberoptic Handles)			
LED	PHAG-0047	LED Laryngoscope Bulb	GTIN: 9330783033799

Standard Spare Parts

S901	Contact Port Assembly	GTIN: 9330783021864
S902	Wiring Assembly Complete	GTIN: 9330783021871
S907	Small Handle Housing	GTIN: 9330783033799
S905	Medium Handle Housing	GTIN: 9330783033812
S908	Small Handle Head Assembly	GTIN: 9330783036707
S909	Medium Handle Head Assembly	GTIN: 9330783036714

Fiberoptic Spare Parts

S7021/4	F'Optic MAC Assembly - Size 1 (4mm Bundle)	GTIN: 9330783021741
S7022/4	F'Optic MAC Assembly - Size 2 (4mm Bundle)	GTIN: 9330783021765
S7023/4	F'Optic MAC Assembly - Size 3 (4mm Bundle)	GTIN: 9330783021789
S7024/4	F'Optic MAC Assembly - Size 4 (4mm Bundle)	GTIN: 9330783021802
O236	Switch Assembly	GTIN: 9330783019168
S703	Handle Housing	GTIN: 9330783034932
B960	Metal Cap for Bulb	GTIN: 9330783009787

Laryngoscope Carry Case (Suits Fiberoptic & Standard Sets)

P902	Plastic Case (1 x Handle & 4 x Blades)	GTIN: 9330783019557
------	---	---------------------



Otoscope

O766	Otoscope Set with Pouch	GTIN: 9330783031900
------	--------------------------------	---------------------



Magill Forceps Stainless Steel

M043	Infant (18cm) Satin Finish Forceps	GTIN: 9330783017423
M044	Child (20cm) Satin Finish Forceps	GTIN: 9330783017430
M045	Adult (23cm) Satin Finish Forceps	GTIN: 9330783017447

Accessory Products



Ideal for sealing containers.
Contains writable contents tag

SECUR-A-TAG Security Ties

- Tamper-Evident Design
- Assortment of Colours
- Easy to Use - Loop, Insert & Pull
- Write-on Identification Tag
- Secures up to 45mm diameter
- Hundreds of Uses

73287-50/A Security Ties - 10 of each colour (Box of 50) GTIN: 9330783027965

73287-50/B Security Ties - BLUE only (Box of 50) GTIN: 9330783027996


73287-50/G Security Ties - GREEN only (Box of 50) GTIN: 9330783027972

73287-50/P Security Ties - PINK only (Box of 50) GTIN: 9330783027941

73287-50/R Security Ties - RED only (Box of 50) GTIN: 9330783027958

73287-50/Y Security Ties - YELLOW only (Box of 50) GTIN: 9330783027989

NOTE: Colours may vary slightly from pictured items



SECUR-A-TAG
The essential travelling companion

Biegler Blood & Infusion Warmer

BIEGLER BW685 S is a high performance, microprocessor-controlled device for warming infusions and transfusions to prevent intra- and postoperative hypothermia.

Temperature can be set freely in increments of 0.5°C from 37°C to 41°C.

The BW685 S encompasses a slim and ergonomic design. The integrated universal mounting clamp allows the device to be placed on an infusion stand as well as on a normal rail.



BW685 S (with TubeFlow connection port)

GTIN: 9330783024919

- Free Temperature Pre-set 37°~ 41°C
- **No Dedicated Consumables - Lowering Costs**
- Microprocessor Controlled
- Optimised Warming Circuit
- High Transfusion Rates
- Minimised Energy Consumption
- Low/High Temperature Alarm
- Ready to Use in 40 seconds
- Patented Designed Heating Surface
- Multiple Extension Set Use at One Time
- Dry Heat System (NO HAZARDOUS WATER)
- Mounted on normal railing and an infusion stand
- Integrated handle
- Surface Temperature Monitor
- No Moving Parts (reduced maintenance)
- Permanent running self-tests
- Easy to Clean
- Slim, ergonomic design

NOTE: Tubeflow Assembly is Optional



Biegler TubeFlow Assembly

LG400TF01 Tubeflow Assembly

GTIN: 9330783024926

What is TubeFlow?

TubeFlow is an actively heated, reusable silicone profile with control unit for temperature optimisation between warmer and patient at low flow rates.

Thermal Grade Extension Sets

Thermal Grade Extension Sets Suitable for BW385, BW485 and BW585 Blood and Infusion Warmers


- To suit dry heat BW385, BW485 & BW585 blood and infusion warmers
- Single Use Only - Sterilised with ethylene oxide
- Luer lock connection



Standard Extension Sets

WB35000 Length = 3.5 m, 3.0mmID/4.1mmOD GTIN: 9330783024155

WB45000 Length = 4.5 m, 3.0mmID/4.1mmOD GTIN: 9330783024193



Needless Port (NP) Extension Sets

WB45000(NP) Length = 4.5 m with needless port, 3.0mmID/4.1mmOD GTIN: 9330783024209

Thermal Grade Extension Sets Suitable for BW685 and BW685 S Blood and Infusion Warmers

- To suit dry heat **BW685 & BW685 S** blood and infusion warmer
- Single Use Only - Sterilised with ethylene oxide
- Luer lock connection



Standard Extension Sets

WB45000 Length = 4.5 m, 3.0mmID/4.1mmOD GTIN: 9330783024193

WB56000 Length = 5.6 m, 3.0mmID/4.1mmOD GTIN: 9330783026005



Needless Port (NP) Extension Sets

WB45000(NP) Length = 4.5 m with needless port, 3.0mmID/4.1mmOD GTIN: 9330783024209

WB56000(NP) Length = 5.6 m with needless port, 3.0mmID/4.1mmOD GTIN: 9330783025992

Blood & Fluid Infusion Warmers

Protherm II Blood and Infusion Warmer



Protherm II

A universal device for a variety of applications, depending on media.

Suitable for:-

- Warming infusions
- Normal and massive infusions
- Heating liquids without risk of contamination
- Simple, fast and user-friendly operation
- Liquid heated gently with three independently controlled heat zones
- Permanently running self-tests and error display
- Easy to service with serial interface
- Ergonomic design makes it easily transportable and mountable on all commercially available infusion stands
- Economical material
- Low temperature alarm
- Automatic adaption of the heating control system to various media

LB1000000 Protherm II Blood and Infusion Warmer

GTIN: 9330783017393



Heating Bag, complete with Bubble Trap

Disposables:

FP1003001 Heating bag large, complete with bubble trap
(Box of 20)

GTIN: 9330783013593

Blood & Fluid Infusion Warmers and Extension Sets



LG400001

Biegler "Autopress" is an important part of our warming philosophy and works directly with either our blood and infusion warmer line to deliver blood and fluids at a high flow rate or works as a pressure stand alone unit to deliver a constant pressure up to 300 mm/Hg wherever needed.

LG400001 Biegler Autopress

GTIN: 9330783017416

- Preset pressure is always maintained
- Precise pressure setting due to sure setting electronic control
- Accomodates 2 x 500ccm, 2 x 1000ccm bags or 1 x 3000 ccm bag
- Fits on any IV pole
- Pressure limited to 300 mm/Hg
- Mains or Battery operated
- Cuffs can be emptied rapidly and easily
- Setup time significantly reduced
- Battery back-up capacity ~ 2h / operation

Spare Bags:

JR0000500 - 500ccm capacity

GTIN: 9330783017201

JR0001000 - 1000ccm capacity

GTIN: 9330783017218

JR0003000 - 3000ccm capacity

GTIN: 9330783017225



LB1000004 (ESH04) Emergency Transport Warmer

Emergency Transport Warmer

LB1000004 (ESH04) Emergency Transport Warmer

GTIN: 9330783017409

The **BIEGLER ESH04** Emergency Transport Warmer is designed to bring one or two flasks of blood or infusion liquids up to body temperature and hold them there (without a mains connection) for up to 90 minutes at the site of use.

The integrated pressure cuff eliminates the need for further pressurised infusion accessories in the emergency medical kit.

- Designed for mobile use
- Versatile, compact size, low weight
- Fits into emergency kit
- Maintains the temperature for up to 90 minutes at the site of use
- 24h continuous operation possible
- Integrated pressure cuff
- Solution dispensed directly from warmer



515800

515820

515824

High / Low Flowmeters

515800	Oxygen Flowmeter with Handwheel 15 LPM	GTIN: 9310574006294
515820	Air Flowmeter with Handwheel 15 LPM	GTIN: 9310574006300
515824	Oxygen Flowmeter with Handwheel 2.5 LPM	GTIN: 9310574006348
551065	Oxygen Flowmeter 0-60 LPM	GTIN: 9310574014985
551064	Air Flowmeter 0-60 LPM	GTIN: 9310574014978

Spare Parts

515830	Nipple/Wingnut Assembly, Oxygen, Ezi-Flow (10/pack)	GTIN: 9330783006311
515837	Service Kit, Ezi Flow	GTIN: 9330783006359
515849	Nipple/Wingnut Assembly, Air, Ezi-Flow (10/pack)	GTIN: 9330783006366
515856	Spare Needle Valve Assembly, Ezi-Flow	GTIN: 9330783006380



4MFA1001AU

6MFA1001AU

Low Flowmeters

6MFA1001AU	Low Flow Oxygen Flowmeter 25-200cc/min	GTIN: 9330783007387
4MFA1001AU	Low Flow Oxygen Flowmeter 0.1-1 LPM	GTIN: 9330783006090

Spare Parts

503213	Flow Tube For 6MFA1001AU	GTIN: 9330783006144
503214	Flow Tube for 4MFA1001AU	GTIN: 9330783006151
0122	White Flowmeter Nipple For 4MF and 6MF	GTIN: 9330783032600



Perflow High-precision Flowmeter

High-precision flowmeter with pre-calibrated flows. The "click-style" positioning is accurate and prevents accidental flow changes. Pressure stabilised, micro regulator ensures flow accuracy in the event of changes to supply pressure.

- Compact & Robust
- Works at ANY angle
- Models available to suit a wide range of applications from Neonatal to Adult
- No more expensive glass tube repairs
- No more "Guessing" of flow rates
- Allows fine titration of flows

4540840	0-1 lpm Precision Oxygen Flowmeter	GTIN: 9330783031443
4540850	0-5 lpm Precision Oxygen Flowmeter	GTIN: 9330783031450
4540860	0-15 lpm Precision Oxygen Flowmeter	GTIN: 9330783031467




Gauge Type Flowmeter

TM17	Gauge Type Flowmeter	GTIN: 9310574006393
------	----------------------	---------------------

Spare Parts


515830	Outlet Nut / Nipple	GTIN: 9330783006311
551041	Spare Flow Gauge	GTIN: 9310574015821

Oxygen Therapy



518804

518809



518808

518810

O-Series Regulators

518804	O-Series Oxygen Regulator with Yoke	GTIN: 9310574007697
518808	O-Series Nitrous Oxide Regulator with Yoke	GTIN: 9310574007710
518809	O-Series Air Regulator with Yoke	GTIN: 9310574007727
518810	Entonox Regulators with Yoke	GTIN: 9310574007734

Spare Parts

519026	Kit, Service, Seat, Series-O	GTIN: 9330783006465
551022	Spare, Bodok Seals, Regulators (10/pack)	GTIN: 9310574015593
551045	Spare, Gauge, Oxygen, 30000kPa, Series-O	GTIN: 9310574015869

Multe-Flow™ Piston Regulators

Oxygen Therapy Model

551155	Piston Regulator with 1, 1.5, 2, 3, 4, 5, 6, 8, 10 and 15 lpm flow control	GTIN: 9310574024564
--------	--	---------------------

The Multe-Flow™ piston style regulator offers superior safety and performance in a compact and lightweight design.

- All brass* design, chrome plated for easy cleaning
 - Complies to safety and performance requirements AS3840.1-1998
 - Suitable for use with Oxygen cylinders filled to pressures up to 200 Bar
 - Oxygen Therapy and Resuscitation models
 - A variety of flow options available
 - Swivel yoke allows preferred positioning when assembled to an Oxygen cylinder
 - Built in pressure relief valve prevents over pressure conditions to occur
 - Easy to read cylinder contents gauge, brass construction and low profile resists damage
 - 2 Self Sealing SIS outlet connections on the resuscitation models allow equipment to be attached or removed without the need to turn off the oxygen
 - For convenience, comes with a permanently attached cylinder key
 - Contoured flow selector with easy to read flow settings
- * Non Brass Yoke, all brass yoke option available on request

Resuscitation Model

551157	Piston Regulator with two (2) SIS outlets and 2, 3, 4, 5, 6, 8, 9, 10, 15 and 25 lpm flow control	GTIN: 9330783006779
551158	Piston Regulator with two (2) SIS outlets and 8 lpm and 15 lpm flow control	GTIN: 9310574024571

Spare Parts

551162	Spare Cylinder Key Spanner	GTIN: 9330783006786
551166	Spare Yoke Wingscrew	GTIN: 9330783006793
551169	Spare Oxygen Gauge	GTIN: 9330783025769
551174	Spare Therapy Outlet Nipple	GTIN: 9330783006809
551175	Spare SIS Oxygen Outlet	GTIN: 9330783036806



Suction Attachments and Regulators



Suction Attachments - Central Vacuum (with On-Off Tap)

- 554006 Suction Attachment - with Gauge (0 to -100kpa) GTIN: 9310574010314
554007 Suction Attachment - without Gauge GTIN: 9310574010307

Spare Parts

- SF2000 Disposable Suction Filter GTIN: 9330783022380
SFH3000 Spare Disposable Outlet Nozzle / Filter Holder GTIN: 9330783014582
SFK3000 Disposable Suction Filter Kit GTIN: 9330783024568
522501 Spare Vacuum Gauge (0 to -100kpa) GTIN: 9330783034963
7964024 Spare O'ring for Filter Kit GTIN: 9330783007448



Suction Attachments - Venturi Suction (without On-Off Tap)

High Suction (0 to -100kpa)

- 554013 Suction Attachment - with Gauge GTIN: 9330783006861
554010 Suction Attachment - without Gauge GTIN: 9330783006854

Low Suction (0 to -25kpa)

- 518535 Suction Attachment - with Gauge GTIN: 9330783008656

Spare Parts

- SF2000 Disposable Suction Filter GTIN: 9330783022380
SFH3000 Spare Disposable Outlet Nozzle / Filter Holder GTIN: 9330783014582
SFK3000 Disposable Suction Filter Kit GTIN: 9330783024568
522501 Spare Vacuum Gauge (0 to -100kpa) GTIN: 9330783006618
518795 Spare Vacuum Gauge (0 to -25kpa) GTIN: 9330783006458
7964024 Spare O'ring for Filter Kit GTIN: 9330783007448



Suction Regulators

- V9000 High Suction Vacuum Regulator
c/w Gauge and Handwheel
0 to -750 mmHg (0 to -100 kPa) GTIN: 9330783023806

- V9200 Low Suction Vacuum Regulator
c/w Gauge and Handwheel
0 to -180 mmHg (0 to -25 kPa) GTIN: 9330783023813

V9000-V9200 Spare Parts

- 3330004 Filter For Vacuum Regulators * V9000/V9200 GTIN: 9330783005406
3330101 Spare Vacuum Gauge 0-100kPa, For V9000 GTIN: 9330783041947
3330102 Spare Vacuum Gauge 0-25kPa, For V9200 GTIN: 9330783005437
4330001 Maintenance Kit For V9000 Vacuum Reg GTIN: 9330783005970
4330150 Watertrap with Security Valve & Particle Filter GTIN: 9330783006007
4330201 Maintenance Kit For V9200 Vacuum Reg GTIN: 9330783006014



Suction Regulators

- 5351000 High Suction Vacuum Regulator
c/w Gauge, Fast ON-OFF Switch and Handwheel
0 to -750 mmHg (0 to -100 kPa) GTIN: 9330783043972

- 5352000 Low Suction Vacuum Regulator
c/w Gauge, Fast ON-OFF Switch and Handwheel
0 to -150 mmHg (0 to -20 kPa) GTIN: 9330783043989

Suction Filter Kits



SF2000 Suction Filter

The SF2000 is a disposable suction filter replacing the common sintered bronze style currently used in all Comweld and Clements suction attachments. Its ability to change colour when in contact with fluid ensures a clear visual indication of a contaminated filter.

- REPLACES WA185 FILTER
- Suits Comweld (C.I.G.) and Clements suction attachments
- Visual indication of contamination
- Economical replacement · Easy to change · Time-saving
- Disposable · NO CLEANING

SF2000 Suction Filter (Pkt. 20 minimum) (HDPE)

GTIN: 9330783022380

FILTER TURNS BLUE WHEN CONTAMINATED



Suction Filter Housing (Disposable)

- REPLACES ORIGINAL HOUSING
- Available for Comweld & Clements suction attachments
- Made from Easy-view clear material
- Economical replacement

Comweld (C.I.G.) Housing

SFH3000 Housing (Pkt. 20) (PSY)

GTIN: 9330783014582

Clements Housing

SFH5000 Housing (Pkt. 20) (PSY)

GTIN: 9330783024599

SFH8000 Suction Tubing Joiner (Pkt. 10) (PSY)

GTIN: 9330783036639



O'ring

7964024 O'ring (N)

GTIN: 9330783007448



SFK3000 Suction Filter Kit

Comweld Suction Filter Kits

Comweld (C.I.G.) Suction Filter Kit

Comprises of: Filter, Housing & O'ring

SFK3000 Comweld Suction Filter Kit (Pkt. 20)
(PSY)(HDPE)(N)

GTIN: 9330783024568



SFK5000 Suction Filter Kit

Clements Suction Filter Kits

Clements Suction Filter Kit

Comprises of: Filter, Housing & O'ring

SFK5000 Clements Suction Filter Kit (Pkt. 20)
(PSY)(HDPE)(N)

GTIN: 9330783024575

NOTE: Above suction products are suitable for use with:

COMWELD - 554006, 554007, 554013, 554010, 518535, 554067, 554068

CLEMENTS - SUC891XX* Series, SUC892XX Series, SUC895XX Series, SUC896XX Series.

Suction



Twin-O-Vac Multi-purpose Suction/Oxygen Device

TM117G Twin-O-Vac Multi-purpose Device (High Vacuum) GTIN: 9310574006270

518632 Twin-O-Vac Multi-purpose Device (Low Vacuum) GTIN: 9310574006522

Spare Parts:

554050 Bacterial Filter (Pkt 100)

GTIN: 9310574019850

522501 Vacuum Gauge, 0 to -100kpa

GTIN: 9330783034963

515754 Service Kit

GTIN: 9330783006267

FOR A COST-EFFECTIVE DISPOSABLE JAR ALTERNATIVE, REFER BELOW



Central Vacuum Systems

518607 Suction Trap Jar with On/Off Valve

GTIN: 9330783006427

FOR A COST-EFFECTIVE DISPOSABLE JAR ALTERNATIVE, REFER BELOW



TV400-03

Assembled on a
Twin-O-Vac Unit

Disposable Twin-O-Vac Suction Jar

Save Time, Money and Resources with the new DISPOSABLE Twin-O-Vac JAR

- No more messy and time-consuming reprocessing.
- Reduce the risk of splash injuries.
- Reduce the risk of cross-infection.
- Now with NEW handy Jar lid to seal Jar contents after removal.

The new DISPOSABLE Twin-O-Vac Jar is the perfect solution for your suctioning needs.

TV400-01 Suction Jar Only GTIN: 9330783032082

TV400-02 Lid Only GTIN: 9330783032075

TV400-03 Suction Jar and Lid Assembly (with SECUR-A-TAG) GTIN: 9330783032068

73287-50/Y SECUR-A-TAG Security Tie (Pkt. 50) (Yellow) GTIN: 9330783027989

Specifications:

Capacity: 400ml

Length: 180mm

Embossed marks every 100ml

Materials: PP



TV400-02



73287-50/Y

(PP) = Polypropylene

Entonox Administration / Spare Parts

3 Way Adaptors



'Y' Body



'T' Body

'Y' Body

551140/Y 3 Way Oxygen Adaptor (White)

GTIN: 9330783006755

551141/Y 3 Way Air Adaptor (Black)

GTIN: 9330783006762

'T' Body

551140 3 Way Oxygen Adaptor (White)

GTIN: 9301574023871

551141 3 Way Air Adaptor (Black)

GTIN: 9301574023888

551142 3 Way Nitrous Oxide Adaptor (Blue)

GTIN: 9301574023895

551139/S 3 Way Suction Adaptor (Yellow)

GTIN: 9330783006748



'Y' Body



'T' Body

3 Way Adaptors with On/Off Valve on Inlet (Handwheel) Side

'Y' Body

G0188-S-OXY 3 Way Oxygen Adaptor c/w On/Off Valve

GTIN: 9330783039708

G0188-S-AIR 3 Way Air Adaptor c/w On/Off Valve

GTIN: 9330783039692

'T' Body

G0187-S-OXY 3 Way Oxygen Adaptor c/w On/Off Valve

GTIN: 9330783039678

G0187-S-AIR 3 Way Air Adaptor c/w On/Off Valve

GTIN: 9330783039685

Entonox Analgesic Administration System

The system is a rugged demand inhalation unit for the administration of 50% Nitrous Oxide / 50% Oxygen analgesic gas mixture (Entonox).

The system is ideal for the relief of pain in a wide range of situations.

- Transport of accident victims
- Casualty areas
- Minor surgical procedures
- Coronary care units with standard ventilators
- Obstetric analgesia



550010 Entonox Analgesia Administration System

GTIN: 9310574011571

Consists of: Series-O-Entonox Regulator, Bleed Pin, 2m Hose Assembly, Demand Valve & Mask

Humidifier



TM11 Humidifier

GTIN: 9310574006171

Spare Parts

551054 Spare Jar

GTIN: 9310574011731

551057 Service Kit - Diffuser, washer, O'ring

GTIN: 9330783006731

Oxygen Therapy Accessories



Handwheel Kits

518405	AIR Handwheel Kit	GTIN: 9301574008526
518406	OXYGEN Handwheel Kit	GTIN: 9301574025790
518407	NITROUS OXIDE Handwheel Kit	GTIN: 9301574025806

Tubing

PH100/OXY	Oxygen Tubing (per metre)	GTIN: 9330783019663
PH100/AIR	Air Tubing (per metre)	GTIN: 9330783019649
PH100/NIT	Nitrous Oxide Tubing (per metre)	GTIN: 9330783019656



Hose Assembly

4m. length Fitted with Colour-coded non-interchangeable handwheels.

MA164	Oxygen Hose Assembly	GTIN: 9330783017713
MA165	Air Hose Assembly	GTIN: 9330783017751
MA166	Nitrous Oxide Hose Assembly	GTIN: 9330783017768

Other lengths available upon request



Gauges

551042	Contents Gauge, Plain, 0-30000kpa (Series-0 Regulator)	GTIN: 9301574015838
551043	Contents Gauge, Nitrous Oxide, 0-20000kpa (Series-0 Regulator)	GTIN: 9301574015845
551044	Contents Gauge, Air, 0-30000kpa (Series-0 Regulator)	GTIN: 9301574015852
551045	Contents Gauge, Oxygen, 0-30000kpa (Series-0 Regulator)	GTIN: 9301574015869
522501	Flow Gauge, Vacuum, 0 to -100kpa (Twin-0-Vac/Suction Attachments)	GTIN: 9330783019649
551041	Flow Gauge, Oxygen, 0-14lpm (TM17 Flowmeter)	GTIN: 9301574015821



Cylinder Key, Wingscrew and Spanner

511956	Cylinder Keywheel	GTIN: 9330783006182
551162	Cylinder Spanner	GTIN: 9330783006786
302920	Black Wingscrew for Medical Yoke Regulator	GTIN: 9330783043194

A keywheel or spanner design product to turn medical gas cylinders on and off



Flowmeter Nipple and Wingnut

515830	White Nipple & Wingnut for Oxygen Flowmeters (Pkt. 10)	GTIN: 9330783006311
515849	Black Nipple & Wingnut for Air Flowmeters (Pkt. 10)	GTIN: 9330783006366



Bodok Seals

551022	Bodok Seal (Pkt. 10)	GTIN: 9310574015593
--------	----------------------	---------------------

Suits Series-0 regulators or any equipment with pin indexed yoke inlet fittings for medical gas cylinders.



Twin-0-Vac Filter Spares

554050	Bacterial Filters (Pkt. 100)	GTIN: 9310574019850
--------	------------------------------	---------------------

Suits all Twin-0-Vac units.

Infrared Thermometers



Infrared No Contact Thermometer

- Comfortable and Easy To Use
- No Probe Filter
- Accurate and Reliable
- No Ready Time

FS100 No Contact Thermometer

GTIN: 9330783039388



Infrared No Contact Thermometer

- No Contact, No Virus
- Comfortable and Easy to Use
- No Probe Filter
- Multi-measurement (Body Temperature/Object Temperature)
- Fever Alarm
- Multi-colour Back Light

FS700 No Contact Thermometer

GTIN: 9330783041992



Infrared Ear Thermometer

- Multiple Use
- Fast Measurement
- 8 Memories
- Auto-Off

TB100 No Contact Thermometer

GTIN: 9330783000128

PF-20 Probe Filter
Suit TB100 Thermometer (Pkt. 20)

GTIN: 9330783019571

Blood Pressure & Glucose Monitors



Compact Blood Pressure Monitor

- Detects and Displays Arrhythmia Symptoms
- 60 Reading Memories with Time and Date
- Easy to Read Digital Display
- Compact Size, Light, portable for Travel

BP100 Blood Pressure Monitor

GTIN: 9330783039401



Desktop Blood Pressure Monitor

- Detects and Displays Arrhythmia Symptoms
- 60 Reading Memories with Time and Date
- Displays 3 Levels of Blood Pressure State
- Simple One Button Operation

BP400 Blood Pressure Monitor

GTIN: 9330783031948



Desktop Blood Pressure Monitor with Cuff Storage

- Detects and Displays Arrhythmia Symptoms
- 60 Reading Memories with Time and Date
- Displays 3 Levels of Blood Pressure State
- Simple One Button Operation
- Cuff Storage Combined

NBP100 Blood Pressure Monitor

GTIN: 9330783039418



Blood Pressure and Glucose Monitor

Blood Pressure

- Detects and Displays Arrhythmia Symptoms
- 180 Memories with Time and Date
- Advanced Accuracy via Fuzzy Logic Technology

Blood Glucose

- Fast and Accurate Results: Just 5 seconds
- Tiny Blood Sample: Just 1 microlitre
- Large memory Capacity with Averaging: 270 test results
- Display Average Value during 14 days / 30 days

HMF100 Blood Pressure and Glucose Monitor

GTIN: 9330783031870

TS-2050 Glucose Monitor Strip (Pkt. 50)

GTIN: 9330783031894



Glucose Monitor

- Fast and Accurate Results: Just 5 seconds
- Tiny Blood Sample: Just 1 microlitre
- Large memory Capacity with Averaging: 250 test results
- Easy To Read LCD and Visual Icons
- User Friendly Functions

GM100 Blood Pressure Monitor

GTIN: 9330783031887

TS-2050 Glucose Monitor Strip (Pkt. 50)

GTIN: 9330783031894

Handheld & Finger Style Pulse Oximeters



Handheld SpO₂/HR Monitor

- Highly Accurate
- Quick SpO₂ and pulse rate reading
- High resolution backlit LCD screen
- Adjustable alarm limits for SpO₂ and PR signals
- Visual and audio alarm signals
- Compact design
- Low power consumption - up to 24 hours of continual use
- Trend view - display and view past readings
- 24 hours record storage

PH-MD300-A Handheld SpO₂ Pulse Oximeter

GTIN: 9330783020218



Handheld SpO₂/HR Monitor

- Battery charge function with adapter or charge stand (optional)
- High resolution 2.8" TFT display
- Adjustable audible and visual alarms
- Protective cover perfect for portable monitoring needs
- Patient ID ranges from 1 to 127
- Auto power-off: Selectable auto power-off time - from 1 minute to 15 minutes

PH-MD300-M Handheld Pulse Oximeter

GTIN: 9330783033164

PH-MD300-MV Handheld Pulse Oximeter

GTIN: 9330783040445

Supplied with Reusable Ear Clip Sensor
(Suitable for Veterinary Applications)



'OxyWatch' Adult SpO₂ Finger Pulse Oximeter

- Highly Accurate
- Special design for children
- Display SpO₂ waveform (pleth or line), SpO₂, PR, bar graph
- LCD display screen with six display modes
- Automatic power-off design

PH-MD300C63 Finger Pulse Oximeter

GTIN: 9330783031542



Paediatric SpO₂ Finger Pulse Oximeter

- Highly Accurate
- Special design for children
- Display SpO₂ waveform (pleth or line), SpO₂, PR, bar graph
- LCD display screen with six display modes
- Automatic power-off design

PH-MD300-C5 Finger Pulse Oximeter

GTIN: 9330783024513

RESPIRATORY CARE PRODUCTS



SpO₂ Sensors

Disposable SpO₂ Sensors



M-50J-A Disposable Sensor (<3kg or >30kgm)
Position: Neonate - Lateral Side of Foot OR
Adult - Index Finger
Material: Medaplast

GTIN: 9330783033188



M-50J-B Disposable Sensor (<3kg or >30kgm)
Position: Neonate - Lateral Side of Foot OR
Adult - Index Finger
Material: Microfoam

GTIN: 9330783033195



M-50J-C Adult Disposable Sensor (>30kgm)
Position: Index Finger
Material: Microfoam

GTIN: 9330783033201

Reusable SpO₂ Sensors



M-50A Reusable Finger Clip Sensor (20-60kgm)
Suit: Small Adult
Position: Index Finger

GTIN: 9330783033003



M-50B Paediatric Reusable Sensor (10-50kgm)
Suit: Paediatric
Position: Index Finger

GTIN: 9330783033010



M-50C Neonate Reusable Wrap Sensor (>10gm)
Suit: 1 - 4kg Position: Lateral Side of Foot
Suit: 3 - 5kg Position: Big Toe
Suit: >10g Position: Big Toe / Thumb

GTIN: 9330783033027



M-50D Adult Reusable Ear Clip Sensor (>40kgm)
Suit: To Adult
Position: Ear

GTIN: 9330783033034



M-50E Adult Reusable Finger Clip Sensor (>40kgm)
Suit: To Adult
Position: Index Finger

GTIN: 9330783033041



M-50G Adult Reusable Soft Sensor (>40kgm)
Suit: Adult
Position: Index Finger
Thickness of Finger: 9 - 20mm

GTIN: 9330783033058



M-50H Paediatric Reusable Soft Sensor (>20kgm)
Suit: Paediatric
Position: Index Finger
Thickness of Finger: 7 - 12mm

GTIN: 9330783033065

Apnoea and Sleep Recorders



SpO₂/PR Wrist Pulse Oximeter

- Professional sleep monitoring device
- Wrist measurement, convenient clinic operation
- Compact and lightweight design
- Application from adult to paediatric
- Low power consumption, working with 1 AAA battery
- Data transmission, storage and replay
- 1-10 patient ID setting
- Low voltage indicator
- Downloadable Sleeping Pattern history

PH-MD300-W Wrist SpO₂ Pulse Oximeter

GTIN: 9330783024476

SpO₂/HR Wrist Pulse Oximeter PC Software

Download to PC, storage information from PH-MD300-W for replay and examination

PH-MD300-W/DOS Wrist SpO₂ Pulse Oximeter
Including Software

GTIN: 9330783024469

NEW



SpO₂/PR Wrist Pulse Oximeter

- Compact and Lightweight
- Four Display Modes
- Low Power Consumption
- Battery Low Indicator
- Rechargeable Lithium-Ion Battery
- Visible, Audible and Vibration Alarms
- Automatic Power-off after Finger Out for 8 secs
- Bluetooth® and USB Modes for Data Transmission

PH-MD300-W2 Wrist SpO₂ Pulse Oximeter

GTIN: 9330783040438

Infant Care Products



NeoShades®

Phototherapy Eye Shades with Tabs

- Attaches to patients with gentle NeoBond Hydrocolloid
- No cumbersome straps around the head
- Available in three sizes
- Special shape provides ideal fit on face
- Cute design encourages parents to photograph baby

N720 NeoShades® Phototherapy Eye Shields (Micro) GTIN: 9330783000067

N721 NeoShades® Phototherapy Eye Shields (Small) GTIN: 9330783000050

N722 NeoShades® Phototherapy Eye Shields (Large) GTIN: 9330783000043

Boxes of 50



Non-slip Surface



NeoShades®

Phototherapy Eye Shades with NeoFoam™ Headstrap


- Attaches to patients with gentle NeoFoam™ headstrap with non-slip surface
- Available in three sizes
- Special shape provides ideal fit on face
- Cute design encourages parents to photograph baby

N727 Phototherapy Eye Shields with headstrap (Micro) GTIN: 9330783028092

N728 Phototherapy Eye Shields with headstrap (Small) GTIN: 9330783028108

N729 Phototherapy Eye Shields with headstrap (Large) GTIN: 9330783028085

Boxes of 50



NeoSmile™

Temperature Probe Cover

- Only temperature probe cover that attaches to patients with gentle NeoBond® hydrocolloid
- Resists skin breakdown from moisture or heat · Leaves no residue on skin
- Doesn't fall off prematurely like gel or acrylic covers

N731 NeoSmile™ - Temperature Probe Cover with NeoBond® Hydrocolloid GTIN: 9330783018772



Hold-A-Line securely holds one or two lines in place

Hold-A-Line


Catheter & Tubing Holder

- Easy to open and close repeatedly
- Holds catheters, IV's and various tubes
- Unique shape allows for larger umbilical stumps
- Clear flap allows for visualisation of markings
- Relieves stress on tubes
- Helps prevent tube dislodgement

N705 Hold-A-Line Catheter & Tubing Holder with Neobond® GTIN: 9330783018673

* All items above available in Malaysia only

Infant Care Products



NeoBar® Neonate

ET Tube Holder with gentle hydrocolloid adhesive tabs

- Attaches to patients with latex free NeoBond® hydrocolloid
- Reduces extubations and helps prevent palate trauma
- Available in 6 colour coded sizes for babies as small as 400 gm
- No tape near nose or mouth, allows for better oral care
- Standardises application procedures and prevents frequent tape changes

N709 NeoBar® Nasal Prong CPAP Holder - Mini - Purple GTIN: 9330783018680


N710 NeoBar® Nasal Prong CPAP Holder - Ultra - Yellow GTIN: 9330783018697

N711 NeoBar® Nasal Prong CPAP Holder - Micro - White GTIN: 9330783018703

N712 NeoBar® Nasal Prong CPAP Holder - Small - Green GTIN: 9330783018710

N713 NeoBar® Nasal Prong CPAP Holder - Large - Blue GTIN: 9330783018727

N714 NeoBar® Nasal Prong CPAP Holder - XL - Burgundy GTIN: 9330783032754



NeoBar® Paediatric

ET Tube Holder with Gentle Hydrocolloid or Foam Adhesive Tabs

Gentle Hydrocolloid Tabs

N715 NeoBar® Nasal Prong CPAP Holder - Macro - Clear GTIN: 9330783018734

N716 NeoBar® Nasal Prong CPAP Holder - Jumbo - Grey GTIN: 9330783025855

N717 NeoBar® Nasal Prong CPAP Holder - Mega - Tan GTIN: 9330783032761

Friendly Foam Adhesive Tabs

N715F NeoBar® Nasal Prong CPAP Holder - Macro - Peach GTIN: 9330783032778

N716F NeoBar® Nasal Prong CPAP Holder - Jumbo - Carbon GTIN: 9330783032785

N717F NeoBar® Nasal Prong CPAP Holder - Mega - Brown GTIN: 9330783032792



NeoBridge®

Umbilical Catheter Holder with Secure Positioning Flaps

- Attaches to patients with latex free NeoBond® hydrocolloid
- Helps prevent dislodgement and standardises procedure
- Available in three sizes
- Flaps can be opened and closed repeatedly
- Accommodates two catheters

N700 NeoBridge® Umbilical Catheter Holder - Micro GTIN: 9330783018635

N701 NeoBridge® Umbilical Catheter Holder - Small GTIN: 9330783018642

N702 NeoBridge® Umbilical Catheter Holder - Large GTIN: 9330783018659



N204



N205



N224



N225



N275

Little Sucker®

Oral & Nasal Suction Devices for Gentle and Effective Suctioning

- Single-handed suctioning
- Thumb port for intermittent suctioning
- Soft, flexible tip like a bulb syringe
- Use in NICU, PICU, respiratory, labour & delivery
- Improve patient care and save nursing time
- Replace up to 3 products with 1 little sucker

N204 Oral & Nasal Suction Adaptor - Preemie GTIN: 9330783018437

N205 Oral & Nasal Suction Adaptor - Standard GTIN: 9330783018444

N224 Nasal Suction Adaptor - Preemie Tip GTIN: 9330783036295

N225 Nasal Suction Adaptor - Standard GTIN: 9330783018543

N275 "Snorkel" Multi-purpose Suction Device GTIN: 9330783018581

* All items above available in Malaysia only



Meconium Aspirator®

Meconium Removal Device that Allows Intermittent Suction

- Eliminates mouth suctioning of meconium
- Attaches to any endotracheal tube
- Thumb port allows for intermittent suctioning
- Latex free

N0101-BK Clear Meconium Aspirator

GTIN: 9330783018369



"Little Sucker"
placed in a
Suction Caddy
and suction
capped off.




Suction Caddy™

Suction Device Holder for Little Sucker

- Easily attaches to any solid surface
- Cap off suction and keep things quiet
- Keep your suction device clean and close at hand
- Ideal for holding the "Little Sucker"
- Low price and easy to use

N245 Suction Caddy

GTIN: 9330783018574



NeoBond®


Medical Grade Hydrocolloid Strips

- Use in place of regular adhesive tapes
- Provides long term attachment but prevents epidermal stripping
- Leaves no residue on skin
- No skin prep required
- Latex free
- Ideal for use with limb boards

N741 NeoBond® Medical Grade Hydrocolloid Strips

GTIN: 9330783018789

Size: 12.7mm x 152.4mm (1/2" x 6")



Standard

Mini

EZ-Hold®

Tubing and Cannula Holder

- Hydrocolloid base is gentle on delicate skin
- Fold-over flap reduces stress on tubes
- Ideal for use with nasal cannulas and feeding tubes
- Easy to use
- Open and close repeatedly

N756 Tubing and Cannula Holder - Mini

GTIN: 9330783032648

N757 Tubing and Cannula Holder - Standard

GTIN: 9330783024711



NeoHold™

Cannula & Tubing Holder

- Attaches to patients with gentle NeoBond hydrocolloid
- Relieves stress on tubes and helps prevent dislodgement
- Easy to apply and remove
- Clear flap allows visualisation of tubes


N751 NeoHold™ Cannula & Tubing Holder

GTIN: 9330783018796

Size: 12.7mm x 31.75mm (1/2" x 1-1/4")

* All items above available in Malaysia only

Infant Care Products



Duotube Silicone Binasal Pharyngeal Airways

- Easy to insert and remove
- Reduces Septal and nares trauma
- Made from soft latex free silicone
- Available in 3 sizes
- Available in 2 lengths
- Stays securely in place

20mm Tube length

95-02504	15mm M to 2.5mm Dia. Duotube	(PS)(S)	GTIN: 9330783026074
95-03004	15mm M to 3.0mm Dia. Duotube	(PS)(S)	GTIN: 9330783026098
83-03504	15mm M to 3.5mm Dia. Duotube	(PS)(S)	GTIN: 9330783031108

35mm Tube length

97-02504	15mm M to 2.5mm Dia. Duotube	(PS)(S)	GTIN: 9330783000487
97-03004	15mm M to 3.0mm Dia. Duotube	(PS)(S)	GTIN: 9330783000470
85-03504	15mm M to 3.5mm Dia. Duotube	(PS)(S)	GTIN: 9330783031115



Lead Saver™ Electrode or Temp Probe Cover Securing Device

- Attaches to patients with gentle NeoBond® hydrocolloid
- Prevents gel oozing out edges
- Saves valuable nursing time


N325 Electrode or Temp Probe Cover Securing Device GTIN: 9330783032655



NeoPulse™ Adjustable and Disposable Pulse Oximeter Wrap

- Made with unique, patented Neofoam™
- Soft and ultra skin friendly
- Lightweight and comfortable
- Provides stable, secure attachment that won't slip
- Wicks away perspiration
- Breathable material resists development of foot odour on baby

N781 NeoPulse™ Pulse Oximeter Wrap - Small GTIN: 9330783032808
N782 NeoPulse™ Pulse Oximeter Wrap - Large GTIN: 9330783032099




NeoLead® Pre-wired Neonatal Electrodes

- Pre-wired
- Hydrocolloid border
- Provides long term attachment
- Produces excellent tracings
- Won't dry out prematurely
- Safe on neonatal skin
- Prevents gel from oozing out the edges
- Hydrogel centre
- Radiolucent Neolead - Doesn't show up on x-rays

N301 Standard NeoLead® - Small GTIN: 9330783032815
N305 Radiolucent NeoLead® - Small GTIN: 9330783032822

Boxes of 20 sets



Micro NeoLead® Pre-wired Neonatal Electrodes

- Pre-wired
- Safely secures to limbs without a strap
- No cumbersome band
- Safe on neonatal skin
- Gentle hydrocolloid ends and Hydrogel adhesives
- Radiolucent materials allow for x-ray wherever they are positioned

N306 Radiolucent Micro NeoLead® - Micro GTIN: 9330783032839

Boxes of 20 sets

* All items above available in Malaysia only

Infant Care Products



NeoGrip™

Multi-Purpose Tubing and Cable Holders

- Leaves no residue on tubes, cable or wires
- Holds suction tubes, no bunching or twisting
- Use in any hospital department
- Keep your equipment with you at all times
- Great for use in NICU
- Ideal for use during transit

N602B NeoGrip™ - Tubing and cable holder - Blue

GTIN: 9330783032853

N602P NeoGrip™ - Tubing and cable holder - Pink

GTIN: 9330783032938



	HN - Hawaiian
	BL - Blue
	CO - Camouflage
	PK - Pink

NOTE: Tracheostomy tube not included

EZCare™

Tracheostomy Tube Holder

- Neoprene material is soft, lightweight and comfortable and won't irritate skin
- Provides stable attachment and reduces chance of accidental dislodgement
- Moisture resistant material allows patient to be bathed without removing holder
- Stretchable material allows for cough reflex
- Colourful straps help resist staining

	<i>Strap Size</i>	<i>Neck Size</i>	
N905 EZCare™	127mm	180-205mm	GTIN: 9330783032877
N906 EZCare™	152mm	203-230mm	GTIN: 9330783032884
N907 EZCare™	178mm	230-254mm	GTIN: 9330783032891
N908 EZCare™	203mm	254-280mm	GTIN: 9330783032907
N909 EZCare™	230mm	280-305mm	GTIN: 9330783032914
N910 EZCare™	254mm	305-330mm	GTIN: 9330783032921



N125



NeoFlex Silicone Gel Adhesive Tape

- Safe on Fragile Skin
- Flexible and Easy to Use
- Secure lines, tubes, IV boards and more
- Ideal for Micro Preemies and Sensitive Skin
- Clear Material Allows Skin Monitoring

N125 Silicone Gel Adhesive Tape

GTIN: 9330783034925

1 Roll per Box (36" (3 feet) / 915mm)



Silicone Nipple Shields

- Easy to apply and use
- Soft silicone material
- Reduce pain from cracked nipples
- Latex free

51410 Nipple Shields (2 per pack) (S)

GTIN: 9330783006236



MFBB.RENS.60



MFBB.FLNS.60

Comepa Neonatal ECG Electrodes

NOW with NEW & IMPROVED water-proof backing - improving adhesion in high humidity conditions.

MFBB.RENS.60 Comepa Pregelled Rectangular Electrode GTIN: 9330783017973

MFBB.FLNS.60 Comepa Pregelled Flower Electrode GTIN: 9330783017959

Legend: (S) = Silicone (PS) = Polysulfone (Foam) = Polyethylene/Polyester (Nylon) = Nylon

* All items above available in Malaysia only

Infant Care Products

BabyIVBoards®

Neonatal, Infant, Paediatric & Adult IV Arm Boards

- Soft and pliable
- Reduce discomfort to limb
- Easy to apply and remove
- Available in 4 sizes
- Aluminium insert to help keep profile

APM602 Neonatal IV Armboards with flexible aluminum insert

Size: 105mm x 30mm (Carry Pack of 50) (Al)(Foam) GTIN: 9330783000029

APM603 Infant IV Armboards with flexible aluminum insert


Size: 130mm x 30mm (Carry Pack of 30) (Al)(Foam) GTIN: 9330783000036

APM604 Paediatric IV Armboards with flexible aluminum insert

Size: 175mm x 50mm (Carry Pack of 15) (Al)(Foam) GTIN: 9330783000081

APM605 Adult IV Armboards with flexible aluminum insert

Size: 250mm x 80mm (Pack of 20) (Al)(Foam) GTIN: 9330783025565



Neonatal SpO₂ Probe Cover / Securing Strap

- Easy to apply and remove
- Soft and gentle
- Velcro™ tie allows repeated opening and closing

B905 Neonatal SPO₂ Probe Cover with Velcro™ tie


(for use with SPO₂ sensor) - 35mm x 135mm (Foam)(Nylon) GTIN: 9330783009732

B907 Neonatal Securing Strap with Velcro™ tie

- suits APM602 & APM603 IV Armboards - 150mm x 15mm (Foam)(Nylon) GTIN: 9330783025978

B908 Securing Strap with Velcro™ tie

- suits APM604 IV Armboards - 210mm x 25mm (Foam)(Nylon) GTIN: 9330783025978



TB2000 Infant Resuscitation Circuit

Infant Resuscitation Circuit


with Adjustable PEEP Flow Controller Function

- Fits NeoPUFF, NeoPIP & NeoGuard Units
- Easy to Use
- Adjustable PEEP Flow Controller
- Double-swivel at Mask and Tubing Ends
- PEEP Controller Indexer

TB2000 Infant Resuscitation Circuit

GTIN: 9330783036929


with Adjustable PEEP Flow Controller Function (PP) (PE) (PC)



Adjustable PEEP Function

Legend: (S) = Silicone (PP) = Polypropylene (PVC) = Polyvinyl Chloride (PC) = Polycarbonate (PE) = Polyethylene

Infant Care Products



BabyIVBoards®

Neonatal, Infant, Paediatric & Adult IV Arm Boards

- Soft and pliable
- Reduce discomfort to limb
- Easy to apply and remove
- Available in 4 sizes
- Aluminium insert to help keep profile

APM602 Neonatal IV Armboards with flexible aluminum insert

Size: 105mm x 30mm (Carry Pack of 50) (Al)(Foam) GTIN: 9330783000029

APM603 Infant IV Armboards with flexible aluminum insert

Size: 130mm x 30mm (Carry Pack of 30) (Al)(Foam) GTIN: 9330783000036

APM604 Paediatric IV Armboards with flexible aluminum insert

Size: 175mm x 50mm (Carry Pack of 15) (Al)(Foam) GTIN: 9330783000081

APM605 Adult IV Armboards with flexible aluminum insert

Size: 250mm x 80mm (Pack of 20) (Al)(Foam) GTIN: 9330783025565

Neonatal SpO₂ Probe Cover / Securing Strap

- Easy to apply and remove
- Soft and gentle
- Velcro™ tie allows repeated opening and closing

B905 Neonatal SPO₂ Probe Cover with Velcro™ tie

GTIN: 9330783009732

(for use with SPO₂ sensor) - 35mm x 135mm (Foam)(Nylon)

B907 Neonatal Securing Strap with Velcro™ tie

GTIN: 9330783025978

- suits APM602 & APM603 IV Armboards - 150mm x 15mm (Foam)(Nylon)

B908 Securing Strap with Velcro™ tie

GTIN: 9330783025978

- suits APM604 IV Armboards - 210mm x 25mm (Foam)(Nylon)



OH3000



Oxydome



Oxypod

Oxygen Headboxes

OH3000 (Oxypod II)

GTIN: 9330783019175

The OH3000 is to be used exclusively as a headbox. The reduced length enables even big babies to use it as a headbox within the restricted space of an incubator. (PC)

NOD-101D (Oxydome)

GTIN: 9330783018970

The NOD-101D is injection moulded from "Food Grade" polycarbonate and has no joints to break or corners to harbour bacteria and defy proper cleaning. The Silicone rubber flap acts as a soft protective buffer minimising chafing and allows tighter control of Oxygen concentrations. (PC)

NOD-102D (Oxydome II)

GTIN: 9330783018987

For larger babies, the NOD-102D is used as a conventional headbox having a larger head aperture to accommodate the shoulders. (PC)

NOD-112D (Oxypod)

GTIN: 9330783019076

The NOD-112DAP is a body box version with a unique, radically shaped Oxygen hood for newborns. Injection moulded from "Food Grade" polycarbonate, it is immensely strong, comfortable for the babies and easy to use and clean. Access is granted through large sliding doors and access slots. (PC)

Legend: (Al) = Aluminium (PC) = Polycarbonate (Foam) = Polyethylene/Polyester

Electric Scooter Range

Extra Large 4 Wheel Scooter with Hard Canopy*



Available in Red or Blue

Features:

- Adjustable Tiller
- 150kg Weight Capacity
- Maximum Speed 10km/hr
- Front and Rear Suspension
- Steering Wheel with LCD Display
- Two Head Lights
- Per Charge Range = 45km (with 75A Battery)
- Seat with Slider
- Double Seating
- Left and Right Steering Wheel
- Front Bumper Design / Trunk
- Hard Canopy / Windshield
- Rearview Mirror

SE4010 Extra Large 4 Wheel Scooter with Hard Canopy* GTIN: 9330783029082

Extra Large 4 Wheel Scooter with Hard Canopy*



Available in Red or Blue

Features:

- Adjustable Tiller
- 180kg Weight Capacity
- Maximum Speed 17km/hr
- Front and Rear Suspension
- Luxurious Seat and Footrest
- Two Head Lights
- Per Charge Range = 42km (with 75A Battery)
- Front Basket
- Rearview Mirror
- Seat with Slider
- Luxurious Hard Roof / Canopy
- Hard Basket (optional accessory)
- Wide Tyre - for a stable and comfortable drive

SX4040 Extra Large 4 Wheel Scooter with Hard Canopy* GTIN: 9330783029068

Extra Large 3 Wheel Scooter with Hard Canopy*



Available in Red or Blue

Features:

- Adjustable Tiller with Hand Brake
- 180kg Weight Capacity
- Maximum Speed 17km/hr
- Luxurious Seat and Arm Rest
- Seat with Slider
- Hard Canopy*
- 360° Rotational Seat - for easy Ascend and Decent
- One Head Light
- Trailer (optional accessory)
- Per Charge Range = 37km (with 75A Battery)

SX3040 Extra Large 3 Wheel Scooter with Hard Canopy* GTIN: 9330783029075

*Canopy not shown

Electric Scooter Range



Large 4 Wheel Scooter with Closed Canopy

Features:

- Adjustable Tiller
- 150kg Weight Capacity
- Maximum Speed 10km/hr
- Front and Rear Suspension
- Steering Wheel with LCD Display
- Per Charge Range = 40km (with 75A Battery)
- Seat with Slider
- Disassemble Doors
- Front Bumper Design
- Storage Bag inside Scooter
- Rearview Mirrors

SC4050 Large 4 Wheel Scooter with Closed Canopy GTIN: 9330783042012



Large 4 Wheel Electric Chair

Features:

- Adjustable Tiller
- 180kg Weight Capacity
- Per Charge Range = 41km
- Easy to Use Control Panel
- Windshield on both sides
- Armrest Width is Adjustable
- Foldable Foot Rest
- Easy to Recline Chair

WL4024 Large 4 Wheel Electronic Wheel Chair GTIN: 9330783029051

Available in Red



Mini Size 4 Wheel Foldable Scooter

Easy as 1, 2, 3

Features:

- Space Saving Foldable Design
- Easy to Carry as Luggage
- Tiller Angle Adjustable
- Automatically Retrievable Anti-tipping Wheel
- Accelerator
- Foldable Foot Rest
- Max. Speed 6kph

SP4010 Mini Size 4 Wheel Foldable Scooter GTIN: 9330783042005



Electric Powered Bicycles

- 2 Sizes Available - 508mm & 700mm
- Attain Speeds of Up to 25kpm
- BM2030 - Up to 35km per Charge
- BM2040 - Up to 30km per Charge
- 200W Motor

BM2030 Elec. Powered Bicycle (508mm Wheel Size) GTIN: 9330783043996



BM2040 Elec. Powered Bicycle (700mm Wheel Size) GTIN: 9330783044009

Index

A			E (cont.)	
Adaptors / Connectors - Disposable	6, 7, 8, 9		Electric Bicycles	99
Adaptors / Connectors - Light Weight Metal, Reusable	10, 11		Electric Scooters	98, 99
Adaptors / Connectors - Polycarbonate Reusable	3, 4, 5		Electric Wheelchairs	99
Adaptors / Connectors - Polysulfone Reusable	1, 2		Expiratory Valve - Metal (DF655)	11
Aerosol Therapy Masks	62		Extension Sets (for Biegler Blood / Fluid Warmer)	76
Air Adaptor - 3 Way	84		EZCare	95
Air Flowmeter Nipples	85		EZ-Head	93
Air Flowmeters	79		F	
Air Handwheel Kits	85		Face Tent Masks	62
Air Regulators	80		Filters (Bacterial / Viral)	24
Anaesthetic Circuits - Adult Disposable	25-43		Filters (FHME)	23
Anaesthetic Circuits - Adult Reusable	46		Filters (for Suction Attachment)	82
Anaesthetic Circuits - Coaxial	43		Filters (for Twin O Vac)	85
Anaesthetic Circuits - Paediatric Disposable	33, 34, 35, 42, 43		Filters (HME Tracheostomy)	24
Anaesthetic Circuits - Paediatric Reusable	46		Filters (HME)	24
Analgesia / Entonox / Mid O Gas Circuits - Disposable	47, 48		Finger Pulse Oximeters	88
Analgesic Administration System (Entonox)	84		Flowmeter Nipples (Oxygen / Air)	85
Armboards - Neonate to Adult	96		Flowmeters (Low Flow)	79
Autopress Pressure Infusor	78		Flowmeters (Oxygen / Air)	79
B			Forceps (Magill, Stainless Steel)	73
Bags - Rebreathing	49		G	
Biegler Autopress (Pressure Infusor)	78		Gauge - Manometer	11
Blood / Fluid Warmer	75		Gauge (Oxygen Air, Suction)	85
Blood / Fluid Warmer "Emergency Transport Warmer"	78		Glucose Monitors	87
Blood / Fluid Warmer "Extension Sets"	76		H	
Blood / Fluid Warmer "Protherm II"	77		Handheld Pulse Oximeters	88
Blood / Fluid Warmer "Tube Flow"	75		Handles (Laryngoscope, Standard, Fibre Optic, LED)	72
Blood Pressure & Glucose Monitors	87		Handwheel Kits (Oxygen, Air, Nitrous Oxide)	85
Blood Pressure Monitors	87		Head Harness (for CPAP Masks)	59
Bodok Seals	85		Head Harness (for Masks)	57
Bulbs (Laryngoscopes)	73		Heating Bags (for Protherm II)	77
C			High Concentration Oxygen Therapy Masks	62
Cable Ties (Secur-A-Tag)	74		HME Filters	23, 24
Catheter Mount - Disposable	13		Hold-A-Line	91
Catheter Mount - Reusable	12		Hook Ring (for Masks)	53
Central Vacuum Suction Attachment	81 & 83		Hose Assembly (Oxygen, Air, Nitrous Oxide)	85
Chamber (Humidification)	18		Humidification Chambers	18
Circuits - Coaxial	43		Humidification Heated & Non-Heated Wire Circuits	19, 20, 21
Circuits - Heated & Non-Heated Humidification	19, 20, 21		Humidification Heated Wire Adaptors	22
Circuits Anaesthetic - Adult Disposable	25-43		Humidification Heater Base	18
Circuits Anaesthetic - Adult Reusable	46		Humidification Temperature Probes	22
Circuits Anaesthetic - Paediatric Disposable	33, 34, 35, 42, 43		Humidifier (Oxygen TM11)	84
Circuits Anaesthetic - Paediatric Reusable	46		I	
Circuits Mapleson / Magill	44, 45, 46		Infrared Thermometers	86
Circuits Mid O Gas / Analgesia	47, 48		IV Armboards / Splints (Neonate to Adult)	96
Circuits Resuscitation - Adult, Paediatric, Neonate	44, 45, 46, 96		J	
CO ₂ Nasal Cannulas	14		Jackson Rees Circuits	44, 45, 46
CO ₂ Sample Lines	16		Jars (for Twin O Vac)	83
Comepa Neonatal	95		K	
Connectors - Disposable	6, 7, 8, 9		Keywheel / Spanner (for Oxygen & Air Bottles)	85
Connectors - Light Weight Metal, Reusable	10, 11		L	
Connectors - Polycarbonate Reusable	3, 4, 5		Laryngoscope Blades (Disposable S/S "Standard Conventional" LED)	69
Connectors - Polysulfone Reusable	1, 2		Laryngoscope Blades (Disposable Stainless Steel F/Optic & LED)	68
CPAP Kits	97		Laryngoscope Blades (Duoscope)	66
Cylinder Key & Spanner	85		Laryngoscope Blades (Onescope)	67
D			Laryngoscope Blades (Reusable Mac, Miller, Kessel, Mac Flexible Tip etc.)	70, 71, 72
Diluter (Adjustable)	65		Laryngoscope Bulbs & Spare Parts	73
Diluter Venturi	65		Laryngoscope Handles (Standard, F/Optic, LED)	72
Diverters (for Peep Valves)	60		Lead Saver	94
Duoscope (Laryngoscope)	66		Little Sucker	92
Duotube	94		Low Flowmeters (Oxygen)	79
E				
ECG Electrodes	94, 95			
Emergency Transport Blood / Fluid Warmer	78			
Entonox Analgesic Administration System	84			

Index

M

Magill Forceps (Stainless Steel)	73
Manometer Gauge	11
Mapleson / Magill Circuits	44, 45, 46
Masks - Aerosol Therapy	62
Masks - CPAP Non-Vented (Full Face & Nasal)	58, 59
Masks - Face Tent	62
Masks - Oxygen Therapy & High Concentration	61
Masks - Oxygen with Concentration Jets	65
Masks - Tracheostomy Masks & Kits	64
Masks - Venturi	65
Masks (Disposable MRI Compatible with Inflate / Deflate port)	56
Masks (Disposable SCENTED with Inflate / Deflate port)	55
Masks (Disposable Twin Port)	57
Masks (Disposable with Inflate / Deflate port)	54
Masks (PVC Disposable)	53
Masks (Reusable Silicone, Adult to Small Infant)	53
MDI Adaptors	8
Meconium Aspirator	93
Micro NeoLead	94
Mid O Gas / Analgesia Circuits	47, 48
Monitors (Blood Pressure and Glucose)	87
Monitors (Handheld & Finger Pulse Oximeters)	88

N

Nasal Cannula CO ₂ / Oxygen (Adult - Paediatric)	14
Nasal Cannula Diagnostic (Adult - Paediatric)	15
Nasal Cannula Oxygen (Adult - Neonate)	14, 15
Nebulizer Bottles	63
Nebulizer Kit	65
NeoBar	92
NeoBond	93
NeoBridge	92
NeoFlex Silicone Gel Adhesive Tape	95
NeoGrip	95
NeoHold	93
NeoLead	94
NeoPulse	94
NeoShades	91
NeoSmile	91
Nipples (for Oxygen & Air Flowmeters)	85
Nitrous Oxide 3 Way Adaptor	84
Nitrous Oxide Handwheel Kits	85
Non-Contact Thermometers	86
Norman Mask Elbow (Disposable)	9
Norman Mask Elbow (Reusable)	10

O

Onescope (Laryngoscope Blade)	67
Otoscope	73
Oxygen / CO ₂ Nasal Cannula (Adult - Paediatric)	14
Oxygen Adaptor - 3 Way	84
Oxygen Flowmeter Nipples	85
Oxygen Flowmeters	79
Oxygen Handwheel Kits	85
Oxygen Hose Assembly	85
Oxygen Humidifier	84
Oxygen Nasal Cannula (Adult - Neonate)	14, 15
Oxygen Nasal Cannula Diagnostic (Adult - Paediatric)	15
Oxygen Regulators	80
Oxygen Therapy Masks	61
Oxygen Therapy Masks with Concentration Jets	65
Oxygen Tubing - Star Lumen Disposable	16

PQ

Peep Valve Diverters (Reusable & Disposable)	60
Peep Valves (Reusable & Disposable)	60
Probe Cover (SpO ₂ Neonate)	35
Protherm II (Blood & Infusion Warmer)	77
Pulse Oximeter Sensors	89

PQ (cont.)

Pulse Oximeters (Handheld, Finger and Wrist Style)	88, 90
--	--------

R

Re-Breathing Bags	49
Regulators (Oxygen, Air, Nitrous Oxide, Entonox)	80
Regulators (Suction)	81
Regulators / Flowmeters (Combined Unit)	80
Resuscitation Circuits (Adult / Paediatric / Neonate)	44, 45, 46, 96

S

Samples Lines - CO ₂	16
Scooters (Electric)	98, 99
Seals 'Bodok'	85
Securitag (Security Cable Ties)	74
Securing Straps (SpO ₂)	96
Sensors (SpO ₂)	89
Silicone Nipple Shields	95
Smoothbor Tubing (Reusable, 10mmID, 15mmID, 22mmID)	52
Splint (IV Armboard - Neonate to Adult)	74
Suction (Twin O Vac High & Low Suction)	83
Suction 3 Way Adaptor	84
Suction Attachments (Central Vacuum / Venturi)	81
Suction Caddy	93
Suction Filter Housing	82
Suction Filters and Kits	82
Suction Regulators	81
Suction Trap Jar	83
Suction Tubing	17
Suction Yankauer Handpiece	17

T

T Adaptors (Oxygen, Air, Nitrous, Suction)	84
T Piece Adaptor / Connector (Disposable)	7
T Piece Adaptor / Connector (Reusable)	5
Teeth Opener	57
Temperature probes (Humidification)	22
Thermometers (Tympanic and No Contact)	86
Tracheostomy Filters	24
Tracheostomy Masks & Kits	64
Tube Flow (for Biegler BW685 S Blood Warmer)	75
Tubing - Disposable Corrugated (10mmID, 15mmID, 22mmID)	50
Tubing - Disposable Expandable (15mmID & 22mmID)	51
Tubing - Extension Sets (for Biegler Blood Warmers)	76
Tubing - Oxygen Star Lumen	16
Tubing - Oxygen, Air, Nitrous (for Blenders & Ventilators)	85
Tubing - Reusable (Catheter Mount Silicone Tube)	5
Tubing - Smoothbor (Reusable 10mmID, 15mmID, 22mmID)	52
Tubing - Suction	17
Twin O Vac - Disposable Jar	83
Twin O Vac - Filters	85
Twin O Vac - High and Low Suction	83

UV

Venturi Diluters & Adjustable	65
Venturi Masks	65
Venturi Suction Attachments	81

W

Washington T Piece (Disposable)	9
Washington T Piece (Reusable)	10
Water Traps	22
Wrist Pulse Oximeters	90

XYZ

Y Adaptors (Oxygen and Air)	84
Y Piece (Disposable)	6
Y Piece (Reusable Polycarbonate)	3, 5
Y Piece (Reusable Polysulfone)	1
Yankauer Suction Handpiece	17



MEDIPLAS RESPIRATORY PRODUCTS SDN BHD

QUALITY POLICY

Mediplas Respiratory Products Sdn Bhd., understands the importance of good healthcare. To support this Mediplas is dedicated to supplying the latest innovative and high quality products to Hospital and Healthcare services worldwide.

The Company is committed to implementing a philosophy together with a system where meticulous quality monitoring is a normal practice at all stages of operation.

As a Company, Mediplas ensures that the relationship between customer and supplier is maintained and that we work together towards continual improvement in the quality of care for Adults through to Neonates. By establishing open communication with our customers we can understand their needs and expectations. This enables Mediplas to develop techniques and equipment to address the changing requirements of our customers, worldwide and promotes the continual improvement of our products, services and Quality Management systems.

Quality objectives are identified consistently throughout our production systems and the Management of Mediplas is committed to monitoring, reviewing and improving the objectives regularly and in a timely manner.

Mediplas Respiratory Products Sdn Bhd is totally committed to operating a quality system throughout the company and fully complies to: ISO 13485.

Managing Director

A handwritten signature in black ink, appearing to read "Ian Bird", is written over a white background.

Mr. Ian Bird

2nd January 2011



TERMS AND CONDITIONS



MEDIPLAS RESPIRATORY PRODUCTS SDN BHD

Head Office & International Department:

No7, Jalan KIP 3, Taman Perindustrian KIP, 52200 Kuala Lumpur, Malaysia

Phone: +603 6272 7625 Facsimile: +603 6273 6452

E-mail: sales@mediplas.com.my

TERMS AND CONDITIONS OF SALE

Mediplas Respiratory Products Sdn Bhd (798640-U) hereafter known as "the Company".

1. All products are supplied at ex-Factory Malaysia prices (unless otherwise agreed in writing).
2. All risks in the products shall pass to the "Purchaser" when the product leaves the "Company" store and from that point, regardless of how delivery is affected, any insurance shall be the purchaser's responsibility. Please state Insurance Open Policy number on confirmation order.
3. Written quotations are firm for 30 days from date of writing except for material changes in exchange rates or increases in duty on imported products.
4. Verbal prices quoted are subject to change without notice and are recommended prices only.
5. All orders must be confirmed in writing.
6. All claims for short delivery will only be recognized if made in writing 48 hours from date of delivery.
7. The "Company" takes no responsibility for incorrect freight methods when freight details are not clearly advised at time of ordering and marked on the confirmation order. Confirmation orders must be received no later than one week after verbal orders are placed.
8. Payment method is to be either; Payment in advance or Irrevocable Letter of Credit, made payable to Mediplas Respiratory Products Sdn Bhd.
9. Title to the products shall not pass to the "Purchaser" until the full invoice price has been received by the "Company".
10. The "Purchaser" shall hold the products solely as bailee, returnable on demand, pending passing of title.
11. A Goods Return Authorisation Number must be requested from the "Company" prior to any item being returned for any reason. This number MUST accompany all documentation and all items returned to the "Company" Stores. The "Company" reserves the right to refuse delivery acceptance of any items returned not bearing this authorisation.
12. A 20% minimum restocking charge will be made on all goods returned. Freight will be the responsibility of the Distributor or customer returning goods. Goods will not be accepted for return when held for more than one (1) month from date of sale.
13. Goods brought into supply special requests or requirements will not be accepted for return under any circumstances.
14. An official Mediplas number must be quoted for payment of any repairs carried out by our Distributors.
15. Disposable or Single use products are sold in box quantities.
16. Insurance is the purchaser's responsibility.



Mediplas Respiratory Products Sdn Bhd
No 7 Jalan KIP 3, Taman Perindustrian KIP, 52200 Kuala Lumpur, Malaysia
Ph: +60 3 6272 7625 Fax: +60 3 6273 6452
E-mail: sales@mediplas.com.my Web: www.mediplas.com.my



World Wide Sales Network

Some products nominated in this catalogue may not be available to other countries due to distribution agreements.

www.mediplas.com.my



MEDIPLAS RESPIRATORY PRODUCTS SDN BHD

International Office & Manufacturing Division:

No 7 Jalan KIP 3, Taman Perindustrian KIP,

52200 Kuala Lumpur, Malaysia

Ph: +60 3 6272 7625 Fax: +60 3 6273 6452

E-mail: sales@mediplas.com.my