


O₂ & CO₂ Nasal Cannulas


With New Softi Smoothflow[®] Cannula - "For Total Patient Comfort"



"Softi" Smoothflow[®] Oxygen Nasal Curved Tips Cannula

- Contoured "Softi" Smoothflow[®] cannula provides total patient comfort
- Star lumen tubing
- Over the ear lit secures position on the patient

Part No.	Cannula Size	Tube Length
PH1161	Adult	2.1 metre
PH1163	Paediatric	2.1 metre



"Softi" Smoothflow[®] Oxygen / CO₂ Sampling Nasal Cannula

- Contoured "Softi" Smoothflow[®] cannula provides total patient comfort
- The design allows delivery of oxygen to one nasal prong and sampling of expired CO₂ gases from the other prong simultaneously
- Four (4) CO₂ connections to select from
- Star lumen tubing
- Over the ear lit secures position on the patient

Part No	Cannula Size	Tube Length	Connection
PHC02-261	Adult	2.1 m	Male Luer
PHC02-261-4	Adult	4.0 m	Male Luer
PHC02-263	Paediatric	2.1 m	Male Luer
PHC02-263-4	Paediatric	4.0 m	Male Luer
PHC02-261/F	Adult	2.1 m	Female Luer
PHC02-261-4/F	Adult	4.0 m	Female Luer
PHC02-263/F	Paed.	2.1 m	Female Luer
PHC02-263-4/F	Paediatric	4.0 m	Female Luer
PHC02-261-ZMF	Adult	2.1 m	Micro Filter
PHC02-261-4-ZMF	Adult	4.0 m	Micro Filter
PHC02-263-ZMF	Paediatric	2.1 m	Micro Filter
PHC02-263-4-ZMF	Paediatric	4.0 m	Micro Filter
PHC02-261-ZMFO	Adult	2.1 m	Micro Filter and Orange Female Luer
PHC02-263-ZMFO	Paediatric	2.1 m	Micro Filter and Orange Female Luer

"Softi" Smoothflow[®] Cannula

Diagnostic Cannula & CO₂ Sample Lines



"Softi" Smoothflow[®] Diagnostic Cannula (PVC)

- Can be used in conjunction with pressure transducers or thermistors in Polysomnography (Sleep Studies)
- Optimizes pressure/snore waveform signals
- Inexpensive
- Contoured "Softi" Smoothflow[®] cannula provides total patient comfort
- Over the ear lit secures position on the patient
- Available in various lengths
- Male and Female luer lock connection

Part No.	50 75 100 130 150 200 Length (cm)	ML Male Luer FL Female Luer Connection	61 Adult 63 Paediatric Size
LD	##	##	##

Example Part Number : **LD50ML63** Paediatric Diagnostic Cannula with Male Luer Lock, Length: 50cm



CO₂ Sample Lines with Filter

- Made of medical grade PVC material
- Tube length 3 meters, customization is available
- Sampling lines for measurement of end-tidal CO₂
- The filter protects monitor from moisture and contamination
- Male & Female Luer locks are compatible with most of the monitoring devices on the market.

Part No.
MP180003 3.0 m (10ft.) CO₂

Sample Line with Filter M/F Luer connection

CO₂ Sample Lines with Inner "Star Lumen" Design

Male to Female Luer Connection

Part No.
MP0300 3.0 m (10ft.) CO₂ Sample Line M/F Luer connection
MP0315 4.5 m (15ft.) CO₂ Sample Line M/F Luer connection
MP0330 9.0 m (30ft.) CO₂ Sample Line M/F Luer connection

Male to Male Luer Connection

Part No.
MP0500 3.0 m (10ft.) CO₂ Sample Line M/M Luer connection
MP0515 4.5 m (15ft.) CO₂ Sample Line M/M Luer connection
MP0530 9.0 m (30ft.) CO₂ Sample Line M/M Luer connection



* Product subject to change without notice



Mediplas Respiratory Products Sdn Bhd
No. 7 Jalan KIP 3, Taman Perindustrian KIP 52200 Kuala Lumpur, Malaysia
Phone: +60 3 6272 7625 Fax: +60 3 6273 6452
E-mail: sales@mediplas.com.my

EC REP Obelis s.a
Boulevard Général Wahis 53, 1030 Brussels, BELGIUM
Phone: +(32) 2 732 5954 Fax: +(32) 2 732 6003 E-Mail: mail@obelis.net

